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EMERADE PRE-FILLED PEN: FAILURE TO ACTIVATE

- **Complaints that Emerade pens (autoinjector) have failed to activate have been received.**
- **If Emerade fails to activate, an additional attempt should be made using an increased force when pressing the pen against the intended injection site**
- **Patients should, in line with the SmPC, be prescribed two Emerade pens which they should always carry.**
- **HCPs should inform patients of this advice and remind of the need to have two Emerade pens always available.**

Background

Emerade 150/300/500 micrograms solutions for injection are approved for treating severe acute allergic reactions (anaphylaxis). It is distributed in two presentations: Single package with one auto-injector (pen) and twin package with two auto-injectors (pens).

Complaints of Emerade pens (auto injector), which have failed to activate, have been received. In currently available data from laboratory testing, 2 out of 625 pens failed to activate initially, and required higher than normal activation force for activation.

Investigation of reports of activation failure is ongoing.

Advice to patients

<p>HAS MY EMERADE PEN ACTIVATED? <u>Fig. 1</u></p>  <p>When Emerade has been used the needle shield will extend and lock</p>	<p>Instructions:</p> <ol style="list-style-type: none">1. After using an Emerade pen, verify that the pen has activated.2. An Emerade pen that has been activated, will have an extended needle shield (Fig. 1)
<p>WHAT DO I DO IF MY EMERADE PEN HAS NOT ACTIVATED? <u>FIG. 2</u></p>  <p>Not Activated</p>	<ol style="list-style-type: none">1. If the needle shield has not extended, the pen has not activated (Fig. 2).2. Perform an additional attempt to activate the pen, using an increased force when pressing the pen against the intended injection site. If not successful, proceed to use your second pen.3. Retain suspected un-activated pen (s) for reporting and investigation purposes.

1. You should call emergency services on 112, ask for an ambulance and state “Anaphylaxis”, even if you start to feel better.
2. If you still feel unwell after the first injection, use your second pen 5 to 15 minutes after the first.
3. Retain any pen that is suspected of non-functioning for reporting and investigation purposes.

For more information, please refer to the PIL.

Call for reporting and Company contact point

The reporting of suspected adverse reactions after approval is of great importance. It allows continuous monitoring of the benefit-risk balance of the drug. Healthcare professionals are encouraged to report any suspected adverse event to www.meldenbivirkning.dk