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1. "AAPM code of practice for radiotherapy accelerators: Report of AAPM Radiation Therapy Committee Task Group." Medical Physics, Vol. 21, No. 7, July 1994
2. "Comprehensive Quality Assurance for radiation oncology: Report of AAPM Radiation Therapy Committee Task Group," Medical Physics, Vol. 21, No. 4, April 1994

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¹ Klein et al.: "QA of Medical Accelerator: Report of Task Group", Medical Physics, Vol.36, No. 9. september 2009,

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Peter J. Coronado,
Acting Director, Global Regulatory Affairs
on behalf of
Mika Miettinen,
Director, Global Quality Assurance and Product Reliability

24-Jan-2012

Dato

Kontaktoplysninger for Varian Oncology Helpdesk

Telefon: USA og Canada: 1.888.VARIAN5 (888.827.4265)
Europa: +41 41 749 8844

E-mail: Nordamerika: support-america@varian.com
Australien/New Zealand: support-anz@varian.com
Europa: support-emea@varian.com
Sydøstasien: support-sea@varian.com
Kina/Asien: support-china@varian.com
Japan: support-japan@varian.com
Latinamerika: supporte.al@varian.com

Internet: Oncology Systems kundewebsted - www.myvarian.com
Varian Medical Systems offentligt websted - www.varian.com

Report Form Field Safety Corrective Action

Medical Devices Vigilance System (MEDDEV 2.12/1 rev 5)

1. Administrative information	
Destination	
Recipient	
Name of national competent authority The Danish Medicines Agency	
Address of national competent authority Axel Heides Gade 1 DK - 2300 – Kobenhavn Denmark	
Date of this report 2012-02-09	
Reference number assigned by the manufacturer CP-06611	
Incidence reference number and name of the coordinating National Competent Authority (if applicable)	
Identify to what other National Competent Authorities this report was also sent Austria, Belgium, Bulgaria, Cyprus, Czech Republic, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Lithuania, Luxembourg, Netherlands, Norway, Poland, Portugal, Romania, Slovakia, Slovenia, Spain, Sweden, Switzerland, Turkey, United Kingdom	

2 Information on submitter of the report
Status of submitter
<input checked="" type="checkbox"/> Manufacturer
<input type="checkbox"/> Authorised representative within EEA
<input type="checkbox"/> Others (identify the role):

3 Manufacturer information	
Manufacturer name Varian Medical Systems, Inc.	
Manufacturer's contact person Peter J. Coronado	
Address 3100 Hansen Way	
Postal code 94304	City Palo Alto, California
Phone 1.650.424.6320	Fax 1.650.842.5113
E-mail peter.coronado@varian.com	Country USA

4 Authorised representative information
Name of the authorised representative Varian Medical Systems UK, Ltd.
The authorised representative's contact person Per Persson
Address Gatwick Road

ANNEX 4 EUROPEAN FIELD SAFETY CORRECTIVE ACTION REPORT FORM V. 05/07

Postal code RH10 9RG	City Crawley
Phone +44 (0)1293 601200	Fax +44 (0)1293 533509
E-mail per.persson@varian.com	Country United Kingdom

5 National contact point information	
National contact point name	
Name of the contact person	
Address	
Postal code	City
Phone	Fax
E-mail	Country

6 Medical device information	
Class	
<input type="checkbox"/> AIMD Active implants	
<input type="checkbox"/> MDD Class III	<input type="checkbox"/> IVD Annex II List A
<input checked="" type="checkbox"/> MDD Class IIb	<input type="checkbox"/> IVD Annex II List B
<input type="checkbox"/> MDD Class IIa	<input type="checkbox"/> IVD Devices for self-testing
<input type="checkbox"/> MDD Class I	<input type="checkbox"/> IVD General
Nomenclature system (preferable GMDN) GMDN	
Nomenclature code 35159	
Nomenclature text Accelerator system, linear	
Commercial name/brand name/make C-Series Clinac or Trilogy	
Model number H14, H18, H27, H29	
Serial number(s) and/or lot/batch number(s) See list attached to Field Safety Notice	
Software version number (if applicable) 2.x through 8.x	
Manufacturing date/expiry date (if applicable) N/A	
Accessories/associated device (if applicable) N/A	
Notified body (NB) ID- number BSi / CE0086	

7 Description of FSCA

Background information and reason for the FSCA

Product Indication :

The Varian High Energy Linear Accelerator is a computer controlled photon and electron beam producing, radiation therapy linear accelerator machine. The major components of the High Energy Linear Accelerator platform are: the Stand, the Gantry, the Couch, the Console with Electronics Cabinet, and the Modulator Cabinet. The Varian High Energy Linear Accelerator is intended to provide stereotactic radiosurgery and precision radiotherapy for lesions, tumors and conditions anywhere in the body when radiation treatment is indicated.

Description of Issue:

An event has been reported to Varian involving a C-Series Clinac® accelerator. The event entails excessive connector resistance, which caused the actual jaw positions to differ from the intended jaw positions without warning the operator. Normal morning daily quality assurance testing and light field verification detected the discrepancy. Varian has not received any report of misadministration as a result of the discrepancy. Investigation found that a similar issue could arise on the gantry or couch axes. This does not affect the MLC.

Varian's investigation revealed excessive connector resistance in the lines providing reference voltages to the potentiometers. This can cause a discrepancy between the actual physical positions compared to the corresponding position readouts, leading to an incorrect setting such as treatment field size.

For versions 2.x – 6.x only: C-Series Clinac® version 3.x through 6.x only provide redundant C-Series Clinac® axis position feedback on the axes intended to move during irradiation, these are the upper jaws and the gantry (this redundancy is not included in version 2.x). On these versions, the axes not intended to move during irradiation are not equipped with secondary feedback of the positions. Thus these axes can be susceptible to failures that can affect the position accuracy of the single feedback source. These single feedback axes are the lower jaws, collimator rotation, and all couch axes.

Three reports has been received related to this issue (CP-00214 – Germany; CP-06495 - USA; CP-06611 - USA). In the reported cases, no misadministration or serious deterioration in the state of health was reported.

Note: There are 2 FSNs associated with this FSCA. The FSN labelled CP-06611 pertains to customers who have v7.x-8.x. The FSN labelled CP-06611v2-6 pertain to customers who have v2.x-6.x. Each FSN has a different consignee list and will therefore be distributed only to customers who are affected.

Root Cause:

Design deficiency: An unknown source of resistance was applied to the connection from the motor interface PCB to the potentiometers causing a discrepancy between the physical vs. readout positions of the gantry, collimator, jaw and couch. The wires carrying both +/- voltage reference are used for both primary and secondary position readout indicators, therefore a change in resistance to the connector would result in the same offset in both readouts.

Risk Analysis:

In this situation, when the reference voltage drops the actual jaw position may differ from the intended jaw position; the digital display will indicate only the intended jaw position and will not accurately show the actual jaw position. Any jaw could be affected and could result in unplanned dose to normal tissue and/or critical structures.

Additionally, the gantry position, collimator rotation and couch axes could be affected and the actual position may differ from the intended position.

The most serious probable scenario is a cranial-spinal treatment plan that would be using a

match line shift to 'feather' the hot and cold spots throughout the treatment course. This type of treatment does not always include additional collimation (MLC or a block) to the area that is being matched. If the physical jaw was not in the correct position, an under dose or over dose would occur. This error could gradually increase due to degradation and cause an increase in harm. It is reasonable to expect that a gradual error would be detected. However, if a sporadic error occurred, the discrepancy may be larger and may not be detected through normal QA processes.

Alternatively, if the light field was used to abut the fields and the field size was correct before performing a double exposure image, a subsequent auto go to the planned field size could result in an incorrect field size used for treatment.

Any organ at risk or critical structure that has reached the associated radiation tolerance limit could cause harm and/or late radiation effects. The degree of harm is based on the degree of error and the number of times the field size is incorrect during the course of treatment.

A risk hazard assessment of the consequences leading to the most probable scenario that could be expected if the anomaly occurred is Critical (Could result in permanent impairment or life-threatening injury). The Probability that harm will come to a patient as a result of this issue is Remote (Unlikely to occur but possible during the life of an individual item; unlikely but reasonably expected to occur in the installed base). The risk hazard concluded with an ALARP.

Description and justification of the action (corrective/preventive)

Varian had decided to further mitigate this issue as a Field Safety Corrective Action and will notify all affected Users with a Field Safety Notice. The FSN provides a description of the problem, and advises Users of the steps Varian is taking to address the issue.

For version 2.x – 6.x only: No technical solution is proposed for devices with v2.x to v6.x. The FSN communication serves as the corrective action mechanism for the installed base. These devices are manually driven and monitored by the user and only the y-jaws and the gantry rotation have secondary position readouts (SPROs).

For versions 7.x – 8.x only: Varian is continuing to evaluate possible technical solutions for this issue, but no technical solution is proposed at this time as mitigation. Varian's external corrective action consists of ensuring affected customer sites receive a copy of the FSN.

FSN distribution to affected users will begin immediately. In addition, the notices, including the customer distribution list will be forwarded to Varian Sales, Marketing, and Service organizations. Customer sites will be notified by appropriate methods, and will be requested to return proof of notification (via completed "Proof of Notification Form" or appropriate Postal system certified receipt of delivery) to Varian. The enclosed FSN will also be posted to the customer support site: <http://www.MyVarian.com>.

Advice on actions to be taken by the distributor and the user

Please refer to the FSN for further details

Attached please find

- Field Safety Notice (FSN) in English
- FSN in national language
- Others (please specify):

Time schedule for the implementation of the different actions

90 days to complete FSN distribution and collection of proofs (April 2012). The timeframe for a technical solution (for v7.x-8.x only) to be available has yet to be determined.

These countries within the EEA and Switzerland are affected by this FSCA

Within EEA and Switzerland:

- | | | | | | | | | | |
|--|--|--|--|--|--|--|--|--|--|
| <input checked="" type="checkbox"/> AT | <input checked="" type="checkbox"/> BE | <input checked="" type="checkbox"/> BG | <input checked="" type="checkbox"/> CH | <input checked="" type="checkbox"/> CY | <input checked="" type="checkbox"/> CZ | <input checked="" type="checkbox"/> DE | <input checked="" type="checkbox"/> DK | <input checked="" type="checkbox"/> EE | <input checked="" type="checkbox"/> ES |
| <input checked="" type="checkbox"/> FI | <input checked="" type="checkbox"/> FR | <input checked="" type="checkbox"/> GB | <input checked="" type="checkbox"/> GR | <input checked="" type="checkbox"/> HU | <input checked="" type="checkbox"/> IE | <input checked="" type="checkbox"/> IS | <input checked="" type="checkbox"/> IT | <input type="checkbox"/> LI | <input checked="" type="checkbox"/> LT |
| <input checked="" type="checkbox"/> LU | <input checked="" type="checkbox"/> LV | <input type="checkbox"/> MT | <input checked="" type="checkbox"/> NL | <input checked="" type="checkbox"/> NO | <input checked="" type="checkbox"/> PL | <input checked="" type="checkbox"/> PT | <input checked="" type="checkbox"/> RO | <input checked="" type="checkbox"/> SE | <input checked="" type="checkbox"/> SI |
| <input checked="" type="checkbox"/> SK | | | | | | | | | |

Candidate Countries:
 HR TR

All EEA, Candidate Countries and Switzerland

Others: French Polynesia, Guadeloupe, Martinique, Reunion

These countries outside the EEA and Switzerland are affected by this FSCA
 Algeria, Angola, Argentina, Australia, Azerbaijan, Bahrain, Bangladesh, Belarus, Bolivia, Bosnia and Herzegovina, Brazil, Canada, Chile, China, Columbia, Costa Rica, Dominican Republic, Ecuador, Egypt, Georgia, Guatemala, Guam, Honduras, Hong Kong, India, Indonesia, Iran, Israel, Jamaica, Japan, Kazakhstan, Kenya, Kuwait, Lebanon, Macau, Macedonia, Malaysia, Mauritania, Mexico, Moldova, Morocco, Nepal, New Zealand, Oman, Pakistan, Panama, Paraguay, Peru, Philippines, Puerto Rico, Qatar, Russia, Saudi Arabia, Serbia, Singapore, South Africa, South Korea, Sri Lanka, Syria, Taiwan, Thailand, Trinidad and Tobago, Tunisia, Ukraine, United Arab Emirates, Uruguay, USA, Venezuela, Vietnam, Virgin Islands (American)

8 Comments

List of customer sites in EEA member states, candidate countries and territories.

V 2.x-6.x:

Austria	3	Ireland	1	Romania	5
Belgium	8	Italy	54	Slovakia	5
Czech Republic	6	Latvia	1	Slovenia	1
Denmark	2	Lithuania	2	Spain	37
Finland	9	Martinique	1	Sweden	10
France	61	Netherlands	5	Switzerland	9
Germany	53	Norway	3	Turkey	17
Greece	1	Poland	10	United Kingdom	34
Iceland	1	Portugal	9		

V 7.x-8.x:

Austria	4	Greece	3	Poland	11
Belgium	15	Guadeloupe	1	Portugal	10
Bulgaria	1	Hungary	1	Reunion	2
Cyprus	1	Ireland	6	Romania	2
Czech Republic	9	Italy	40	Slovakia	3
Denmark	6	Latvia	3	Slovenia	1
Estonia	1	Lithuania	3	Spain	31
Finland	11	Luxembourg	1	Sweden	9
France	88	Martinique	1	Switzerland	20
French Polynesia	1	Netherlands	8	Turkey	38
Germany	96	Norway	4	United Kingdom	39

I affirm that the information given above is correct to the best of my knowledge.

 Signature

Peter J. Coronado
 Name

Palo Alto, California, USA
 City

2012-02-09
 Date

Submission of this report does not, in itself, represent a conclusion by the manufacturer and/or authorized representative or the national competent authority that the content of this report is complete or accurate, that the medical device(s) listed failed in any manner and/or that the medical device(s) caused or contributed to the alleged death or deterioration in the state of the health of any person.