

URGENT FIELD SAFETY NOTICE

ATTENTION	Operating Room Director and Materials Management
TYPE OF ACTION	URGENT DEVICE RECALL
REF / DATE	ENSEALG2TS-2012-11 / 27-November-2012
PRODUCT	Ethicon Endo-Surgery (EES) is initiating a voluntary recall of the ENSEAL® G2 Curved and Straight Tissue Sealer devices due to a potential to self-activate
DEVICE DETAILS	The recall involves the following product codes:
	NSLG2C14 ENSEAL 5 mm Diameter Tissue Sealer G2 14 cm Length Curved Jaw
	NSLG2C25 ENSEAL 5 mm Diameter Tissue Sealer G2 25 cm Length Curved Jaw
	NSLG2C35 ENSEAL Laparoscopic 5 mm Diameter Tissue Sealer G2 35 cm Length Curved Jaw
	NSLG2C45 ENSEAL Laparoscopic 5 mm Diameter Tissue Sealer G2 45 cm Length Curved Jaw
	NSLG2S14 ENSEAL 5 mm Diameter Tissue Sealer G2 14 cm Length Straight Jaw
	NSLG2S25 ENSEAL 5 mm Diameter Tissue Sealer G2 25 cm Length Straight Jaw
	NSLG2S35 ENSEAL Laparoscopic 5 mm Diameter Tissue Sealer G2 35 cm Length Straight Jaw
	NSLG2S45 ENSEAL Laparoscopic 5 mm Diameter Tissue Sealer G2 45 cm Length Straight Jaw
	<p>Affected products can be identified by expiration date.</p> <p style="color: red;">NOTE: ALL LOT NUMBERS WITH EXPIRATION DATE FROM OCTOBER 2013 THROUGH SEPTEMBER 2017 ARE IMPACTED.</p> <p>Please use the product identification tool in Attachment A to identify the correct expiration date.</p> <p>The recall <i>does not include</i> the ENSEAL® G2 Super Jaw, ENSEAL® Trio or ENSEAL® Round Tip products.</p> <p>See Attachment B f or a detailed list of product codes that can be used for substitutions.</p>
REASON	<p>There is a potential for self activation while the device is in use, which could result in injury to patients or users (thermal damage). To date, we have not received any reports of patient or user injuries.</p> <p>The root cause has been determined to be damage to the insulation on the activation wire during manufacturing, which in rare cases may lead to self activation during usage.</p>
ACTION	<p>We need your help in ensuring that all affected products are located, accounted for, and returned to [Affiliate Name].</p> <p style="color: red;">EFFECTIVE IMMEDIATELY – DO NOT USE EES ENSEAL® G2 Curved and Straight Tissue Sealers.</p> <ol style="list-style-type: none"> 1) Examine your inventory immediately to determine if you have affected product on hand and remove the affected product. 2) Fill out the Business Reply Form and return it back to [Affiliate Name] within 3 business days, even if you do not have affected product. If you have product to be returned, keep a copy of this form for your records. 3) To return affected product, enclose a copy of the Business Reply Form with the product, and use the pre-paid shipping label to return to:

URGENT FIELD SAFETY NOTICE

	<p>[Affiliate Name / Affiliate Address]</p> <p>Your Sales Representative is available to provide assistance in the completion of this voluntary recall if you should request help.</p>
TRANSMISSION	<p>Please share this information with all of the appropriate staff at your facility and any other organization where the product has been transferred.</p>
CONTACT	<p>[Affiliate Name] will process your product return and issue a credit/replacement upon return of the product and the Business Reply Form.</p> <p>If you have additional questions about this action, please contact your Sales Representative or call [Affiliate Name].</p> <p>We apologize for any inconvenience this will cause you, but rest assured it is our utmost intent to make this process as easy for you as possible.</p>
CONFIRMATION	<p>The Field Safety Action is being conducted with the full knowledge of the U.S. Food and Drug Administration (FDA) and EU National Competent Authorities</p>