Danish Medicines Agency (DMA)

Axel Heides Gade 1, DK - 2300 - Kobenhavn,
DENMARK

Tel:+45 44 88 9595, fax:+45 44 88 9599

E-mail: med-udstyr@dkma.dk,

Web sites: www.medicinskudstyr.dk, <http://medicaldevices.dk/>

Date: 24-Jan-2013

SCOPE: Notification of Field Safety Corrective Action (FSCA) - Resolve Locking Biliary Drainage Catheters

Merit Reference: 1721504-1/25/13-001-R

Dear Sir/Madam,

Merit Medical Systems, Inc. is voluntarily recalling all Resolve Locking Biliary Drainage Catheters (RBC & RBDC). Merit has become aware of one catheter failure during a routine 90-day catheter replacement. The catheter was observed to be in several pieces under fluoroscopic visualization prior to removal. Unused product must be quarantined immediately and returned to Merit. Merit recommends that medical judgment be used in order to assess the risk-benefit of removing catheters from patients. Merit also advises users to discontinue placement of any further devices. Merit has received no reports of patient harm or injury resulting from this issue.

The products involved in the product retrieval are classified under the Medical Device Directive as:

BRAND NAME:	ReSolve® Drainage Catheters	
CLASSIFICATION:	Class IIa sterile, rule 7, according to Annex IX of the Medical Device Directive 93/42/EEC	
CONFORMITY ASSESSMENT ROUTE:	Annex II, Section 3.2 of the EC Directive 93/42/EEC as amended by Directive 2007/47/EC	
NOTIFIED BODY:	BSI Product Services, Number 0086	
EC CERTIFICATE:	CE 541900	
INTENDED USE:	The RBC and RBDC with locking pigtail and hydrophilic coating are used for drainage of bile within the biliary system.	
AUTHORIZED REPRESENTATIVE:	Merit Medical Ireland, Ltd. Parkmore Business Park West Galway, Ireland	Contact: Mark Mullaney Phone: +353-91-703-761 Fax: +353-91-771-888

Root Cause and Corrective Action

- The root cause has not been determined. Corrective action beyond the field action has not been determined at this time.

Removal Actions

- Merit initiated the recall on January 16, 2013. On January 16, sales representatives began contacting each physician and user facility that had been provided product.
- Written communications were drafted and distribution began on 24 January 2013.
- Merit also assessed its on-hand inventory for affected products. The affected products have been quarantined at Merit's facilities.

Complaint Details

- To date, Merit has received one (1) product complaint from a physician in Denmark. There have not been any reports of injury or illness reported with the complaint.
- The complaint stated the following: The user reported that during a routine drainage catheter exchange in the treatment for intrahepatic strictures, the physician noticed that the device had fractured into several pieces under fluoroscopic visualization. A guide wire was inserted and used to straighten the pig tail. The guide wire exited a fracture point of the catheter near the tip. The majority of the device was removed over the wire while holding the suture. The suture near the pigtail was cut which allowed the distal tip of the catheter to remain in the patient. The tip was left in the patient to pass normally through the bowel. The physician visualized the tip of the catheter under fluoroscopy to be moving distally via peristalsis. The device was exchanged without further incident. The patient was discharged as normal. No reports of harm or injury.

Distribution Details:

- See Attachment A.

Field Safety Notice

- Final Field Safety Notice (FSN) proposal is enclosed per Attachment B.

Product Labeling

- In addition, product labeling is enclosed per Attachment C.

It is Merit's objective to finalize recall activities within 30 days. Please do not hesitate to contact me if you have questions or need additional information.

Sincerely,

Merit Medical Ireland Ltd.

Mark Mullaney, Regulatory Affairs Manager

Parkmore Business Park West, Galway, Ireland

Phone #: +353 91 703 761; Fax #: +353 91 771 888; E-mail: mmullaney@merit.com

Attachment A

Distribution Details

Customer	Part #	Lot #	Ship Date	Qty	Qty Used	Quantity Returned	TOTAL
Novotel Cote De Nacre Attn: Anne-Laure Devos Caen, France 14000	RBC-8-038	H399931S1	09/19/12	5	4	36	40
	RBC-10-038	H399938S1	09/19/12	5			
	RBC-12-038	H399940S1	09/19/12	5			
	RBC-14-038	H399942S1	09/19/12	5			
	RBDC-8-038	H399948S1	09/19/12	5			
	RBDC-12-038	H399952S1	09/19/12	5			
	RBDC-10-038	H399950S1	09/20/12	5			
	RBDC-14-038	H399953S1	09/20/12	5			
Rigshospitalet University Hospital of Copenhagen Blegdamsvej 9 2100 Copenhagen East Denmark	RBC-8-038	H399931S1	09/19/12	5	13	In process	40
	RBC-10-038	H399938S1	09/19/12	5			
	RBC-12-038	H399940S1	09/19/12	5			
	RBC-14-038	H399942S1	09/19/12	5			
	RBDC-8-038	H399948S1	09/19/12	5			
	RBDC-10-038	H399950S1	09/19/12	5			
	RBDC-12-038	H399952S1	09/19/12	5			
	RBDC-14-038	H399953S1	09/19/12	5			
Jonas Pedersen Skolelodsvej 1, Gundsoemagle Roskilde, DK-4000 Denmark	RBDC-10-038	H399950S1	23/10/2012	15	0	15	15

Confidential Information

Attachment B

Final Field Safety Notice (FSN) proposal

January xx, 2013

URGENT FIELD SAFETY NOTICE

Address and contact details

Product: Merit ReSolve[®] Locking Biliary Drainage Catheter

Product Numbers: See attached Product Retrieval Form

Lot Number: See attached Product Retrieval Form

Dear contact name,

This letter is to inform you that Merit Medical Systems, Inc. is recalling all ReSolve Locking Biliary Drainage Catheters (RBC & RBDC). Merit has become aware of one catheter failure during a routine 90-day catheter replacement. The catheter tip was observed to be in several pieces under fluoroscopic visualization prior to removal. Unused product must be quarantined immediately and returned to Merit. Do not distribute these devices further. Merit recommends that medical judgment be used in order to assess the risk-benefit of removing catheters from patients. Merit has received no reports of patient harm or injury resulting from this issue.

Our records indicate that you have received ReSolve Locking Biliary Drainage Catheters (RBC & RBDC) as noted on the attached Product Retrieval Form.

Actions required of you:

1. Please immediately quarantine them and do not distribute.
2. Review and complete the attached Product Retrieval Form. Make one copy of the completed form for returning product.
3. Fax a completed, signed copy of the form to Jackie Preece at 801-208-3378.
4. Mail the original completed form using the postage paid, self addressed envelope to:
Jackie Preece
Merit Medical Systems, Inc.
1600 West Merit Parkway, South Jordan, UT 84095-2416
5. Place the copy of the signed completed Product Retrieval Form with the product to be returned to Merit. The form must accompany all products being returned.
6. Return affected products by shipping them back to Merit via Federal Express 2nd Day Air (Acct# 1129-5275-6), include the assigned RGA number from the table above on the outside of the box and ship to:
RGA Department, RGA # (as applicable from Retrieval Form)
Merit Medical Systems, Inc.
Attn: Richard Collard
1600 West Merit Parkway, South Jordan, UT 84095-2416
7. If you have any questions concerning this communication, please don't hesitate to contact Jesse Hansen at 801-826-4034.

Your prompt attention and cooperation is greatly appreciated.

Yours sincerely,

Merit Medical Ireland Ltd.
Mark Mullaney,
Manager, Regulatory Affairs
Phone #: +353 91 703 761; Fax #: +353 91 771 888
E-mail: mmullaney@merit.com

Merit Medical Systems, Inc.
Casey Hughes
Manager, Regulatory Affairs
Phone #: 801-316-4932 Fax #: 801-316-4860

Urgent Field Safety Notice Product Retrieval Form

Merit Medical Systems, Inc.

Address	Site Representative _____
	Title _____
	Phone # _____
	Date _____

THIS IS AN URGENT PRODUCT RECALL / FIELD SAFETY NOTICE:

Reason for product retrieval:

Merit Medical Systems, Inc. is recalling all ReSolve® Locking Biliary Drainage Catheters (RBC & RBDC). Merit has become aware of one catheter failure during a routine 90-day catheter replacement. The catheter tip was observed to be in several pieces under fluoroscopic visualization prior to removal. Unused product must be quarantined immediately and returned to Merit. Do not distribute these devices further. Merit recommends that medical judgment be used in order to assess the risk-benefit of removing catheters from patients. Merit has received no reports of patient harm or injury resulting from this issue.

Please provide status on the following:

Lot #	Part #	Qty	Mfg. Date	Exp. Date	Customer PO #	Merit Order #	RGA	Qty Used	Qty On Hand & Unused	Qty Being Returned	Date Returned	Shipment Tracking #

Qty Discarded	Issue Credit	Send Replacements	Qty	** Sales Order #	** Ship Via

Attention: _____

Department: _____

** Merit Medical Customer Service Use Only

Please fill out and sign the Product Retrieval Form and complete the following steps. It is very important that you complete these steps in order to assist Merit in complying with applicable government regulations.

1. Fax the completed Product Retrieval Form to: Jackie Preece 801-208-3378
2. Mail the completed Product Retrieval Form from (#1) above using the provided postage-paid, self-addressed envelope to Merit.
3. Place a copy of the signed completed Product Retrieval Form with the products to be returned as below. The form must accompany all products being returned to Merit.

Product Return Instructions

Return the affected products by shipping them back to Merit via Federal Express 2nd Day Air (Acct# 1129-5275-6), include the assigned RGA number from the table above on the outside of the box and ship to:

RGA Department, RGA # (as applicable from table above)
 Merit Medical Systems, Inc.
 Attn: Richard Collard
 1600 West Merit Parkway, South Jordan, UT 84095-2416

If you have further questions, please contact: Jesse Hansen 801-826-4034

I certify that I received and understood this notice. I certify that the above listed products have been either: used, removed from inventory, discarded, or returned to Merit Medical Systems, Inc. according to the notification instructions. Furthermore, if I have further distributed product listed on this form, I certify that a copy of this notice was provided to said consignee(s).

 Signature of Site Representative

 Date

Confidential Information

Attachment C

Product Labelling

ReSolve®

HYDROPHILIC COATED


Locking Biliary Drainage Catheter

10F (Standard)

(3.3 mm)

Drainage Location

12 cm

REF RBDC-10-038

MAX GUIDE WIRE
0.038"
(0.97 mm)

LOT RBDC-10-038

Use By: 1980-01

Caution: Consult accompanying documents

U.S. Patent Nos. 7,740,608 and 7,909,814 and Others Pending

ReSolve®

HYDROPHILIC COATED


Locking Biliary Drainage Catheter

10F (Standard)

(3.3 mm)

Drainage Location

12 cm

REF RBDC-10-038

MAX GUIDE WIRE
0.038"
(0.97 mm)

LOT RBDC-10-038

Use By: 1980-01

Caution: Consult accompanying documents

U.S. Patent Nos. 7,740,608 and 7,909,814 and Others Pending

ReSolve®

HYDROPHILIC COATED

Locking Biliary Drainage Catheter

- 1 ReSolve Biliary Drainage Catheter with locking pigtail and straightener
- 1 Dead end cap
- 1 Flexible Stiffening Cannula
- 1 Metal Stiffening Cannula
- 1 Repositioning tool

Drainage Location

12 cr

REF RBDC-10-038

LOT RBDC-10-038

Use By: 1980-01

MAX GUIDE WIRE
0.038" (0.97 mm)

STERILE EO Single Use Caution: accomp docume

Rx ONLY U.S. Patent Nos. 7,740,608 and 7,909,814 and Others Pending. Sterile if package is unopened or undamaged. Non-pyrogenic. Merit Medical Systems, Inc. 1600 West Merit Parkway Sol 1-801-253-1600 U.S.A. Customer Service 1-800-356-3748 EC REP Merit Medical Ireland Ltd, Parkmore Business Park West, European Customer Service +31 43 358 82 22 www.merit.com

ReSolve®

HYDROPHILIC COATED

REF RBDC-10-038

LOT

10F ReSolve® 40 cm /0.038" (0.97 mm) Locking Biliary Drainage Catheter

Description: 10F ReSolve® 40 cm /0.038" (0.97 mm) Locking Biliary Drainage Catheter RBDC-10-038

Description: 10F ReSolve® 40 cm /0.038" (0.97 mm) Locking Biliary Drainage Catheter RBDC-10-038

Non-pyrogenic Sterile if package is unopened or undamaged.

STERILE EO

CE 0086

Single Use

1 Unit

10/08/32 AM 01/17/13 LFCAZ0094 REV 004 08/12

Merit Medical Systems, Inc. 1600 W. Merit Parkway South Jordan, Utah 84095 U.S.A. Customer Service 1-800-356-3748 EC REP Merit Medical Ireland Ltd Parkmore Business Park West, Galway, Ireland European Customer Service +31 43 358 82 22 www.merit.com Made in U.S.A. Rx ONLY ID 033012 402304009/A

Description: 10F ReSolve® 40 cm /0.038" (0.97 mm) Locking Biliary Drainage Catheter RBDC-10-038

Description: 10F ReSolve® 40 cm /0.038" (0.97 mm) Locking Biliary Drainage Catheter RBDC-10-038

Non-pyrogenic Sterile if package is unopened or undamaged.

STERILE EO

CE 0086

Single Use

1 Unit

10/08/32 AM 01/17/13 LFCAZ0094 REV 004 08/12

Merit Medical Systems, Inc. 1600 W. Merit Parkway South Jordan, Utah 84095 U.S.A. Customer Service 1-800-356-3748 EC REP Merit Medical Ireland Ltd Parkmore Business Park West, Galway, Ireland European Customer Service +31 43 358 82 22 www.merit.com Made in U.S.A. Rx ONLY ID 033012 402304009/A

Description: 10F ReSolve® 40 cm /0.038" (0.97 mm) Locking Biliary Drainage Catheter RBDC-10-038

Description: 10F ReSolve® 40 cm /0.038" (0.97 mm) Locking Biliary Drainage Catheter RBDC-10-038

ReSolve®

HYDROPHILIC COATED


Locking Biliary Drainage Catheter

10F (Long)

(3.3 mm)

Drainage Location

15 cm

REF RBC-10-038

MAX GUIDE WIRE
0.038"
(0.97 mm)

LOT RBC-10-038

Use By: 1980-01

Caution: Consult accompanying documents

U.S. Patent Nos. 7,740,608 and 7,909,814 and Others Pending.

Description:
10F ReSolve® 40 cm /0.038" (0.97 mm)
Locking Biliary Drainage Catheter
RBC-10-038

Description:
10F ReSolve® 40 cm /0.038" (0.97 mm)
Locking Biliary Drainage Catheter
RBC-10-038

Non-pyrogenic
Sterile if package is unopened or undamaged.

STERILE EO

CE 0086

Single Use

1 Unit

10:08:45 AM 01/17/13 LFCAZ0094 REV 004 08/12

Merit Medical Systems, Inc.
1600 W. Merit Parkway South Jordan, Utah 84095 U.S.A.
Customer Service 1-800-356-3748
EC REP Merit Medical Ireland Ltd Parkmore Business Park
West, Galway, Ireland European Customer Service
+31 43 358 82 22 www.merit.com
Made in U.S.A.
R_X ONLY ID 033012 402304009/A

Description:
10F ReSolve® 40 cm /0.038" (0.97 mm)
Locking Biliary Drainage Catheter
RBC-10-038

ReSolve®

HYDROPHILIC COATED


Locking Biliary Drainage Catheter

10F (Long)

(3.3 mm)

Drainage Location

15 cm

REF RBC-10-038

MAX GUIDE WIRE
0.038"
(0.97 mm)

LOT RBC-10-038

Use By: 1980-01

Caution: Consult accompanying documents

U.S. Patent Nos. 7,740,608 and 7,909,814 and Others Pending.

Description:
10F ReSolve® 40 cm /0.038" (0.97 mm)
Locking Biliary Drainage Catheter
RBC-10-038

Description:
10F ReSolve® 40 cm /0.038" (0.97 mm)
Locking Biliary Drainage Catheter
RBC-10-038

Non-pyrogenic
Sterile if package is unopened or undamaged.

STERILE EO

CE 0086

Single Use

1 Unit

10:08:45 AM 01/17/13 LFCAZ0094 REV 004 08/12

Merit Medical Systems, Inc.
1600 W. Merit Parkway South Jordan, Utah 84095 U.S.A.
Customer Service 1-800-356-3748
EC REP Merit Medical Ireland Ltd Parkmore Business Park
West, Galway, Ireland European Customer Service
+31 43 358 82 22 www.merit.com
Made in U.S.A.
R_X ONLY ID 033012 402304009/A

Description:
10F ReSolve® 40 cm /0.038" (0.97 mm)
Locking Biliary Drainage Catheter
RBC-10-038

ReSolve®

HYDROPHILIC COATED

Locking Biliary Drainage Catheter

- 1 ReSolve Biliary Drainage Catheter with locking pigtail and straightener
- 1 Dead end cap
- 1 Flexible Stiffening Cannula
- 1 Metal Stiffening Cannula
- 1 Repositioning tool

Drainage Location

15 cm

REF RBC-10-038

LOT RBC-10-038

Use By: 1980-01

MAX GUIDE WIRE
0.038" (0.97 mm)

(01)10884450012809(17)800100(10)RBC-10-038

STERILE EO

Single Use

Caution

R_X ONLY

U.S. Patent Nos. 7,740,608 and 7,909,814 and Others Pending. Sterile if package is unopened or undamaged. Non-pyrogenic

Merit Medical Systems, Inc. 1600 West Merit Parkway
1-801-253-1600 U.S.A. Customer Service 1-800-356-3748

EC REP Merit Medical Ireland Ltd, Parkmore Business Park West
European Customer Service +31 43 358 82 22
www.merit.com

ReSolve®

HYDROPHILIC COATED

REF RBC-10-038

LOT

10F ReSolve® 40 cm /0.038" (0.97 mm)
Locking Biliary Drainage Catheter

0002 01/1

ReSolve[®]

Biliary Catheter

INSTRUCTIONS FOR USE
MODE D'EMPLOI
ISTRUZIONI PER L'USO
GEBRUIKSAANWIJZING
INSTRUCCIONES DE USO
INSTRUÇÕES DE UTILIZAÇÃO
GEBRUIKSAANWIJZING
BRUKSANVISNING
BRUGERVEJLEDNING
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ

ReSolve[®] Biliary Catheter

English


INSTRUCTIONS FOR USE

PRODUCT DESCRIPTION

The ReSolve Biliary Catheter with locking pigtail and hydrophilic coating is a multi-side holed, radiopaque catheter used for percutaneous biliary drainage. The components of the catheter allow for introduction and placement using the over-the-wire method.

The product may be packaged with the following components:

- One (1) ReSolve Biliary Drainage Catheter with locking pigtail (A), hydrophilic coating and pigtail straightener (B) drainage location: Standard or Long (G)
- One (1) Metal stiffening cannula (C)
- One (1) Flexible stiffening cannula (D)
- One (1) Repositioning Tool (E)
- One (1) Dead end cap (F)

INTENDED USE OF PRODUCT

INDICATIONS:

The ReSolve Biliary Drainage Catheter with locking pigtail and hydrophilic coating is used for drainage of bile within the biliary system.

CONTRAINDICATIONS:

The ReSolve Biliary Drainage Catheter is contraindicated for use where percutaneous biliary drainage catheterization is unacceptable.

The ReSolve Biliary Drainage Catheter is contraindicated for intravascular use.

CAUTIONS: Where long-term use is indicated, it is recommended that indwelling time not exceed 90 days. This catheter should be evaluated by the physician on or before 90 days post-placement.

To maximize the advantages of the hydrophilic coating on the surface of the distal portion of the catheter, wet the catheter prior to use with sterile water or saline. Keep catheter wet during placement.

WARNINGS:

- RX Only Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.
- Contents of unopened, undamaged package are sterile and non-pyrogenic
- This device is intended for single use only
- DO NOT wipe catheter with dry gauze or any solvents because it may damage the catheter coating
- DO NOT over tighten the connection between the drainage catheter and drainage tubing or dead end cap
- If using alcohol to clean the catheter hub, allow sufficient time for alcohol to dry before connecting the drainage tubing or dead end cap

REUSE PRECAUTION STATEMENT

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or crossinfection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.


POTENTIAL COMPLICATIONS:

- Hemobilia
- Pneumothorax
- Sepsis
- Catheter leak
- Hemorrhage
- Perforation of the bile ducts, liver and/or duodenum
- Cholangitis
- Hemothorax
- Biloma
- Catheter occlusion
- Bile peritonitis
- Pancreatitis
- Bilithorax
- Fever
- Catheter dislodgement
- Skin Infection

INSTRUCTIONS FOR USE:

Seldinger Entry Technique or Guide Wire Exchange

1. Remove the stiffening cannula from catheter.
2. Ensure that the distal portion of the catheter is wet prior to placement. See Cautions, Page 1.
3. Flush catheter prior to use.
4. Ensure suture locking mechanism is in the proximal position.
5. Slide pigtail straightener along distal portion of catheter to straighten curve prior to placing stiffening cannula into the catheter. Place the stiffening cannula into the catheter and tighten the luer lock fittings. See Figure 1.

6. Remove pigtail straightener from catheter prior to insertion.
7. Place catheter/cannula assembly over appropriate guide wire and advance into the biliary system. Disconnect stiffening cannula from catheter and continue to advance catheter over guide wire into desired position. The catheter accommodates a 0.038" (0.97mm) wire. See Figure 2. Placement should be confirmed with diagnostic imaging.


CAUTION: Failure to follow could cause stiffening cannula to lodge in catheter.

8. After placement is confirmed, remove cannula and guide wire.


Engaging the Suture Locking Mechanism

9. Pull the suture until desired pigtail is formed. Do not overly tighten the pigtail as over tightening may damage the catheter. Rotate the suture locking mechanism distally to hold the suture in place. See Figure 3.

NOTE: If the catheter needs to be repositioned, unlock the suture locking mechanism by rotating the arm proximally to the point of resistance. Do not rotate the suture locking mechanism beyond the point of resistance.

10. Once placement is confirmed, and the suture locking mechanism has been rotated to the most distal position, press the suture locking mechanism into the hub to secure it. The suture locking mechanism is now locked into position. See Figure 4.

11. The catheter is now ready to be connected to appropriate drainage bag or tubing.
12. A flush regimen should be designed for the circumstances of each patient and the protocol of the physician.

CATHETER EXCHANGE, REPOSITIONING OR REMOVAL

1. Disconnect catheter from drainage bag or tubing.
2. To release the pigtail loop choose one of the following options:


Option 1: Using the Repositioning Tool, align the opening of the round section of the Repositioning Tool in line with the handle of the suture locking mechanism which will release the suture and allow the pigtail to straighten upon removal.

Option 2: For exchange or removal only, advance 0.038" (0.97mm) heavy duty floppy tip guide wire past distal tip of catheter cut the hub off the drainage catheter and sever suture. This will release the suture and the pigtail loop.

CAUTION: The suture will no longer be secured to the catheter. Take care to remove both the suture and catheter.

3. For catheter exchange or if access is to be maintained, advance appropriate guide wire through catheter; use diagnostic imaging to confirm wire placement. Guide wire will maintain access to drainage site. To ease guide wire placement, the flexible stiffening cannula may be used.

4. Carefully remove the catheter. Proceed with either catheter exchange or skin closure.

ATTENTION ATTENDING PHYSICIAN: IF PATIENT WILL NOT BE FOLLOWED UP BY YOU, IT IS RECOMMENDED THAT THE "INSTRUCTIONS FOR USE" OR THE SECTION ON HOW TO REMOVE THE CATHETER BE ATTACHED TO THE PATIENT'S CHART.

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Manufacturer: Merit Medical Systems, Inc. South Jordan, Utah 84095 U.S.A. 1-801-253-1600
Authorized Representative: Merit Medical Ireland Ltd, Parkmore Business Park West, Galway, Ireland