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ATT.: Monica R. Roos

SSI Diagnostica hereby recalls Dermatophyt PCR kit, lot no. DERM S

The reason for the recall is that the internal plasmid control in the primer mix in some PCR test is not amplified. Such tests are inconclusive.

PCR test of the plasmid stock, used for internal control in the primer mix in DERM-S shows that the plasmid was unstable when stored at -20 degree Celcius. The instability was not observed in our QC check.

As the dermatophyte PCR is a multiplex PCR there is competition between the three possible fragments amplified (internal control, pan-derm and *T. rubrum* fragment). A decrease in the amounts of amplified internal control fragment will lead to a higher specificity for the amplification of *T. rubrum* and pan-derm fragments, and vice versa.

The combination between the unstable plasmid stock and the aim of getting the highest specificity resulted in production of a lot where the plasmid concentration was below the limit, resulting in fine amplification of the internal plasmid in some cases but absolutely no internal control in other cases.

The positive results obtained with this lot are valid, but negative samples may be inconclusive due to the lack the internal control fragment.

SSI Diagnostica would like you to dispose of the faulty Dermatophyte PCR lot DERM S.

You will receive a replacement lot or a credit note as agreed upon in the earlier phone or mail contact.

This notice needs to be passed on to all those who need to be aware within your organization where the potentially affected devices have been transferred.

Please transfer this notice to other organizations on which this action has an impact.

Yours sincerely,

Jytte Jørgensen

Jytte Jørgensen, QA/RA manager

Tel: +45 48299124

jjj@ssi.dk

MEDIQ EESTI OÜ
Kungla Street 2
76505 Saue
Estonia

SSI Diagnostica
Our ref:

13 October 2014

Urgent Field Safety Notice

Dermatophyt PCR-kit
FSCA 13-10-2014 (01-2014)
Recall of Dermatophyt PCR-kit lot DERM-S

ATT.: To whom it may concern

SSI Diagnostica hereby recalls Dermatophyt PCR kit, lot no. DERM S
The reason for the recall is that the internal plasmide control in the primer mix in some PCR test is not amplified. Such tests are inconclusive.

PCR test of the plasmid stock, used for internal control in the primer mix in DERM-S shows that the plasmid was unstable when stored at -20 degree Celcius. The instability was not observed in our QC check.

As the dermatophyte PCR is a multiplex PCR there is competition between the three possible fragments amplified (internal control, pan-derm and *T. rubrum* fragment). A decrease in the amounts of amplified internal control fragment will lead to a higher specificity for the amplification of *T. rubrum* and pan-derm fragments, and vice versa.

The combination between the unstable plasmid stock and the aim of getting the highest specificity resulted in production of a lot where the plasmid concentration was below the limit, resulting in fine amplification of the internal plasmid in some cases but absolutely no internal control in other cases.

The positive results obtained with this lot are valid, but negative samples may be inconclusive due to the lack of the internal control fragment.

SSI Diagnostica would like you to dispose of the faulty Dermatophyte PCR lot DERM S.
You will receive a replacement lot or a credit note as agreed upon in the earlier phone or mail contact.

This notice needs to be passed on to all those who need to be aware within your organization where the potentially affected devices have been transferred.

Please transfer this notice to other organizations on which this action has an impact.

Yours sincerely,

Jytte Jørgensen

Jytte Jørgensen, QARA manager

Tel: +45 48299124

jjj@ssi.dk