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**Vedlagt:**

Bilag 1: Mycobacterium Data Sheet

Bilag 2: Disinfection Excerpt from the Operating Instructions for 3T Heater Cooler

Bilag 3: Customer Response Form

## **Bilag 1**

### **Mycobacterium Data Sheet:**

#### **IMPORTANT INFORMATION**

Sorin Heater Cooler Devices - Mycobacterium  
Reference # IIS 9611109-07-14-14

#### **FOR GENERAL HEALTHCARE STAFF - INTRODUCTION**

- Mycobacteria are widely distributed in the ecosystem, are present in water, and even survive in chlorinated drinking water. Some species are classic human pathogens, for example *M. tuberculosis* that causes TB. However, most mycobacteria are not considered harmful to humans. These environmental mycobacteria are termed nontuberculous mycobacteria (NTM). These do occasionally cause opportunistic infections. (Van Ingen, J Med Microbiol, September 2012 vol. 61 no. Pt 9, 1234-1239).
- For example, *M. chimaera* was identified as the causative agent of a respiratory tract infection in patients with cystic fibrosis. (Cohen-Bacrie et al. Journal of Medical Case Reports 2011, 5:473).
- The various NTM species have greatly differing microbiological culture requirements and infective characteristics. Special methods must be employed to determine presence or absence and to obtain accurate identification.
- In order to control these organisms in medical environments and equipment, it is critical to identify effective chemicals against NTM and follow Instructions for Use provided by the manufacturer of disinfectants and medical equipment/devices.
- Particular attention **must** be paid to the use of water in critical health care situations because water that is not sterilized, for example by filtration, often carries these kinds of bacteria.

#### **MYCOBACTERIA ARE INHERENTLY RESISTANT TO CHEMICAL DISINFECTANTS AND ANTIBIOTICS**

- Mycobacteria have a natural resistance to many disinfectants. All mycobacteria share a characteristic cell wall, which is thicker than in most bacteria, which is hydrophobic, waxy, and rich in mycolic acids. This cell wall makes a substantial contribution to the environmental and chemical tolerance of this group.  
[http://www.cdc.gov/hicpac/Disinfection\\_Sterilization/4\\_0efficacyDS.html](http://www.cdc.gov/hicpac/Disinfection_Sterilization/4_0efficacyDS.html)

#### **FOR TECHNICAL & MEDICAL STAFF - MYCOBACTERIA REQUIRE SPECIALIZED MICROBIOLOGY**

- Mycobacteria require specialized growth media and techniques for successful culture. Therefore, they may be present in the samples and not detected due to the inadequate culture media or incubation periods. Extension of incubation for up to 60 days or more may be required to recover these types of isolates.
- Sampling for these organisms can be made challenging by their highly hydrophobic nature. For sampling surfaces directly, or via a liquid medium, swabbing vigorously to recover these organisms may be needed. Also, using a surfactant (e.g. Tween) can assist in recovering during sampling and manipulation during culture/testing.
- It is likely that the presence of these NTM can be overlooked by routine medical microbiological and environmental testing.
- Sampling of the environment (air and surfaces), equipment, and clinical specimens, should be considered whenever mycobacterial infection is a risk. This can include water filled equipment in health care settings. For general mycobacteriology see:  
[http://www.health.gov.au/internet/main/publishing.nsf/Content/cda-cdi3001-pdf-cnt.htm/\\$FILE/cdi3001f.pdf](http://www.health.gov.au/internet/main/publishing.nsf/Content/cda-cdi3001-pdf-cnt.htm/$FILE/cdi3001f.pdf)  
For detailed on NTM information see:  
<https://www.thoracic.org/statements/resources/mtpi/nontuberculous-mycobacterial-diseases.pdf>

**Bilag 2**  
**Operating Instructions Excerpt:**  
Cleaning and Disinfection of the Heater Cooler Device

**IMPORTANT INFORMATION**

Sorin Heater Cooler Devices - Mycobacterium  
Reference # IIS 9611109-07-14-14

The following chapters are attached from the current Operating Manual GA-16-XX-XX ENG (rev 11) of the 3T heater cooler device.

5.2 Filling the water tanks

6.2 Cleaning and disinfection of the housing

6.2.1 Disinfection of the water circuits

6.2.2 Protecting the water circuits from microbial growth

## Bilag 3 Customer Response Form

### IMPORTANT INFORMATION

Sorin Heater Cooler Devices - Mycobacterium

Reference # IIS 9611109-07-14-14

According to our records you have the following impacted products:

<< Fill in the customer related codes and serial numbers only- Use Product Trace list (Excel File)>>

Product Code	Product description	Serial Number

**Please return this completed form to:**

Sorin Site/ Distributor Name: <<Print Your Company name here>>  
 Country: <<Print Your Country here>>  
 Contact Name: << Print Your Contact Name here>>  
 E-mail: <<Print Your E-mail address here>>  
 Fax No.: <<Print Your Fax No. here>>  
 Phone Number: <<Print Your Phone No. here>>

Section 1 - Please Complete:

1. We HAVE reviewed and understand the attached Important Information
2. Yes- We do have the listed affected products and we will follow the instruction
3. We DO NOT understand the attached Important Information and request more information

Please contact us:

Christian Peis, Director QA, Sorin Group Deutschland GmbH at +49 89 323 01 152, via fax at +49 89 323 01 333 or via e-mail at [SGD.fsca@sorin.com](mailto:SGD.fsca@sorin.com)

Customer Name: <<Print Your Company name here>>  
 Country: <<Print Your Country here>>  
 Contact Name: << Print Your Contact Name here>>  
 E-mail: <<Print Your E-mail address here>>  
 Fax No.: <<Print Your Fax No. here>>  
 Phone Number: <<Print Your Phone No. here>>

Submitted by .....

Signature .....

Date ...../...../.....