

**To the attention of:
Distributor of products from Synthes GmbH**

9 December 2015

**URGENT NOTICE:
MEDICAL DEVICE FIELD SAFETY NOTIFICATION – FSN2015163
Kirschner Wire 1.6mm w/Drill Tip, L 200mm (sterile)**

Please note that this is a Medical Device Field Safety Notification only; it is not required to return any of the affected products.

Part Description, Part and Lot Numbers

Part Description	Part Numbers	Lot Numbers
Kirschner Wire Ø 1.6 mm, with drill tip, length 200 mm	02.113.001	All
Kirschner Wire Ø 1.6 mm, with drill tip, length 200 mm, sterile	02.113.001S	All

Dear Sir/Madam,

Synthes GmbH is initiating a Field Safety Notification for the Kirschner Wire for the above mentioned Part and Lot Numbers. Kirschner Wire implants are indicated for a wide range of orthopedic trauma applications including stand-alone device for fracture fixation and fracture fixation in conjunction with other fixation systems.

Our records indicate that your facility may have the affected product(s) subject to this Field Safety Notification or has been using the affected product(s) from a loan set.

Synthes asks that you review the information contained in this Field Safety Notification and complete the actions on page two of this notification.

Reason for the Field Safety Notification

The Kirschner wire is made from Co35Ni35Cr20Mo10 (MP35N) which is a cobalt based alloy containing 35% nickel. However, the label and technique guide incorrectly indicate that the devices are made of stainless steel.

The product is shipped with an insert which references the potential risk of allergy/hypersensitivity reactions, but there is no specific information on the label or in the technique guide with regards to the percentage of nickel content in the devices.

Potential hazard

For patients with known nickel sensitivity, products labeled as containing stainless steel should be avoided since stainless steel contains nickel and an alternate product is generally chosen.

Implanted wires or unintentionally retained fragments that are left in the patient pose a potential risk of an Adverse Tissue Reaction which could be worse in a population sensitive to nickel. As MP35N is a more corrosion-resistant alloy and more biocompatible than stainless steel, the likelihood of a patient being exposed and sensitized to a higher amount of nickel or other alloys is low compared to a product manufactured from stainless steel. Therefore, the risk for an Adverse Tissue Reaction is lower than if the product was made of stainless steel. MP35N is considered implant-grade quality.

Customer immediate actions

Please take the following actions:

- Review the instructions listed in this notification.
- Forward this notice to anyone in your facility that needs to be informed.
- If any of the affected products has been forwarded to another facility, contact that facility and provide them with a copy of this letter.
- Identify all customers who have bought the affected products, provide them with a copy of the Accounts Letter and ask these customers to complete and return the Verification Section (page 3 of the Account Letter).
- Report the following information to Synthes Field Action in Zuchwil (fieldaction.int@synthes.com): total number of hospitals/accounts affected, total number of signed acknowledgements received, for not responding customers the total number of reminders sent to those clients who did not sign the Verification Section.
- Keep a copy of this notice.

We sincerely apologize for any inconvenience that this issue may create. Should you have any questions please do not hesitate to contact your DePuy Synthes sales consultant or contact person.

Thank you for your attention and cooperation.

Synthes GmbH

Pierre van Iwaarden
Field Action Manager

David Carvin
Sr. QA Manager, Product Safety and
Performance

Cc:

Account Name: _____

**URGENT NOTICE:
MEDICAL DEVICE FIELD SAFETY NOTIFICATION – FSN2015163
Kirschner Wire 1.6mm w/Drill Tip, L 200mm (sterile)**

Please note that this is a Medical Device Field Safety Notification only; it is not required to return any of the affected products.

Part Description, Part and Lot Numbers

Part Description	Part Numbers	Lot Numbers
Kirschner Wire Ø 1.6 mm, with drill tip, length 200 mm	02.113.001	All
Kirschner Wire Ø 1.6 mm, with drill tip, length 200 mm, sterile	02.113.001S	All

____ We have located the Kirschner Wire Ø 1.6 mm, with drill tip, length 200 mm (sterile) within our stock, and acknowledge receipt of the updated information.

____ We acknowledge receipt of this information, but do not have Kirschner Wire Ø 1.6 mm, with drill tip, length 200 mm (sterile) in stock.

Name/Title (please print): _____

Address: _____

Country: _____

Phone Number: _____

Signature and Date: _____

Please complete and return this page your local DePuy Synthes contact person

Note: If the Verification Section is answered on behalf of more than one facility and/or individual, please clearly indicate the name and address of the facility and/or individual on this page of the notification.