Viktig produktinformation

Produkt: Seronorm™ Paediatric
Varenr.: 208005
FSN ID: 2016-03
Handling: Information

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Att.: SERO AS-grossist /SERO AS-kunde

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
Med venlig hilsen

Silje Solheim Johnsen
Head of Product Management

Åse Davanger
Quality Manager

Seronorm™ Paediatric

REF 208005

LOT 1412555

2019-01

6 x 3 mL

IVD

CE

Collaborating laboratories

The analytical data of Seronorm™ Paediatric have been elaborated in collaboration with the following independent laboratories:

Collaborating laboratories:

- The Central Laboratory of Vestfold Hospital, Tønsberg, Norway
- Laboratory of Akershus University Hospital, Lørenskog, Norway
- The Laboratory of SERO AS, Billingstad, Norway

EN

Storage and stability

This product will be stable until the expiry date when stored unopened at 2-8 °C.

After reconstitution, all analytes are stable for 7 days when stored tightly capped at 2-8 °C.

Vials aliquoted and frozen within 30 minutes at ≤ -20 °C in screw cap tubes with O-ring are stable for 1 month.

Limitations

All stability data require that bacterial contamination is avoided. Increased turbidity may indicate bacterial growth.

Assignment of values

Seronorm™ Paediatric is intended for use as an assayed quality control material/serum to monitor the precision and trueness of laboratory measurement procedures. Analytical values of components which can be used to monitor trueness are given with an uncertainty value (U) and traceability in

accordance with ISO 17511¹⁾ and ISO 18153²⁾. All the analyses have been performed by independent collaborative laboratories using SERO's protocol. The reference values are based on reference methods or transferred values in accordance to certified reference materials. At the time of establishment of this package insert, this reference method or material constitutes the top of the calibration hierarchy for the given components. These values are presented with shaded background.

The mean analytical values and acceptable ranges given in the table below are derived from replicate analyses obtained through collaboration with a number of independent laboratories and they are specific for this lot. The uncertainty of the analytical value is presented as a single number, U (with a coverage factor $k = 2$), which takes into account various factors including imprecision of the method applied for the assignment as well as the stability of the given component. Individual laboratory means should fall within the acceptable range. However, variations in the reagent and/or calibrator lots as well as the instrument and/or method used. It is therefore recommended that each laboratory establishes its own analytical mean and acceptable range and only use the provided values as a guide. For information, the relevant Rilibak range (Richtlinien der Bundesärztekammer zur Qualitätssicherung in medizinischen Laboratorien (2008)) is provided.

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- 1) ISO 17511. In Vitro Diagnostic Medical Devices- Measurement of quantities in samples of biological origin- Metrological traceability of values assigned to calibrators and control materials. Geneva, International Organisation for Standardization 2003.
 2) ISO 18153 In vitro diagnostic medical devices -- Measurement of quantities in biological samples -- Metrological traceability of values for catalytic concentration of enzymes assigned calibrators and control materials.
 3) This is not a Rilibäk range, but SERO's assessment