Viktig produktinformation

Produkt: Pathonorm™ H
Varenr.: 100705
FSN ID: 2016-02
Handling: Information

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Att.: SERO AS-grossist /SERO AS-kunde

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Med venlig hilsen

Silje Solheim Johnsen
Head of Product Management

Åse Davanger
Quality Manager

Pathonorm™ H

REF 100705

LOT 1403120

2018-04

10 x 5 mL

IVD

CE

Collaborating laboratories

The analytical data of Pathonorm™ H have been elaborated in collaboration with the following independent laboratories:

- Department of Clinical Chemistry, Asker og Bærum Hospital, Bærum, Norway
- Furst Medical Laboratory, Oslo, Norway

- The Laboratory of SERO AS, Billingstad, Norway
- Department of Medical Biochemistry, St.Olav's Hospital, Trondheim, Norway

EN

Storage and stability

This product will be stable until the expiry date when stored unopened at 2-8 °C.

After reconstitution, all analytes are stable for 7 days when stored tightly capped at 2-8 °C.

Vials to be aliquoted and frozen within 30 minutes are stable 1 month at ≤ -20 °C.

Limitations

All stability data require that bacterial contamination is avoided. Increased turbidity may indicate bacterial growth.

CK is sensitive to temperature changes.
CK and Bilirubin are sensitive to light.

Assignment of values

The mean analytical values and acceptable ranges given in the table below are derived from replicate analyses obtained through collaboration with a number of independent laboratories and they are specific for this lot. The uncertainty of the analytical value is presented as a single number, U (with a coverage factor $k = 2$), which takes into account various factors including imprecision of the method applied for the assignment as well as the stability of the given component. Individual laboratory means should fall within the acceptable range. However, variations in the analytical value and acceptable range for each component may vary over time due to differences in reagent and/or calibrator lots as well as the instrument and/or method used. It is therefore recommended that each laboratory establishes its own analytical mean and acceptable range and only use the provided values as a guide. If applicable, the Rilibäk range (Richtlinien der Bundesärztekammer zur Qualitätssicherung in medizinischen Laboratorien (2008)) is provided.

The value assignment has been established in accordance with the Essential Requirements (Annex 1 of the IVD Directive) 98/79/EC, and the ISO 17511² International standard.

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1) Council Directive 98/79/EC of the European Parliament and the Council of 27 October 1998 on In Vitro Diagnostic Medical Devices (Annex 1: Essential requirements).

2) ISO 17511. In Vitro Diagnostic Medical Devices- Measurement of quantities in samples of biological origin- Metrological traceability of values assigned to calibrators and control materials. Geneva, International Organisation for Standardization 2003.

3) This is not a Rilibäk range, but SERO's assessment

4) To avoid Zinc contamination from vial and rubber stopper, transfer the serum to plastic tubes after reconstitution