### Viktig produktinformation

Produkt: Seronorm™  
Varenr.: 100105  
FSN ID: 2016-01  
Handling: Information

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**Att.:** SERO AS-grossist /SERO AS-kunde

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Med venlig hilsen

Silje Solheim Johnsen  
Head of Product Management

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- Department of Clinical Chemistry, Asker og Bærum Hospital, Bærum, Norway
- Furst Medical Laboratory, Oslo, Norway
- The Laboratory of SERO AS, Billingstad, Norway

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#### Assignment of values

Seronorm™ is intended for use as an assayed quality control material/serum to monitor the precision and trueness of laboratory measurement procedures. Analytical values of components which can be used to monitor trueness are given with an uncertainty value (U) and traceability in accordance with ISO 17511<sup>1)</sup> and ISO 18153<sup>2)</sup>. All the analyses have been performed by independent collaborative laboratories using SERO's protocol.

The mean analytical values and acceptable ranges given in the table below are derived from replicate analyses obtained through collaboration with a number of independent laboratories and they are specific for this lot. The uncertainty of the analytical value is presented as a single number, U (with a coverage factor  $k = 2$ ), which takes into account various factors including imprecision of the method applied for the assignment as well as the stability of the given component. Individual laboratory means should fall within the acceptable range. However, variations in the analytical value and acceptable range for each component may vary over time due to differences in reagent and/or calibrator lots as well as the instrument and/or method used. It is therefore recommended that each laboratory establishes its own analytical mean and acceptable range and only use the provided values as a guide. For information, the relevant Rilibak range (Richtlinien der Bundesärztekammer zur Qualitätssicherung in medizinischen Laboratorien (2008)) is provided.

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- 1) ISO 17511. In Vitro Diagnostic Medical Devices – Measurement of quantities in samples of biological origin – Metrological traceability of values assigned to calibrators and control materials. Geneva, International Organisation for Standardization 2003.
- 2) ISO 18153 In vitro diagnostic medical devices – Measurement of quantities in biological samples – Metrological traceability of values for catalytic concentration of enzymes assigned calibrators and control materials.
- 3) This is not a Rilibäk range, but SERO's assessment.

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- The Laboratory of SERO AS, Billingstad, Norway
- Department of Medical Biochemistry, St.Olav's Hospital, Trondheim, Norway

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Vials to be aliquoted and frozen within 30 minutes are stable 1 month at  $\leq -20$  °C.

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#### Assignment of values

The mean analytical values and acceptable ranges given in the table below are derived from replicate analyses obtained through collaboration with a number of independent laboratories and they are specific for this lot. The uncertainty of the analytical value is presented as a single number, U (with a coverage factor  $k = 2$ ), which takes into account various factors including imprecision of the method applied for the assignment as well as the stability of the given component. Individual laboratory means should fall within the acceptable range. However, variations in the analytical value and acceptable range for each component may vary over time due to differences in reagent and/or calibrator lots as well as the instrument and/or method used. It is therefore recommended that each laboratory establishes its own analytical mean and acceptable range and only use the provided values as a guide. If applicable, the Rilibak range (Richtlinien der Bundesärztekammer zur Qualitätssicherung in medizinischen Laboratorien (2008)) is provided.

The value assignment has been established in accordance with the Essential Requirements (Annex 1 of the IVD Directive) 98/79/EC, and the ISO 17511<sup>2)</sup> International standard.

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1) Council Directive 98/79/EC of the European Parliament and the Council of 27 October 1998 on In Vitro Diagnostic Medical Devices (Annex 1: Essential requirements).  
 2) ISO 17511. In Vitro Diagnostic Medical Devices- Measurement of quantities in samples of biological origin- Metrological traceability of values assigned to calibrators and control materials. Geneva, International Organisation for Standardization 2003.  
 3) This is not a Rilibäk range, but SERO's assessment  
 4) To avoid Zinc contamination from vial and rubber stopper, transfer the serum to plastic tubes after reconstitution

REF 100105

LOT 1502086

 2019-03

10 x 5 mL

IVD

CE

### Collaborating laboratories

The analytical data of Seronorm™ have been elaborated in collaboration with the following independent laboratories:

- Department of Clinical Chemistry, Asker og Bærum Hospital, Bærum, Norway
- Furst Medical Laboratory, Oslo, Norway
- The Laboratory of SERO AS, Billingstad, Norway
- Department of Medical Biochemistry, St.Olav's Hospital, Trondheim, Norway

### EN

#### Storage and stability

This product will be stable until the expiry date when stored unopened at 2-8 °C.

After reconstitution, all analytes are stable for 7 days when stored tightly capped at 2-8 °C.

Vials aliquoted and frozen within 30 minutes at  $\leq -20$  °C in screw cap tubes with O-ring are stable for 1 month.  
Discard the remaining material after each use.

#### Limitations

All stability data require that bacterial contamination is avoided. Increased turbidity may indicate bacterial growth.

CK is sensitive to temperature changes.  
CK and Bilirubin are sensitive to light.

#### Assignment of values

Seronorm™ is intended for use as an assayed quality control material/serum to monitor the precision and trueness of laboratory measurement procedures. Analytical values of components which can be used to monitor trueness are given with an uncertainty value (U) and traceability in accordance with ISO 17511<sup>®</sup> and ISO 18153<sup>®</sup>. All the analyses have been performed by independent collaborative laboratories using SERO's protocol.

The mean analytical values and acceptable ranges given in the table below are derived from replicate analyses obtained through collaboration with a number of independent laboratories and they are specific for this lot. The uncertainty of the analytical value is presented as a single number, U (with a coverage factor  $k = 2$ ), which takes into account various factors including imprecision of the method applied for the assignment as well as the stability of the given component. Individual laboratory means should fall within the acceptable range. However, variations in the analytical value and acceptable range for each component may vary over time due to differences in reagent and/or calibrator lots as well as the instrument and/or method used. It is therefore recommended that each laboratory establishes its own analytical mean and acceptable range and only use the provided values as a guide. For information, the relevant RilibÄK range (Richtlinien der Bundesärztekammer zur Qualitätssicherung in medizinischen Laboratorien (2008)) is provided.

To avoid Zinc contamination from vial and rubber stopper, transfer the serum to plastic tubes after reconstitution.

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1) ISO 17511. In Vitro Diagnostic Medical Devices – Measurement of quantities in samples of biological origin – Metrological traceability of values assigned to calibrators and control materials. Geneva, International Organisation for Standardization 2003.  
 2) ISO 18153 In vitro diagnostic medical devices – Measurement of quantities in biological samples – Metrological traceability of values for catalytic concentration of enzymes assigned calibrators and control materials.  
 3) This is not a Rilibak range, but SERO's assessment.