



Important Safety Notice

Withdrawal of possible failed assembly of castor fork, model S3.

Date: 16 of June 2016

FSCA-identity: 6.6.2-2016-28983

Type of action: Exchange of castor fork.

Important safety information

Affected component.



Castor fork S3 with Loctite.

Panthera AB has recommended all their customers to exchange this article towards replacement article.

Affected art.no:

2100010,2100030,2100060,2100090
2100020,2100040,2100070,2100100
2100050,2100080,2100110

Description of problem

We have discovered at random sample that a few castor forks were poorly assembled on our new wheelchair Panthera S3. They lacked thread locking, Loctite, between fork and axle. This could in worst case, cause the axle to come loose from the fork and then eventually cause a fall out the wheelchair for the user.

Advice of action that should be taken by the owner/user.

We have contacted all receivers of wheelchairs with the affected component mounted and asked them to exchange it with the attached replacement article.

We urge all our users to change their castor fork towards the offered replacement castor fork.

Forward this safety message

Spread this information to affected persons and to all customers and suppliers.

Panthera AB apologize for the inconvenience this causes and stresses that these actions has been taken with the users safety in mind.

Contact information

Please contact Panthera if you have questions:

Panthera AB

Gunnebogatan 26

163 53 Spåna

Sweden

panthera@panthera.se

www.panthera.se

+46-8-761 40 50

Undersigned verifies that this message has been sent to the supervisory authority.

Nils Sallnäs

Mikael Beskow

VD /CEO

Contact person



.....

.....