30. november 2016

Til: Risikomanagere og kirurger

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#### **Repair Service**

Brasseler USA recommends that all powered devices and accessories be returned for preventative maintenance every six (6) months. Brasseler USA warrants any service or repair work performed will be free from defects in material or workmanship for the period of six (6) months from date of service or repair. This warranty applies to the actual work performed.

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#### **Cleaning Procedure**

1. Leave hose attached or insert the cleaning cap (optional) into the hose connector of the instrument. Remove bur prior to cleaning.
2. Scrub, using a nylon brush, the instrument thoroughly with mild soap and water. Remove all traces of blood and debris.
3. With the air/nitrogen hose or cleaning cap(optional) still attached, rinse all traces of contaminants and detergent under running water.

*Do NOT immerse. If accidental immersion occurs please see Accidental Immersion.*

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
Med venlig hilsen

Kevin W. Escapule  
Post Market Surveillance & Regulatory Compliance Director

# BILAG 1

## Uddrag af kirurgisk teknik

SECTION  
**6**

### Mill Femoral Trochlea


#### 6.2 Perform the Milling Operation

- Open Collet on Milling Handpiece: Unlock the burr locking mechanism by rotating the locking collar into the "open" position (a click can be felt as the collet fully opens).
- Install PFJ Mill Burr (00-5927-050-00): Insert burr fully into collet until flush (Fig. 32).

**TECHNIQUE TIP**

**6.B**


*If burr does not sit flush initially, make sure collet is unlocked and rotate burr while inserting.*

**Fig. 32**

Insert burr into collet.

- Lock Collet: Rotate the locking collar into the "lock" position (a click can be felt when fully secure, and the "red" dots will be aligned) (Fig. 33).
- Tug slightly on burr to ensure it is fully locked.

**Fig. 33**

Rotate locking collar into "Lock" position.

**N** **NOTE:** If the Milling Handpiece is inoperable, refer to the Milling Handpiece operating instructions for trouble shooting (packaging insert or [brasselerusamedical.com](http://brasselerusamedical.com)). If after trouble shooting the mill remains inoperable, use a second Milling Handpiece. Refer to Appendix A for the backup milling procedure if a second Milling Handpiece is not available.

- When operating the Milling Handpiece, take care to place the slotted end (burr guard) into the track before initiating power on the handpiece (Fig. 34).

**Fig. 34**

Place burr guard into track before initiating power.

**! WARNING:** The Milling Handpiece should only be operated when engaged in the Milling Guide. Serious injury could occur if instructions are not followed.

## BILAG 2

### Uddrag af kirurgisk teknik

#### Appendix A: Backup Milling Procedure

##### Rationale:

The PFJ Drill Milling Adapter can be used to complete the milling operation in the event the Milling Handpiece is inoperable at the time of surgery. The Drill Milling Adapter joins the PFJ Mill Burr to a standard surgical drill via a Jacob's Chuck and utilizes a PFJ Milling Guide for accurate bone preparation. *The milling operation is more challenging when using this approach and should only be used as a backup to the Milling Handpiece.*

##### 1. Assembly:

Insert new PFJ Mill Burr (00-5927-050-00) into the PFJ Drill Milling Adapter as shown in (Fig. 1).

Fig. 1 Insert burr into adapter

With a 1/4 in. or smaller Jacob's Chuck facing upright, insert the shank of the PFJ Mill Burr with PFJ Drill Milling Adapter into the Jacob's Chuck until they both bottom out (Fig. 2). Tighten securely with chuck key (Fig. 3). Rotate assembly to visually verify that the burr is clamped straight and that the PFJ Drill Milling Adapter rotates independently from the PFJ Mill Burr.

Fig. 2 Burr and adapter in Jacob's Chuck

Fig. 3 Tighten Jacob's Chuck securely

**Caution:** To minimize the potential for over-resection, hold Jacob's Chuck vertically during tightening and ensure there is no axial play in the assembly (Figs. 4 & 5).

Fig. 4 Backup Milling Assembly (CORRECT)

Fig. 5 Backup Milling Assembly (INCORRECT)

Place Backup Milling Assembly into a standard surgical drill.

**Caution:** The use of a reamer/wire driver is not recommended. Ensure the standard surgical drill/Jacob's Chuck can achieve 750-25,000 rpm. Too slow/fast may result in poor bone preparation and/or damage to the instruments.

##### 2. Secure the Milling Guide: Additional Fixation

In addition to the three (3) anterior fixation screws, it is recommended to utilize a 48mm MIS headless screw (00-5983-041-48) in the Milling Guide's posterior fixation hole when using this approach. This will help ensure adequate fixation (Fig. 6).

Fig. 6 Posterior screw insertion

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