**Advanced Medical Solutions**  
(Plymouth) Ltd.

Advanced Medical Solutions (Plymouth) Ltd  
Western Wood Way,  
Langage Science Park,  
Plympton, Plymouth,  
Devon, PL7 5BG, UK  
**Web** : [www.admedsol.com](http://www.admedsol.com)

## Field Safety Notice

**Date issued:** 3<sup>rd</sup> March 2017

### Affected products

Product name: LiquiBand® FIX8™  
Product code: FX001

### Description of the situation

Dear valued LiquiBand® FIX8™ user,

Following further clinical evaluation by Advanced Medical Solutions (**AMS**), with effect from May 2015 the indications for use of the LiquiBand® FIX8™ device (**Product**) were extended from “laparoscopic surgical repair of inguinal hernia, achieved through the fixation of polypropylene or polypropylene / polyester combination hernia mesh to the abdominal wall” to “laparoscopic repair of abdominal hernias, achieved through the fixation of prosthetic mesh to the abdominal wall and the approximation of peritoneum” (**Extended Claims**). During a recent routine technical file audit by our Notified Body (BSI), they confirmed that more data was required to support these Extended Claims.

### Risk to health

As Patient Safety is a high priority for AMS, additional clinical data, a review and statement by the field expert in respect of that data, and risk assessments by both AMS and the field expert have already been provided to BSI for review. The following information is consistent with the additional data provided to BSI in respect of the Extended Claims:

- Post Market Surveillance data collected by AMS shows no evidence of adverse events relating to the Extended Claims; and
- the field expert’s review of the current clinical data and our risk assessment concluded that the Product is suitable for the Extended Claims and the risk to the patient and user is low.

Please note that until BSI has completed a review of this additional data, we need to communicate to our users of the Product in all European CE-Mark territories that the Extended Claims are under review and, whilst this review takes place, we will only release Products to customers under the original indications for use.

### Description of the temporary change

AMS will be temporarily removing revision 03 of the Instructions for Use (**IFU**) from use and re-introducing revision 01. As a result, the following sections of the IFU will be affected:

**Advanced Medical Solutions**  
(Plymouth) Ltd.

Advanced Medical Solutions (Plymouth) Ltd  
Western Wood Way,  
Langage Science Park,  
Plympton, Plymouth,  
Devon, PL7 5BG, UK  
**Web** : [www.admedsol.com](http://www.admedsol.com)

	Extended indications under review	Original indications
IFU	DRM 06 0528 revision 03	DRM 06 0528 revision 01
Indications for use	<ul style="list-style-type: none"> <li>The LiquiBand® FIX8™ device is intended for use in laparoscopic repair of abdominal hernias, achieved through the fixation of prosthetic mesh to the abdominal wall and the approximation of peritoneum.</li> </ul>	<ul style="list-style-type: none"> <li>The LiquiBand® FIX8™ HMF device is intended for use in the laparoscopic surgical repair of inguinal hernia, achieved through the fixation of polypropylene or polypropylene / polyester combination hernia mesh to the abdominal wall.</li> </ul>
Precautions	<ul style="list-style-type: none"> <li>There is only limited clinical data for mesh other than polypropylene and titanium coated mesh. Bench tests have demonstrated suitable compatibility with meshes made from polyester and polyvinylidene fluoride (PVDF).</li> </ul>	<ul style="list-style-type: none"> <li>There is only limited medium to long term clinical data available for the safe use of this device for the closure of the peritoneal sac.</li> <li>There is only limited clinical data to support the safe use of this device for the fixation of polyester hernia repair meshes. Bench tests have demonstrated suitable compatibility and performance.</li> </ul>

#### Required actions regarding the use of the Product

- No corrective action needs to take place for patients who received a treatment according to the Extended Claims.
- There is no physical action required regarding any Products you may have in stock. Please disregard revision 03 of the IFU found in the Product packaging and refer to revision 01 of the IFU which is attached to this notice.
- No action is required for any future shipments of the Products as they will be shipped with revision 01 of the IFU until further notice.

#### Required actions regarding this Field Safety Notice

- Please ensure that all those who need to be aware of this notice within your organisation receive a copy of this notice.
- Please transfer this notice to any additional people or organisations who also need to be notified.
- As soon as possible, and **no later than 14 days after receipt of this notice**, you should complete the attached form and return it to AMS either by post or by email to the addresses stated on the form.

#### Contacts

AMS sincerely apologises for any inconvenience caused by this Field Safety Notice. We assure you that we will contact all customers directly and in a timely manner to communicate an update to the

**Advanced Medical Solutions**  
**(Plymouth) Ltd.**

Advanced Medical Solutions (Plymouth) Ltd  
Western Wood Way,  
Langage Science Park,  
Plympton, Plymouth,  
Devon, PL7 5BG, UK  
**Web** : [www.admedsol.com](http://www.admedsol.com)

IFU in respect of the Extended Claims once the review by our Notified Body is complete. In the meantime, please do not hesitate to call your current AMS representative or contact us by email ([Regulatory.Plymouth@admedsol.com](mailto:Regulatory.Plymouth@admedsol.com)) if you have any further questions regarding this notification.

The undersigned confirms that this notice has been notified to the appropriate Regulatory Agencies.

.....  
**Rose Guang**  
**Group QA/RA Director**  
**For and on behalf of Advanced Medical Solutions (Plymouth) Limited**

**Advanced Medical Solutions**  
(Plymouth) Ltd.

Advanced Medical Solutions (Plymouth) Ltd  
Western Wood Way,  
Langage Science Park,  
Plympton, Plymouth,  
Devon, PL7 5BG, UK  
Web : [www.admedsol.com](http://www.admedsol.com)

## Field Safety Notice: Confirmation of receipt form

**Date issued:** 3<sup>rd</sup> March 2017

### Affected products

Product name: LiquiBand® FIX8™  
Product code: FX001

### Customer information

I confirm receipt of this Field Safety Notice and that the information has been communicated to all relevant parties.

Distributor Name  
and Address

---

Hospital / Location

---

Name and Position

---

Contact Email

---

Contact Number

---

Signature

---

Date

---

This form is to be returned **no later than 14 days after receipt of this notice.**

By post to: Advanced Medical Solutions, Western Wood Way, Langage Science Park, Plympton, Plymouth, Devon, PL7 5BG, UK. For the attention of: Rose Guang – Group QA/RA Director.

By e-mail to: [Regulatory.Plymouth@admedsol.com](mailto:Regulatory.Plymouth@admedsol.com).

Registered in England 3603261

VAT No. GB 636 15551 27

Registered Office: Premier Park, 33 Road One, Winsford Industrial Estate, Winsford, Cheshire, CW7 3RT UK

Certificate No. MD 585407

## LiquiBand® FIX8™ Hernia Mesh Adhesive Fixation Device

### Single use Instrument

Read all directions, precautions and warnings prior to use.

These instructions for use provide directions for the proper use of the Advanced Medical Solutions (Plymouth) Ltd. LiquiBand® FIX8™ Hernia Mesh Adhesive Fixation (HMF) device. They are not intended to be a comprehensive instruction manual for the performance of a laparoscopic surgical procedure. Prior to the use of the device, surgeons should be fully proficient in minimally invasive/laparoscopic and open hernia repair surgical techniques.

### Device Description

LiquiBand® FIX8™ HMF device is designed for the application of n-butyl-2-cyanoacrylate adhesive to an implanted hernia repair mesh, in order to fix the mesh to the underlying tissue. The device consists of:

- n-butyl-2-cyanoacrylate adhesive monomer, in liquid form, supplied in a thin-walled, sealed glass vial; and,
- a surgically invasive, laparoscopic delivery instrument comprising a 5mm diameter cannula, with a handle at the proximal end incorporating a loading chamber, filter, piston chamber and trigger. The distal tip of the device is open to allow the adhesive to be dispensed from it.

The device is designed for introduction and use through a 5 mm diameter laparoscopic port sleeve. A larger diameter sleeve will require the use of a converter.

Both the cyanoacrylate adhesive in the glass vial and the surgically invasive delivery device are supplied sterile, for single use only.

### Principle of Operation

The glass vial containing the liquid cyanoacrylate adhesive monomer comes preloaded into, and is subsequently broken in, a loading chamber in the handle of the delivery instrument. The adhesive is drawn through a filter into a piston chamber in the delivery instrument handle. After priming, each press of the trigger on the handle of the delivery instrument dispenses 12.5 mg of adhesive from the distal tip of the cannula. A single device contains sufficient adhesive for at least 33 individual applications of the adhesive. The gauge on the side of the delivery instrument gives an indication of the amount of adhesive delivered and the approximate amount remaining.

### Mechanism of Action

When applied to the proximal surface of the hernia mesh, the liquid adhesive monomer permeates through the perforations in the mesh to the surface of the underlying tissue, where it polymerises on contact with moisture on the tissue surface. This process of chemical polymerisation fixes the mesh to the surface of the tissue at the site of adhesive contact, maintaining the mesh in position while it is incorporated into the abdominal wall through the normal process of tissue fibrosis.

The adhesive completes its polymerisation reaction within approximately 10 seconds. Once completely set, the adhesive no longer possesses adhesive properties such that tissues or surgical instruments may be placed in contact with it without risk of unwanted adhesion.

### Indications for use

The LiquiBand® FIX8™ HMF device is intended for use in the laparoscopic surgical repair of inguinal hernia, achieved through the fixation of polypropylene or polypropylene/polyester combination hernia mesh to the abdominal wall.

### Contraindications

- The device is not intended for use when prosthetic material fixation is contraindicated.
- Do not use on patients with a hypersensitivity to cyanoacrylate adhesives or formaldehyde.
- Do not use for the fixation of meshes constructed with polytetrafluoroethylene (PTFE) or absorbable materials.
- Do not use device for closure or fixation of cerebral tissues, blood vessels or peripheral nerves.

### Warnings

- Only physicians having adequate training and familiarity with endoscopic techniques should use this device. A thorough understanding of the operating principles, risks versus benefits, and hazards involved in utilizing an endoscopic approach is necessary to avoid possible hazards to the user and/or patient.
- This device is provided sterile and is intended for use in a single patient. Do not reuse, reprocess, clean, disinfect or resterilize this device as this may compromise the sterility and performance of the device.
- Ensure the device is properly activated and primed before use according to the instructions in the Directions for Use section.
- This device does not contain any user-serviceable parts. Do not attempt to repair or dismantle the device. If at any point the device appears to be damaged or not functioning correctly, discard and replace it with another device.
- Do not dilute or mix the adhesive with other substances.
- Accidental bonding of unwanted tissue may occur due to misapplication of adhesive. Separation of tissues after accidental bonding should only be performed if deemed necessary.

### Precautions

- The fixation method for any prosthetic device should be determined on the basis of accepted surgical techniques, procedural requirements, and the instruction for use of the prosthetic material.
- Ensure the mesh is held in contact with the underlying tissue during each application of the adhesive and for approximately 10 seconds afterwards, to allow for polymerization.
- The viscosity of the adhesive is only slightly greater than that of water, so adhesive should be applied very carefully to prevent its spread to unwanted areas.
- The adhesive should always be applied in minimal amounts, i.e. avoid multiple applications of adhesive in any given location. A second application of adhesive can be applied over the first only after full polymerization.
- The application of an excessive amount of adhesive in a single location prolongs setting time and can prevent adherence. After polymerization, any excess adhesive may lead to detachment of the adhesive film and/or give rise to the formation of small fragments of polymerised adhesive.


- The adhesive will readily adhere to most substrates. Care should be taken to avoid unwanted contact with the adhesive during polymerization. Polymerized adhesive can be removed from metal instruments with acetone.
- There is only limited medium to long term clinical data available for the safe use of this device for the closure of the peritoneal sac.
- There is only limited clinical data to support the safe use of this device for the fixation of polyester hernia repair meshes. Bench tests have demonstrated suitable compatibility and performance.

### Potential Adverse Reactions

As with the majority of implanted devices, adverse reactions associated with the use of this device may include transient local irritation at the implant site and a transitory inflammatory foreign body response.

### Directions for Use

#### Priming the device

The device must be properly activated and primed to ensure that it operates correctly.

- Remove the device from its packaging. Hold the device with the handle vertical and the distal tip of the cannula tilted downwards.

Lift the blue lid firmly upward until a cracking sound of the glass ampoule breaking inside the lid is heard.

- Wait for 5 seconds then, whilst the tip pointing downwards and handle held vertically slowly pull out red tab no. 2, observing that the adhesive is being drawn into the transparent barrel. This should take approximately 20 seconds. At the end of piston travel gently disengage the red tab and dispose of.

Warning: A fast transfer stage could result in excessive air being drawn into the dispensing chamber.

- Gently close the blue lid until it is completely aligned with the body and clicks into place.
- Remove the red pull tab no. 4 and dispose of.
- Rotate the blue knob clockwise approximately 320° until blue dispenser trigger alongside the handle is released. Care should be taken to avoid impeding the release of the trigger. The distal tip should be directed into the packing tray well no. 6 before the blue knob is rotated.
- With the cannula tip placed into packing tray well no. 6, prime the device by actuating the trigger, approximately three times, until a droplet of adhesive freely emerges.
- Wipe any excess adhesive from the tip of the cannula onto the inside of the packing tray or on a sterile wipe. The device is now ready for surgical use.

### Using the Device

- Adhesive should be applied in discrete locations distributed across the mesh to fix the mesh in place.
- Upon complete fixation of the mesh, surgical graspers should be used to test for the adhesion of the mesh and for any areas left unadhered.

### Disposal

After use, dispose of the device according to local procedures and guidelines.

### Storage

This device should always be stored in its original packaging. Store only at temperatures between 5°C (41°F) and 25°C (77°F). Do not use the device after the expiration date shown on the blister pack.

### Sterility

This device is packaged for single patient use and is sterilized by irradiation and ethylene oxide gas. Do not re-sterilize. Do not use if package is open or damaged.

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