ADVIA®-kemi systemer

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ADVIA is a trademark of Siemens Healthcare Diagnostics.

Appendix:

Maximum % bias observed at 300 mg/L of sulfasalazine and sulfapyridine

Assay	Concentration of analyte in common unit (SI Unit)	Maximum% bias observed at 300 mg/L Sulfasalazine	Maximum% bias observed at 300 mg/L Sulfapyridine
Ammonia (AMM)	~60 µg/dL (35 umol/L)	75.9%	-18%
Salicylate (SAL)	~25 mg/dL (1.8 mmol/L)	25%	-24.8%
Alanine Aminotransferase (ALT) and equivalent assay ALT_c	~50 U/L	-25.6%	Interference ≤ 10% observed
ALTP5P , and equivalent assay ALTP_c	~50 U/L	-34.6%	Interference ≤ 10% observed

FIELD CORRECTION EFFECTIVENESS CHECK

Sulfasalazine and Sulfapyridine Interference with NADH and/or NADPH reaction Assays

This response form is to confirm receipt of the enclosed Siemens Healthcare Diagnostics Urgent Field Safety Notice CHC17-06.A.OUS.CHC dated June, 2017 regarding Sulfasalazine and Sulfapyridine Interference with NADH and/or NADPH reaction Assays. Please read the question and indicate the appropriate answer. Fax this completed form to Siemens Healthcare Diagnostics at the fax number provided at the bottom of this page.

1. I have read and understood the Urgent Field Safety Notice Correction instructions provided in this letter. Yes No

Name of person completing questionnaire:

Title:

Institution:

Instrument Serial Number:

Street:

City:

State:

Phone:

Country:

Customer Sold To #:

Customer Ship To #:

Please fax this completed form to the Customer Care Center at (###) ###-####. If you have any questions, contact your local Siemens technical support representative.