Dimension® klinisk kemi system og Dimension Vista® System

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Dimension and Dimension Vista are trademarks of Siemens Healthcare Diagnostics.

Sulfasalazine and Sulfapyridine Interference with NADH and/or NADPH reaction Assays

Appendix:

Maximum bias observed for sulfasalazine and sulfapyridine with Dimension Assays

Assay	Concentration of analyte	Maximum% bias observed at 300 mg/L [0.75 mmol/L] Sulfasalazine	Maximum% bias observed at 300 mg/L [1.2 mmol/L] Sulfapyridine
Ammonia (AMM)	~426 µg/dL [~250 µmol/L]	<10%	-19%
Alanine Aminotransferase (ALTI)	~55 U/L [~0.92 µkat/L]	-29%	<10%
Aspartate Aminotransferase (AST)	~37 U/L	-10%	<10%
Glucose (GLUC)	~126 mg/dL [~7.0 mmol/L]	-17%	11%
Creatine Kinase MB (MBI)	~20 U/L [0.33 µkat/L]	<10%	-11%
Thyroxine (T4)	~8 µg/dL [~103 nmol/L]	15%	<10%

Maximum bias observed for sulfasalazine and sulfapyridine with Dimension Vista Assays

Assay	Concentration of analyte	Maximum% bias observed at 300 mg/L [0.75 mmol/L] Sulfasalazine	Maximum% bias observed at 300 mg/L [1.2 mmol/L] Sulfapyridine
Ammonia (AMM)	~426 µg/dL [~250 µmol/L]	66%	-19%
Alanine Aminotransferase (ALTI)	~55 U/L [~0.92 µkat/L]	-72%	-19%
Aspartate Aminotransferase (AST)	~37 U/L	-19%	-12%
Glucose (GLU)	~126 mg/dL [~7.0 mmol/L]	-21%	11%
Creatine Kinase MB (MBI)	~20 U/L [0.33 µkat/L]	22%	<10%
Thyroxine (T4)	~8 µg/dL [~103 nmol/L]	18%	<10%

FIELD CORRECTION EFFECTIVENESS CHECK

Sulfasalazine and Sulfapyridine Interference with NADH and/or NADPH reaction Assays

This response form is to confirm receipt of the enclosed Siemens Healthcare Diagnostics Urgent Field Safety Notice CHC-17-06.A.OUS.DM dated June 2017 regarding Sulfasalazine and Sulfapyridine Interference with NADH and/or NADPH reaction Assays. Please read the question and indicate the appropriate answer. Fax this completed form to Siemens Healthcare Diagnostics at the fax number provided at the bottom of this page.

1. I have read and understood the Urgent Field Safety Notice instructions provided in this letter. Yes No

Name of person completing questionnaire:

Title:

Institution:

Instrument Serial Number:

Street:

City:

State:

Phone:

Country:

Customer Sold To #:

Customer Ship To #:

Please fax this completed form to your local Siemens technical support representative. If you have any questions, contact your local Siemens technical support representative.