Urgent Field Safety Notice (FSN) NephroCam (Gamma Camera) Risk of patient data mix-up

Language: English / Dansk

2017-12-13

Dear valued customer

The purpose of this letter is to notify users of **NephroCam** systems that under specific conditions there is a risk of patient data mix-up on PACS or workstation.

The purpose of this document is to inform you of an issue which has been brought to our attention, the immediate actions required to mitigate the issue, and the corrective actions which DDD-Diagnostic will take to solve the issue.

This document contains important information regarding the correct use of your equipment. This notice needs to be passed on to everyone who needs to be aware of the issue within your organization and in any other organization to which this issue may apply.

Please keep a copy of this notice together with the Instructions for Use.

Summary

- **Issue:** DDD has identified that the DICOM Study Instance UIDs used by NephroCam systems are not globally unique as required by the DICOM standard. If acquisition data from multiple NephroCam systems is transferred to the same PACS or workstation, this can lead to incorrect storage of acquisition data under another patient's name.
- **Resolution:** DDD will promptly release a software update (version 1.2.2) to ensure that DICOM Study Instance UIDs are globally unique.
- **Precautions:** NephroCam is still considered safe to use as long as patient data is not transferred to PACS or a workstation holding patient data from other NephroCam systems. If data from multiple **NephroCams** has been transferred to the same PACS or workstation, we recommend reviewing of previous patient data to reduce the risk of data mix-up caused by this issue.

- **Next Step:** Our DDD service engineer will contact you in the near future to arrange for the software update and closure of the FSCA (Field Service Corrective Action)

Transmission of this Field Safety Notice:

- This notice needs to be passed on to everyone who needs to be aware of the issue concerned within your organization and any other organization to which this issue may apply.
- Please maintain awareness on this notice and the resulting action for an appropriate period of time to ensure effectiveness of the corrective action.

Affected Products: all NephroCam Systems

NephroCam Systems (Model Number: 9KID2275-B02) with the following serial numbers:

13112201, 13112202, 13112203, 14111001, 15061701, 15090701, 15121001, 16062901, 16070401, 16120601, 16122101

Contact Information

If you have any questions regarding this Field Safety Notice, please feel free to contact your Service and Support Manager Albert Radoor or our QA / RA Manager Stephane Labrosse:

<p>Service and Support Manager Albert Rømer Radoor Tel.: +45 2386 7654 Email: support@ddd-diagnostic.dk</p>	<p>QA / RA Manager Stephane Labrosse Tel.: +45 6058 8020 Email: sl@ddd-diagnostic.dk</p>
---	---

We sincerely apologize for any inconvenience this may cause. Please contact us if you have any questions or concerns.

With kind regards,

Stephane Labrosse
Quality Assurance and Regulatory Affairs Manager

Vigtig sikkerhedsmeddelelse (Field Service Notice) NephroCam (Gammakamera) Risiko for sammenblanding af patientdata

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