

**COOK®**

**Cook Medical Europe**  
 O'Halloran Road,  
 National Technological Park,  
 Limerick, Ireland.  
 Phone: + 353 61 334440  
 Fax: + 353 61 334441

**Urgent Field Safety Notice**

**Commercial name of the affected product:** Entuit® Secure Gastrointestinal Suture Anchor Set  
**Manufacturer :** Cook Incorporated,  
**Cook Reference Number:** 2018FA0010  
**Type of action:** Field Safety Corrective Action (FSCA)

-----  
 Date: 13 Dec 2018

Attention: Chief Executive / Risk Management / Purchasing

**Details on affected devices:**

PRODUCT BRAND NAME	REFERENCE PART NUMBER (RPN)	ORDER NUMBER	LOT NUMBER
Entuit® Secure Gastrointestinal Suture Anchor Set	GIAS-SRM-2	G35569	8144390, 8159665, 8173295, 8204594, 8261263, 8294836, 8308856, 8331166, 8395981, 8459751, 8471244, 8511162, 8535303, 8541938, 8574041, 8582863, 8597539, 8541938X, 8582863X
	GIAS-SRM-3	G35570	8056941, 8135292, 8165522, 8228910, 8261265, 8293447, 8331165, 8357394, 8389427, 8395982, 8445405, 8471250, 8483543, 8511203, 8556989, 8557006, 8564659, 8622758, 8628912, 8663134

**Description of the problem:**

Cook has received complaints for difficulty in sliding down the retention mechanism on Entuit® Secure Gastrointestinal Suture Anchor Set manufactured with a specific extension spring lot. Therefore, Cook is initiating a voluntary recall of the 39 lots of Entuit® Secure Gastrointestinal Suture Anchor Set that were manufactured with the affected extension spring lot.

Potential adverse events that may occur if an affected product is used include delayed or prolonged procedure, additional intervention to place more devices, and prolonged hospitalization.

**Advise on action to be taken by the user:**

1. Immediately collect all remaining affected products as per the specified lot listing from your inventory.
2. Please complete the enclosed Customer Response Form. Where product is indicated as being returned, our Customer Services department will contact you to organize the return and issue you

with the relevant Returns Authorization number. Please include contact details on the Customer Response form.

Product should be addressed to:

Cook Medical EUDC  
Robert-Koch-Straße, 2  
52499 Baesweiler  
GERMANY

Credit will be provided for the returned affected products where applicable.

3. Send the Customer Response Form via email to [European.FieldAction@CookMedical.com](mailto:European.FieldAction@CookMedical.com) or alternatively by fax to Cook Medical marked for the attention of European Customer Quality Assurance (fax number +353 61239294). Do not enclose the response form with the returned product.
4. Please report any adverse events to Cook Medical by contacting our Customer Support Department.

**Transmission of this Field Safety Notice:**

This notice needs to be passed on all those who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred. (If appropriate)

Please transfer this notice to other organisations on which this action has an impact. (If appropriate)

Please maintain awareness on this notice and resulting action for an appropriate period to ensure effectiveness of the corrective action. (if appropriate)

**Contact reference person:**

Larry Pool  
Post Market Director  
Cook Incorporated  
50 Daniels Way, PO Box 489, Bloomington, IN 47402, United States

The undersign confirms that this notice has been notified to the appropriate Regulatory Agency.

Should you have any questions, please feel free to contact us for more information (e-mail: [European.FieldAction@CookMedical.com](mailto:European.FieldAction@CookMedical.com), phone +353 61 334440).



Larry Pool  
Post Market Director