

January 30, 2019

To: Hospital

Subject: **URGENT MEDICAL DEVICE FIELD SAFETY NOTICE - REMOVAL**

Reference: ZFA2019-00348

Affected Product: Zimmer® Sterile Disposable Tourniquet Cuff, Dual Bladder

Item Number	Description	Lot Number
60-7080-152-00	Zimmer® Sterile Disposable Tourniquet Cuff with Protective Sleeve and PLC, Dual Port, Dual Bladder	27968701



Zimmer Surgical Inc. is conducting a medical device Field Safety Corrective Action (removal) for one lot of the Zimmer® Sterile Disposable Tourniquet Cuff, Dual Bladder due to a potential leak between the two bladders.

The issue can be recognized easily when inflating/deflating a single bladder, both bladders could potentially inflate/deflate due to the leak. Additionally, alarms would sound in the event of an unintended pressure change in either cuff bladder. To date no adverse effects are reported.

Risks		
Describe immediate health consequences (injuries or illness) that may result from use of or exposure to the product issue.	Most Probable	Highest Severity
	<i>None</i>	<i>Exposure to Anesthesia, Critical - Local Anesthetic Systemic Toxicity (L.A.S.T)</i>
Describe long range health consequences (injuries or illness) that may result from use of or exposure to the product issue.	Most Probable	Highest Severity
	<i>None</i>	<i>None</i>

Our records indicate that you may have received one or more of the potentially affected units. The potentially affected units were distributed between March 2019 and May 2019 (Local deployment may differ).

Hospital Responsibilities:

1. Review this Field Safety Notice and ensure that affected personnel are aware of the contents.
2. If you have any potentially affected units at your facility, assist your Zimmer Biomet sales representative and quarantine all potentially affected units. Your Zimmer Biomet sales representative will remove the potentially affected units from your facility.
3. Complete **Attachment 1 – Certificate of Acknowledgement** and send to fieldaction.emea@zimmerbiomet.com. This form must be returned even if you do not have potentially affected units at your facility.
4. Retain a copy of the acknowledgement form with your field action records in the event of a compliance audit of your facility's documentation.
5. If you have further questions or concerns after reviewing this Field Safety Notice, please contact your Zimmer Biomet representative.

Other Information


This medical device Field Safety Notice was reported to all relevant Competent Authorities and the related Notified Body as required under the applicable regulations for Medical Devices as per MEDDEV 2.12-1 in Europe.

Please keep Zimmer Biomet informed of any adverse events associated with this units or any other Zimmer Biomet product by emailing winterthur.per@zimmerbiomet.com or to your local Zimmer Biomet contact.

Please be aware that the names of user facilities notified are routinely provided to the Competent Authorities for audit purposes. The undersigned confirms that this notice has been delivered to the appropriate Regulatory Agencies.

We would like to thank you for your co-operation in advance and regret any inconveniences caused by this field action.

Sincerely,



Kevin W. Escapule
Post Market Surveillance & Regulatory Compliance Director



ATTACHMENT 1 Certificate of Acknowledgement

IMMEDIATE RESPONSE REQUIRED – TIME SENSITIVE ACTION NEEDED

Affected Product: Zimmer® Sterile Disposable Tourniquet Cuff, Dual Bladder

Field Action Reference: ZFA 2019-00348

Please return the completed form to your Zimmer Biomet contact person or by e-mail fieldaction.emea@zimmerbiomet.com

I received and understood the Field Safety Notice.

Regarding the parts:

All inventories for the potentially affected units have been checked and following parts are to be returned:

Item Reference	Lot Number	Number of parts returned

OR

The potentially affected units which are unavailable for return have been used

By signing below, I acknowledge that the required actions have been taken in accordance with the Field Safety Notice.

Hospital Facility **Surgeon** *(Please check one as applicable)*

Printed Name: _____ **Signature:** _____ **Date:** ____/____/____

Title: _____ **Telephone:** () ____ - _____

Facility Name: _____ **Facility Address:** _____

City: _____ **ZIP:** _____ **Country:** _____