«Hospital_Name»

«Users_Name»

«Department»

«Customer_Address»

«Zip_Code» «City»

«Country_name»

<Reference: 92629317-FA>

1 December 2020

Urgent Field Safety Notice - Urgent Medical Device Recall AXIOS Stent and Electrocautery - Enhanced Delivery System, and Hot AXIOS Stent and Electrocautery - Enhanced Delivery System

Dear «Users_Name»,

Boston Scientific Corporation (BSC) is initiating a voluntary removal of specific lots/batches of the AXIOS Stent and Electrocautery Enhanced Delivery System, and Hot AXIOS Stent and Electrocautery Enhanced Delivery System. During manufacturing BSC identified that the rotating luer fitting at the distal end of the delivery system handle may become detached from the nose of the delivery system handle, which can occur during device preparation, stent delivery, or delivery system removal. This issue only affects the delivery system of the device and there is no effect on implanted stents. There have been no complaints reported to date related to this issue.

The **most common** potential consequences due to delivery system handle separation during AXIOS stent placement for a pseudocyst or walled-off necrosis are a prolongation of procedure in order to exchange out for another device or an additional intervention to treat the patient and/or complete the procedure. In addition, there is a risk of infection as a result of the opening created by the AXIOS delivery system or from additional manipulation of the pseudocyst, walled-off necrosis or biliary tract that is not adherent to the gastric or bowel wall. The **most severe** potential consequences due to delivery system handle separation during AXIOS stent placement are peritonitis and need for surgical intervention. The likelihood of the most severe potential consequence occurring is remote.

Our records indicate that your facility received some of the concerned product. **The table below provides a complete list of all affected products**, including Product Description, Material Number (UPN), GTIN? Lot/Batch numbers and expiry date. Please note that **only the devices listed below are affected. No other Boston Scientific product is involved in this Field Safety Notice.**

Further distribution or use of any remaining product affected by this action should cease immediately.

| Product Description | US or Outside US | UPN # | GTIN | Lot/Batch # | Expiration Date Range |
|--|------------------|-----------|----------------|--|------------------------|
| Hot AXIOS™ Stent and Electrocautery – Enhanced Delivery System | Outside US | M00553540 | 08714729904564 | 26130517 | 10/5/2022 |
| | Outside US | M00553550 | 08714729904571 | 26156059, 26167275, 26167276, 26156742, 26156740 | 10/9/2022 – 10/12/2022 |
| | Outside US | M00553560 | 08714729951100 | 26167277, 26171099, 26152171, 26171383, 26171380, 26152172 | 10/8/2021 – 10/12/2021 |
| AXIOS™ Stent and Electrocautery - Enhanced Delivery System | US | M00553640 | 08714729904588 | 26134666 | 10/6/2022 |
| | US | M00553650 | 08714729904595 | 26152465, 26152466, 26152467 | 10/8/2022 |
| | US | M00553660 | 08714729951179 | 26167271, 26152176, 26167274 | 10/8/2021 – 10/12/2021 |

INSTRUCTIONS:

1- **Please immediately discontinue use of the Boston Scientific product reported in the list and remove all of the affected units from your inventory**, regardless of where these units are stored in your facility. Segregate the units in a secure place, pending return to Boston Scientific.

2- **Please complete the attached Verification Form** even if you do not have any product to return.

3- **When completed, please return the Verification Form to your local Boston Scientific office** for the attention of «Customer_Service_Fax_Number» on or before **24 December 2020.**

4- **If you have products to return**, please package them in an appropriate shipping box and **contact «Customer_Service_Tel» of your local Boston Scientific office**, to arrange return.

5- Please pass this notice to any healthcare professional from your organization that needs to be aware and to any organization where the potentially affected devices have been transferred (if appropriate). Please provide Boston Scientific with details of any affected devices that have been transferred to other organizations (if appropriate).

Your Competent Authority is being notified of this Field Safety Notice.

We regret any inconvenience that this action may cause, and we appreciate your understanding as we act to ensure patient safety and customer satisfaction.

If you have any questions or would like assistance with this Field Safety Notice, please contact your local Sales Representative.

Yours sincerely,

Marie Pierre Barlangua
Quality Department
Boston Scientific International S.A.

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