FSN-reference: FSN 20-001
FSCA-reference: FSCA 20-001

URGENT: FIELD SAFETY NOTICE

WiSE CRT Transmitter Model 4100 - potential premature Battery depletion

03 March, 2021

Dear _____

EBR Systems, Inc is issuing a **Field Safety Notice** on **WiSE CRT Transmitter Model 4100**. **Competent Authorities are aware of this action.**

Our records indicate that you have implanted the WiSE CRT Transmitter Model 4100 and may have received the product serial numbers subject to this Field Safety Notice identified in **Appendix B and C**.

This Field Safety Notice is **EFFECTIVE IMMEDIATELY – DO NOT USE** identified affected products.

PLEASE DISTRIBUTE THIS INFORMATION TO ALL PERSONNEL RESPONSIBLE FOR or WHO MAY USE the WiSE CRT Transmitter Model 4100 in your facility.

EBR Systems, Inc has initiated a voluntary Field Safety Corrective Action (FSCA) for WiSE CRT Transmitter Model 4100 because of increasing rate of premature Battery depletion associated with the WiSE CRT Transmitter Model 4100. If this occurs, the System will continue to function normally and deliver biventricular pacing until the Battery is depleted. Other WiSE CRT system components are unaffected by this action.

Product intended use: The WiSE CRT System is an implantable wireless cardiac pacing system intended for the leadless, endocardial stimulation of the left ventricle in patients indicated for cardiac resynchronization therapy. The WiSE CRT Transmitter Model 4100 is a subcutaneously implanted ultrasound transmitter that initiates an ultrasonic energy pulse that travels through the tissue to intersect an ultrasound receiver implanted in the heart.

A total of 8 Transmitters were explanted, returned for analysis and replaced with new devices. One Battery pocket infection was associated with these revision surgeries. There were no other clinical sequelae reported. Analysis of these 8 devices confirmed the failure mode was an insulation breach in the Transmitter feedthrough, resulting in a failure rate of 2.8% (8 out of 288 devices implanted globally). The time to loss of therapy for these 8 patients ranged from 7.1 to 22.6 months. An additional 8 Transmitters are suspected of having a similar issue that could raise the rate to as high as 5.6%. Clinical impact (patient risk) is provided as **Appendix D**.

For these reasons, EBR Systems will pause further shipments of the Transmitter Model 4100 for patients until product which do not have propensity for premature Battery depletion are available.

Although investigation to date does not indicate all units are prone to premature Battery depletions, users of WiSE CRT Transmitter Model 4100 who have un-implanted units are advised to not implant these units as a precaution. EBR Systems, Inc will continue to monitor any adverse events related to the issue.

FSN-reference: FSN 20-001
FSCA-reference: FSCA 20-001

Hence, all healthcare professionals, affected staff, service and/or facilities working with/using the WiSE CRT System should refer to the attached Field Safety Notice, which includes:

- i. Patient management recommendations (Appendix A)
- ii. A list of potentially affected devices implanted by your hospital (Appendix B)
- iii. A list of devices currently at your site to be returned (Appendix C)
- iv. Clinical impact (Appendix D) and
- v. An Acknowledgement Form (Appendix E) that will need to be completed and returned via email to compliance@ebrsystemsinc.com or to your local EBR representative.

ACTION REQUIRED: NEXT STEPS and IDENTIFICATION OF THE PRODUCT SERIAL NUMBERS SUBJECT TO THIS FIELD SAFETY NOTICE:

1. Refer to **Appendix C** for assistance in identifying the product serial numbers subject to this Field Safety Notice. Examine your inventory immediately to determine if you have any products subject to this Field Safety Notice in stores or at points of use and quarantine such product(s).
2. Remove the products subject to this Field Safety Notice from your inventory and communicate the issue to all relevant cath lab or materials management personnel, or anyone else in your facility who needs to be informed.
3. Complete relevant **Acknowledgement Form (Appendix E)** confirming receipt of this notice and return it to EBR Systems, Inc immediately upon receipt. Please return Acknowledgement Form (Appendix E) even if you do not have the product serial numbers subject to this Field Safety Notice. Please feel free to contact EBR Systems, Inc, local EBR representative, EU Authorised Representative (for affected products in the EU) or UK Responsible Person (for affected products in the UK) if you require any assistance completing the **Acknowledgement Form (Appendix E)**.
4. Customers are required to immediately return all unused WiSE CRT Transmitter Model 4100 subject to this Field Safety Notice that are in their inventory. EBR Systems, Inc will provide detailed instructions for product return once **Acknowledgement Form (Appendix E)** received by customer service. Additionally, EBR Systems, Inc will supply product for any patients requiring device replacement based on your clinical assessment.

If you have additional questions regarding this Field Safety Notice, please contact EBR Systems, Inc via the contact details below:

Manufacturer:

EBR Systems, Inc

E-mail: support@ebrsystemsinc.com

Telephone: +1 408.720.1906

Local Representative

Name: Tobias Wilde

E-mail: tobias@ebrsystemsinc.com

Telephone: +49 172 263 165 7

Alternatively, please contact UK Responsible Person (affected product in the UK) and EU Authorised Representative (affected product in the EU):

FSN-reference: FSN 20-001
FSCA-reference: FSCA 20-001

EU Authorised Representative

Name: Advena Ltd

Address: Tower Business Centre, 2nd Flr., Tower Street, Swatar, BRK 4013, Malta

E-mail: info@advenamedical.com

Telephone: +356 2546 6689

UK Responsible Person

Name: Advena Limited

Address: Pure Offices, Plato Close, Tachbrook Park, Warwick, CV34 6WE, UK

E-mail: info@advenamedical.com

Telephone: +44 1926 800 153

We regret any inconvenience that this action may cause, but we appreciate your understanding as we take action to ensure patient safety and customer satisfaction.

Sincerely,

Madhuri Bhat

Madhuri Bhat
Senior Vice-President
Clinical, Quality, Regulatory & Compliance

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