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Vi beklager den ulempe dette måtte medføre.

Venlig hilsen  
SINOVO health solutions GmbH/Ypsomed AG

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## Important notice

### to mylife™ App users with a connected Dexcom G6

We would like to inform you of a software defect in our product called “mylife™ App” that might have an impact on patient safety. You are contacted because you may be an active user of the mylife™ App and with a Dexcom G6 Continuous Glucose Monitoring (CGM) system connected to it.

The following notice is sent out in coordination with responsible health authorities and on behalf of SINOVO health solutions GmbH (legal manufacturer of the product) and Ypsomed AG (distributor of the product).

#### Affected users

The identified defect concerns persons with diabetes who use the mylife™ App as follows. The conditions are cumulative.

- Users on pump therapy (CSII) or pen therapy (MDI)
- Users with a Dexcom G6 CGM system that is connected with the mylife™ App
- Users with the country setting of Germany, Switzerland, Austria or Denmark

#### Affected installations

The affected installations of the mylife™ App are limited as follows:

- mylife™ App version **1.8.0** (UDI: 111652811196.V1.8.0; you can view the installed version of the app in the About menu of the app)
- Installations on Android smartphones only (any **Android** OS version)

#### Description of the safety-relevant defect

Due to the defect and under certain conditions, the mylife™ App is not able to consistently recognise a loss of the Bluetooth® connection with the Dexcom G6 transmitter. Instead of displaying a signal loss alert on the Last Values screen (home screen), the app may continue to show the last received sensor reading and trend arrow from the Dexcom G6. This last sensor reading and trend arrow, even when outdated, may keep showing over an extended period of time without changing.

Note that the bolus calculator in the mylife™ App is **not** affected of this defect. As per design, the bolus calculator will **not** accept any outdated CGM values as input for a bolus calculation.

#### How can you recognise this situation?

- You observe the exact same sensor reading (CGM value) and trend arrow over an extended period of time.
- The CGM trend graph is not filled with sensor readings anymore. Where normally you can see the last 3 hours of CGM values, you now see a growing gap, starting from the right-hand side of the graph.

### How should you react in this situation?

- Use a blood glucose meter and fingersticks to monitor your blood glucose levels and to make treatment decisions.
- Shut down the mylife™ App completely (not just in the background) and restart it after 5 minutes.
- After the restart, accurately observe the incoming sensor readings from the Dexcom G6. Pay attention if the values are changing every 5 to 10 minutes and if the trend graph is filled with the readings.

### Countermeasures and next steps

We are working with full focus on the correction of this software defect. The corrected mylife™ App is targeted for release in the Google Play Store within the next 2–3 weeks and for the affected countries of Germany, Switzerland, Austria and Denmark.

Once the new app release is available, you will receive a notification in the app to update your current installation. Perform the update immediately when you are notified.

### Contact details

If you have any question regarding this notice, do not hesitate to contact your local mylife™ Diabetescare representative.

**Danmark** Ypsomed ApS // Smedeland 7 // 2600 Glostrup //  
info@ypsomed.dk // www.mylife-diabetescare.dk // +45 48 24 00 45

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We apologise for any inconvenience that this might cause you.

Yours sincerely  
SINOVO health solutions GmbH/Ypsomed AG

#### Disclaimer:

Dexcom and Dexcom G6 are registered trademarks of Dexcom, Inc. in the United States and/or other countries.  
The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc., and any use of such marks by Ypsomed is under license.  
Other trademarks and trade names are those of their respective owners.

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**Deutschland** Ypsomed GmbH // Höchster Straße 70 // 65835 Liederbach //
**Österreich** Ypsomed GmbH // Am Euro Platz 2 // 1120 Wien //
**Schweiz** Ypsomed AG // Markt Schweiz // Weissensteinstrasse 26 // 4503 Solothurn //