«Hospital_Name»

«Users_Name»

«Department»

«Customer_Address»

«Zip_Code» «City»

«Country_name»

<Reference: 97090504-FA>

3 October 2023

Urgent Field Safety Notice - Urgent Medical Device Recall EXALT™ Model D Single-Use Duodenoscope

Dear «Users_Name»,

Boston Scientific is recalling certain batches of the EXALT Model D Single-Use Duodenoscope following an increase in reports of poor image quality due to fluid ingress in the lens. This is a known, uncommon malfunction that has been reported as occurring at a higher rate than anticipated in certain batches of the EXALT Model D scope. In an effort to address this issue, Boston Scientific has implemented corrective changes and will make replacement EXALT Model D scopes available as soon as possible.

If the user does not have a clear, unobstructed, live endoscopic view, it is anticipated that, in accordance with the warnings in the Instructions For Use, the user will not insert, advance, or further operate the device, as doing so could cause patient injury, such as perforation, bleeding, or tissue damage.

The most severe health consequences reasonably foreseeable to occur if a user experiences poor image quality with the inability to identify anatomical structures are bleeding, tissue damage, and removal of the scope and accessories when accessories have been advanced into the biliary and pancreatic ducts, requiring additional interventions and starting the procedure from the beginning.

The most common adverse health consequence reasonably foreseeable to occur if a user experiences poor image quality with the inability to identify anatomical structures is a minor prolongation of the procedure while the device is exchanged for another, or the procedure is continued with poor image quality.

The issue cannot be reliably detected prior to use and cannot be corrected following detection.

Our records indicate that your facility received some of the concerned product. **The table below provides a complete list of all affected products, including Product Description, Material Number (UPN), GTIN, Lot/Batch numbers.** Please note that **only the devices listed below are affected. No other Boston Scientific product is involved in this Field Safety Notice.**

Further distribution or use of any remaining product affected by this action should cease immediately.

PLEASE NOTE: We are aware that hospitals often remove products from the outer carton and store on the shelves in the inner-pouch only. If this is a practice at your facility, **it is very important that you carefully use the product table and consider both the inner and outer packaging UPN codes when searching for affected product, as the UPN numbers on the inner and outer labelling may be different. The product information listed on your specific Verification Form (enclosed with this letter) provides outer package product coding only** and should be utilized when reporting product to return.

Verify by product batch/lot number in the product table to determine if the batch within your inventory is affected. If so, indicate on your Verification Form the quantity of units from each batch that you will be returning. **As the product within these batches are sold as Single packs and 2-packs, it is important that all reported quantities represent the actual number of single unit being returned and not the number of cartons/boxes or multi-packs.**

Box of 1 unit: UPN and Batches

Product Description	Outer Box UPN #	GTIN	Lot / Batch #		Expiration Date Range
EXALT™ Model D Single-Use Duodenoscope	M0054242CE0	8714729995746	31199348	31712358	23 February 2025 - 24 July 2025
	M0054242CE0	8714729995746	31574944		

Box of 2 units: UPN and Batches

Product Description	Outer Box UPN # (2 Pack)	Inner Box UPN # (Single)	GTIN	Lot / Batch #		Expiration Date Range
EXALT™ Model D Single-Use Duodenoscope	M0054242CE1	M0054242CE0	8714729995753	31118201	31331654	23 February 2025 - 13 June 2025
	M0054242CE1	M0054242CE0	8714729995753	31129503	31817723	
	M0054242CE1	M0054242CE0	8714729995753	31265893		

INSTRUCTIONS:

1- **Please immediately discontinue use of the Boston Scientific product reported in the list and remove all of the affected units from your inventory**, regardless of where these units are stored in your facility. Segregate the units in a secure place, pending return to Boston Scientific.

2- **Please complete the attached Verification Form even if you do not have any product to return.**

3- **When completed, please return the Verification Form to your local Boston Scientific office** for the attention of «Customer_Service_Fax_Number», on or before **23 October 2023**.

4- **If you have products to return, please package them in an appropriate shipping box. After receipt of the Verification Form, Boston Scientific will contact you to arrange return.**

5- Please pass this notice to any healthcare professional from your organization that needs to be aware and to any organization where the potentially affected devices have been transferred (if appropriate). Please provide Boston Scientific with details of any affected devices that have been transferred to other organizations (if appropriate).

Your Competent Authority is being notified of this Field Safety Notice.

We regret any inconvenience that this action may cause, and we appreciate your understanding as we act to ensure patient safety and customer satisfaction.

If you have any questions or would like assistance with this Field Safety Notice, please contact your local Sales Representative.

Yours sincerely,

Conor Dolan
Vice President, Global Quality

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Boston Scientific International S.A.