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Hello there,

We hope you are doing well.

We are reaching out to you to inform you about a relevant topic related to the mySugr Bolus Calculator. We have recently identified a bug affecting the mySugr Bolus Calculator on iOS devices. **Good news:** The bug is already fixed!

Please read below about the steps we kindly ask you to take, like updating the mySugr app to its latest version.

What was the issue and who was affected?

During internal testing on July 16, 2024, we discovered that under rare circumstances an Insulin Correction Factor (ICF) value entered in the mySugr app on iOS in mmol/L unit may be incorrectly applied to the mySugr Bolus Calculator. This may lead to insulin recommendations issued by the mySugr Bolus Calculator to be higher than intended. Using an incorrect ICF (in this case a very low value) for a bolus calculation may lead to a serious incident, if the recommended insulin amount, which would be too high, is injected.

The bug was available since mySugr Bolus Calculator version 3.0.0 (UDI:(01)9120084260044(10)IV3.0.0).

What we ask you to do

Please immediately update the mySugr app to the latest version (3.99.1). After the update, you need to carefully review your mySugr Bolus Calculator settings and pay special attention to the Insulin Correction Factor (ICF) to ensure it is set as you intend. After modifying your mySugr Bolus Calculator settings, you should always carefully review them before administering insulin based on a bolus recommendation for the first time and periodically thereafter. If the mySugr Bolus Calculator is not available to you, please reach out to our customer support.

What does mySugr do to deliver safe products?

mySugr as a medical device manufacturer applies international standards such as ISO 13485 for Quality Management. The company and its medical devices are regularly audited by TÜV SÜD to demonstrate compliance to these international standards and regulations. We have strict, multi-stage development and testing processes to prevent errors of this kind. Nevertheless, an error has occurred and we will do everything we can to further optimize our quality assurance to prevent something like this from happening in the future.

We take this incident very seriously and apologise for any inconvenience this may have caused.

If you have any questions or problems with your mySugr app, our Customer Support team at support@mysugr.com will be happy to help you!

Sincerely,

Michael Forisch
Person Responsible for Regulatory Compliance (PRRC)
support@mysugr.com

Disclaimer: This is a required notice about the mySugr Bolus Calculator in the mySugr app. It is not a marketing or promotional email. This email therefore does not contain an unsubscribe link.

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mySugr GmbH,