Vascutek FSN Ref: FSN2024_01
Dato: 10. okt. 2024
ATT.: Alle implantationer af Thoraflex Hybrid,
hospitalsrisiko
EU-producent SRN: GB-MF-000003643
UDI: diverse

Vascutek Ltd
Newmains Avenue, Inchinnan
Renfrewshire PA4 9RR. Storbritannien

Tlf. +44 (0) 141 812 5555
terumo-aortic.com

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Adrienne Day
Vasutek Ltd

APPENDIX 2 EXTRACT FROM IFU (REF 301-216)

SECTION 8 REFERENCE TO THROMBOSIS AND DEVICE KINKING

8 POTENTIAL ADVERSE EVENTS

Apart from risks associated with (general) open surgeries, and anaesthesia, risks related to the use of Thoraflex Hybrid include, but are not limited to:

- Aneurysm sac or false lumen diameter enlargement
- Aortic rupture
- Consequences of exposure to radiation
- Endoleaks
- Hypersensitivity
- Infection due to device contamination
- Migration (proximal migration of the distal end)
- Patency issues (e.g., stenosis, kinking, thrombosis or incomplete expansion)

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

SECTION 14 DEVICE SIZING

Figure 6 Device sizing

Stent ring oversizing and landing zone guidelines are applicable to all designs. Some movement of the distal ring of the Thoraflex Hybrid implant may occur following re-perfusion of the thoracic aorta. Excessive aortic tortuosity may result in inability to properly position the stent-graft, or stent-graft kinking with thrombus formation. If balloon modelling is desired (e.g., for endoleak, stent-graft kinking or stenosis), use a compliant balloon equal in size to the largest target vessel's diameter. Balloon inflation should not exceed 1 atm.

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SECTION 14.1 DEVICE SIZING AND SELECTION

14 DEVICE SIZING AND SELECTION

14.1 ANEURYSM SIZING (SINGLE STAGE)

This is the recommended sizing for the Thoraflex Hybrid implant for aneurysm treatment, when landed distally in healthy vessel of the descending thoracic aorta. The Thoraflex Hybrid aneurysm sizing chart incorporates a suitable oversize of ring stent diameter to aortic diameter. Aortic diameter is based on inner vessel diameter (ID) measurements therefore no further oversize is required. If outside vessel diameters (OD) are measured, then an allowance for the vessel wall thickness must be made before using the sizing chart for device selection.

Based on testing that has been performed it is recommended that a 40mm distal landing zone length is used and will provide optimum sealing within healthy vessel (Table 1, Figure 4, Figure 5, Figure 6).

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SECTION 14.2

Once a suitable Relay NBS device has been selected to treat D3 (e.g. a 34mm device) then a compatible Thoraflex Hybrid device can be chosen with a relevant D2 (e.g. 32mm)

In these cases where the Thoraflex Hybrid does not create a complete distal seal, using larger devices than required will increase the complexity of sizing the extension device and may additionally increase the risk of thrombus generation until completion of the therapy

Figure 8 Thoraflex Hybrid Extended with a Relay NBS Stent-Graft

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SECTION 15 ANGULATION

The splitter should be positioned in the distal aorta so that when the device is deployed, the collar is in the correct position (Figure 13). For the Thoraflex Hybrid Plexus version, the delivery system should be orientated so that the device branches and aortic arch vessels are aligned.

Note: Excessive aortic tortuosity may result in inability to properly position the stent-graft, or stent-graft kinking with thrombus formation. If balloon modelling is desired (i.e. for endoleak, stent-graft kinking or stenosis), use a compliant balloon equal in size to the largest target vessel's diameter. Balloon inflation should not exceed 1 atm.

Figure 13 Positioning the Thoraflex Hybrid delivery system

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