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✉ support@hedia.com

Urgent Field Safety Notice

Important Update for Hedia Diabetes Assistant Users

We want to inform you about an update to our bolus calculator. Recently, we identified a bug related to unit conversion in the app which required Hedia to restrict access to the app. The bug was internally assessed and a new software update is already in place.

What was the issue, and who was affected?

On February 13th, during a user test, we discovered a conversion error in the latest version of Hedia Diabetes Assistant (version 3.1.2).

This issue affects users who selected mg/dL as their blood glucose unit when setting up the bolus calculator in version 3.1.0 or later versions. The bug caused:

- Incorrect insulin sensitivity values being logged, or
- An inability to complete the calculator setup.

This error happened when the app converted values from mg/dL to mmol/L (since mmol/L is the default unit in the app and server).

If incorrect insulin sensitivity was saved, it could result in higher-than-intended insulin dose recommendations.

What has Hedia done?

To ensure the safety of all users (even if the issue could only potentially affect a small proportion of users) Hedia restricted access to the app.

Hedia deleted the insulin sensitivity for users who could potentially be affected. If your insulin sensitivity has been deleted, it does not mean that you were affected.

A new update has been released to fix the unit conversion issue, so the problem is fixed.

What you need to do now

To ensure everything works correctly:

- ✓ **All users need to update the app** to the latest version (**3.1.3**).
- ✓ **Some users are requested to go through the onboarding flow** to set up the app again.
- ✓ **Double-check your settings** before using the bolus calculator for insulin recommendations—especially the first time after the update and periodically afterwards.

How does Hedia ensure product safety?

At Hedia, ensuring the safety of our products is our top priority. As a medical device company, we adhere to internationally recognized standards, including ISO 13485 for Quality Management and ISO 14971 for Risk Management. Our products undergo regular audits by our Notified Body, TÜV SÜD, to maintain compliance with these strict regulations.

Despite our rigorous testing processes, an issue was identified. We are fully dedicated to learning from this experience and continuously enhancing our quality assurance measures.

We sincerely regret any inconvenience this may have caused and greatly appreciate your trust in Hedia.

Need help?

If you have any questions or experience any issues, our Customer Support team is ready to assist you. Contact us at:

✉ support@hedia.com