

June 17, 2025

## URGENT FIELD SAFETY NOTICE

### Potential for False Positive Results with QuickVue® Dipstick Strep A Test

Dear Valued Customer,

The purpose of this notification is to inform you of the potential for false positive results with QuickVue® Dipstick Strep A Test.

Affected Product Name	Product Code (Unique Device Identifier)	Affected Lot(s)	Expiry
QuickVue® Dipstick Strep A Test, 50T	<b>20108</b> (30014613201083)	709673	05-FEB-2026
		709676	05-FEB-2026
		709484	09-JAN-2026
		709501	09-JAN-2026
		709520	16-JAN-2026
		709535	16-JAN-2026
QuickVue® Dipstick Strep A Test, 25T	<b>20125</b> (30014613201250)	709487	09-JAN-2026
		709503	16-JAN-2026
		709509	16-JAN-2026
		709684	12-FEB-2026
QuickVue® Dipstick Strep A Test, 50T, SCN	<b>20108SC</b> (0130014613330073)	709485	09-JAN-2026
		709502	16-JAN-2026

#### Summary

QuidelOrtho™ received complaints regarding false positive results from customers testing with the QuickVue® Dipstick Strep A Test. False positives may include faint pink lines as well as pink-to-red Test lines.

#### Impact to Results

QuidelOrtho confirmed there is an increase in false positive results, as low as 10% and as high as 40% for certain lots of QuickVue Dipstick Strep A.

Positively biased Strep A results may impact clinical judgement regarding or continued provision of antibiotic therapy in patients with suspected infection, potentially contributing to inappropriate antibiotic use.

The QuickVue Dipstick Strep A test is used in acute patient management; QuidelOrtho does not recommend a review of previous patient results as such review would not be able to detect impacted results without repeat testing of the patients who are unlikely to be ill at

the time of review. Such a review is unlikely to impact patients' management. The results from any diagnostic test should be evaluated in conjunction with a patient's history, risk factors, clinical presentations, signs, and symptoms, and other test results. Please consult your Medical Director regarding any additional actions.

As of 09-JUN-2025, QuidelOrtho has received 31 complaints relating to this issue with no reports of adverse effects.

### REQUIRED ACTION

- Discontinue using, render unusable, and discard your remaining inventory of the above listed lots of QuickVue Dipstick Strep A Test.
- Review the content of this communication with your Medical Director and retain this letter for your laboratory.
- Complete and return the enclosed Confirmation of Receipt form no later than **June 25, 2025**. Upon receipt of your completed Confirmation of Receipt form, QuidelOrtho will replace your discarded inventory, as needed.
- Please forward this notification if the affected product was distributed outside of your facility.

### Contact Information

We apologize for the inconvenience this will cause your laboratory. If you have further questions, please contact Global Services Organization at **1-800-874-1517** or [technicalsupport@quidelortho.com](mailto:technicalsupport@quidelortho.com).

*Enclosures: Confirmation of Receipt form (Ref. CL2025-159a\_EU\_CofR)*

© 2025 QuidelOrtho Corporation. All rights reserved. All trademarks are the property of QuidelOrtho Corporation or its subsidiaries.

Quidel Corporation (Quidel) and Ortho Clinical Diagnostics (Ortho), wholly owned subsidiaries of QuidelOrtho Corporation, are transitioning to our new logo and brand. Due to legal and regulatory requirements for diagnostic products, you may continue to see the names and brands of Quidel and Ortho in addition to QuidelOrtho on our packaging, contracts, and marketing materials