

**URGENT:**  
**MEDICAL DEVICE RECALL – REMOVAL**  
**Recall Number: RA2026-4297377**

## InZone Detachment System

**Attn: Supply Chain Management/ Recall Coordinator/  
Inventory Manager**  
**Date: XXX XX, 2026**



Stryker Neurovascular has initiated a Voluntary Medical Device Recall – Removal on **certain lots** of InZone Detachment System devices. Our records indicate that you have been supplied with the subject devices. We therefore request that you read this notice carefully and complete the actions requested.

**Product Affected** Affected product is listed in Attachment A (InZone Detachment System)

**Product Description** The InZone Detachment System is a sterile, handheld, single-patient-use device designed for use with Stryker Neurovascular detachable embolization coils. The device consists of an enclosure with a detachment button, four LED indicator lamps, and a funnel inset at its distal end. The device comes pre-loaded with two AAAA batteries. Stryker Neurovascular Detachable Coils are sold separately. For the specific intended use and indications for use of each type of Stryker Neurovascular Detachable Coil, consult the device Instructions for Use packaged with each device.

**Product Issue** Stryker Neurovascular has identified an issue affecting certain lots of InZone Detachment System devices. Stryker Neurovascular has observed that some devices may experience premature battery drain. When this occurs, the devices may:

- Not power on;
- Power on with faint audible and visual indicators;
- Be unable to detach a coil as intended.

**Potential Harms** The potential harms associated with the failure to detach the coil include:

- Increased time under anesthesia due to replacement of a defective InZone Detachment System with a backup device.
- Additional measures to safely complete the procedure if a replacement InZone Detachment System is not available, such as removing the coil or using an alternative endovascular or surgical approach.

To date, no patient injuries or deaths have been reported.

**Type of Action** Stryker Neurovascular has investigated and identified the issue which contributed to the premature battery drain noted with the InZone Detachment System devices. A corrective manufacturing step has been implemented. Replacement devices are available.



# Business Reply Form

**FSCA identifier: RA 2026-4297377**

Account number:  
Account name:  
Account Address:

## Product: InZone Detachment System

Please choose the most appropriate box below and fill out the corresponding table(s):

- I no longer have any affected product listed on Attachment A.
- I am returning the following unused devices from Attachment A.
- I further distributed the devices listed and will forward a copy of this notice to the new users. I will also notify Stryker of the devices' new location.

Catalog Number	GTIN	Product Description	Quantity on Hand	Affected Lot Numbers
M00345100950	04546540697950	INZONE DETACHMENT SYSTEM		

*\*Stryker assumes that you are sending all unused devices, and other devices are deemed to be used or thrown away at your account.*

I have received the notification from Stryker Neurovascular stating that they have initiated a product field action for the above referenced product, and I acknowledge receipt and review of this URGENT MEDICAL DEVICE Recall - Removal.

Please return this signed and dated form to **your local Stryker representative**.

**Note:** Your signature indicates that you have received and understand the enclosed notification.

\_\_\_\_\_  
Printed name

\_\_\_\_\_  
Title

\_\_\_\_\_  
Contact phone number

\_\_\_\_\_  
Signature

\_\_\_\_\_  
Date

\_\_\_\_\_  
Email address

\_\_\_\_\_  
Phone Number

## Attachment A Recall –

Catalog number	GTIN	Product description
M00345100950	04546540697950	INZONE DETACHMENT SYSTEM

Material Number	Material Description	Batch	Date of Manufacture	Expiration Date
M00345100950	INZONE DETACHMENT SYSTEM	WMP133638	03-Oct-24	02-Oct-26
M00345100950	INZONE DETACHMENT SYSTEM	WMP133871	18-Oct-24	17-Oct-26
M00345100950	INZONE DETACHMENT SYSTEM	WMP133872	18-Oct-24	17-Oct-26
M00345100950	INZONE DETACHMENT SYSTEM	WMP133988	25-Oct-24	24-Oct-26
M00345100950	INZONE DETACHMENT SYSTEM	WMP133989	25-Oct-24	24-Oct-26
M00345100950	INZONE DETACHMENT SYSTEM	WMP134016	28-Oct-24	27-Oct-26
M00345100950	INZONE DETACHMENT SYSTEM	WMP134132	04-Nov-24	03-Nov-26
M00345100950	INZONE DETACHMENT SYSTEM	WMP134133	04-Nov-24	03-Nov-26
M00345100950	INZONE DETACHMENT SYSTEM	WMP134234	08-Nov-24	07-Nov-26
M00345100950	INZONE DETACHMENT SYSTEM	WMP134235	08-Nov-24	07-Nov-26
M00345100950	INZONE DETACHMENT SYSTEM	WMP134255	11-Nov-24	10-Nov-26
M00345100950	INZONE DETACHMENT SYSTEM	WMP134256	11-Nov-24	10-Nov-26
M00345100950	INZONE DETACHMENT SYSTEM	WMP134380	14-Nov-24	13-Nov-26
M00345100950	INZONE DETACHMENT SYSTEM	WMP134381	14-Nov-24	13-Nov-26
M00345100950	INZONE DETACHMENT SYSTEM	WMP134432	18-Nov-24	17-Nov-26
M00345100950	INZONE DETACHMENT SYSTEM	WMP134433	18-Nov-24	17-Nov-26
M00345100950	INZONE DETACHMENT SYSTEM	WMP134518	21-Nov-24	20-Nov-26
M00345100950	INZONE DETACHMENT SYSTEM	WMP134519	21-Nov-24	20-Nov-26
M00345100950	INZONE DETACHMENT SYSTEM	WMP134609	27-Nov-24	26-Nov-26
M00345100950	INZONE DETACHMENT SYSTEM	WMP134610	27-Nov-24	26-Nov-26
M00345100950	INZONE DETACHMENT SYSTEM	WMP134794	06-Dec-24	05-Dec-26
M00345100950	INZONE DETACHMENT SYSTEM	WMP134795	11-Dec-24	10-Dec-26
M00345100950	INZONE DETACHMENT SYSTEM	WMP134939	12-Dec-24	11-Dec-26
M00345100950	INZONE DETACHMENT SYSTEM	WMP134940	12-Dec-24	11-Dec-26
M00345100950	INZONE DETACHMENT SYSTEM	WMP135187	31-Dec-24	30-Dec-26
M00345100950	INZONE DETACHMENT SYSTEM	WMP135243	02-Jan-25	01-Jan-27
M00345100950	INZONE DETACHMENT SYSTEM	WMP135302	07-Jan-25	06-Jan-27
M00345100950	INZONE DETACHMENT SYSTEM	WMP135314	08-Jan-25	07-Jan-27
M00345100950	INZONE DETACHMENT SYSTEM	WMP135354	10-Jan-25	09-Jan-27
M00345100950	INZONE DETACHMENT SYSTEM	WMP135355	10-Jan-25	09-Jan-27
M00345100950	INZONE DETACHMENT SYSTEM	WMP135387	13-Jan-25	12-Jan-27
M00345100950	INZONE DETACHMENT SYSTEM	WMP135417	14-Jan-25	13-Jan-27
M00345100950	INZONE DETACHMENT SYSTEM	WMP135473	16-Jan-25	15-Jan-27
M00345100950	INZONE DETACHMENT SYSTEM	WMP135475	16-Jan-25	15-Jan-27
M00345100950	INZONE DETACHMENT SYSTEM	WMP135491	17-Jan-25	16-Jan-27

M00345100950	INZONE DETACHMENT SYSTEM	WMP135492	17-Jan-25	16-Jan-27
M00345100950	INZONE DETACHMENT SYSTEM	WMP135623	27-Jan-25	26-Jan-27
M00345100950	INZONE DETACHMENT SYSTEM	WMP135720	30-Jan-25	29-Jan-27
M00345100950	INZONE DETACHMENT SYSTEM	WMP135727	31-Jan-25	30-Jan-27
M00345100950	INZONE DETACHMENT SYSTEM	WMP135807	05-Feb-25	04-Feb-27
M00345100950	INZONE DETACHMENT SYSTEM	WMP135865	07-Feb-25	06-Feb-27
M00345100950	INZONE DETACHMENT SYSTEM	WMP135897	10-Feb-25	09-Feb-27
M00345100950	INZONE DETACHMENT SYSTEM	WMP136084	20-Feb-25	19-Feb-27
M00345100950	INZONE DETACHMENT SYSTEM	WMP136107	24-Feb-25	23-Feb-27
M00345100950	INZONE DETACHMENT SYSTEM	WMP136190	28-Feb-25	27-Feb-27
M00345100950	INZONE DETACHMENT SYSTEM	WMP136200	28-Feb-25	27-Feb-27
M00345100950	INZONE DETACHMENT SYSTEM	WMP136611	24-Mar-25	23-Mar-27
M00345100950	INZONE DETACHMENT SYSTEM	WMP136642	25-Mar-25	24-Mar-27
M00345100950	INZONE DETACHMENT SYSTEM	WMP136736	28-Mar-25	27-Mar-27
M00345100950	INZONE DETACHMENT SYSTEM	WMP136748	28-Mar-25	27-Mar-27
M00345100950	INZONE DETACHMENT SYSTEM	WMP136852	03-Apr-25	02-Apr-27
M00345100950	INZONE DETACHMENT SYSTEM	WMP136894	07-Apr-25	06-Apr-27
M00345100950	INZONE DETACHMENT SYSTEM	WMP136985	11-Apr-25	10-Apr-27