Januar 2019

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Impedances at the time of artifact detection	
Sensor RA Ring > Can	2257 Ω
Sensor RA Tip >> Can	2257 Ω
Sensor RV Ring >> Can	387 Ω
Sensor RV Tip >> Can	430 Ω
A Pace Impedance	> 3000 Ω
RV Pace Impedance	536 Ω

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Brady			
Settings			
Mode	DDDR	Output	
RYTHMIQ™	Off	●A	Trend 3.5 V @ 0.4ms
Lower Rate Limit	60 ppm	■V	Trend 3.5 V @ 0.4ms
Maximum Tracking Rate	130 ppm	Sensitivity	
Maximum Sensor Rate	130 ppm	●A	Fixed 0.75 mV
Paced AV Delay	80 - 180 ms	■V	Fixed 2.5 mV
Sensed AV Delay	65 - 150 ms	Leads	
A-Refractory (PVARP)	240 - 280 ms	●A	
V-Refractory (VRP)	230 - 250 ms	Pace	Bipolar
PVARP after PVC	400 ms	Sense	Bipolar
AV Search +	Off	Safety Switch	On
Blanking		■V	
A-Blank after V-Pace	125 ms	Pace	Bipolar
A-Blank after V-Sense	45 ms	Sense	Bipolar
V-Blank after A-Pace	65 ms	Safety Switch	On
Magnet Response	Pace Async	Rate Adaptive Pacing	
Noise Response	DOO	Minute Ventilation	On
Rate Enhancements		Response Factor	8
Rate Smoothing		Fitness Level	Active
Up	Off %	Ventilatory Threshold	120 ppm
Down	Off %	Ventilatory Thresh. Response	70 %
Sudden Brady Response	Off	Accelerometer	Passive
		Minute Ventilation Sensor Settings	
		Minute Ventilation Sensor	On
		Excitation Current	320 µA
		Vector Selection	Auto Select
		Signal Artifact Monitor	On

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
Model 3120 ZOOM™-programmer

Vælg knappen About (om), og kontroller softwaremodel og version

The screenshot shows the 'ABOUT' window of the ZOOM software. The window title is 'ABOUT' and it has a 'Close' button in the top right corner. The window contains the following information:

Institution: Boston Scientific
Programmer Model: 3120
Programmer Serial Number: 1234567

There are two tabs: 'System Information' and 'Acknowledgements'. The 'System Information' tab is active and displays a table of installed software components:

Component	Model	Version
CONFIDENT	2888	1.08
CONTAK CD/VENTAK CHF	2848	4.1
CONTAK RENEWAL	2845	4.34
CONTAK RENEWAL AVT	2893	3.01
CONTAK RENEWAL TR	2865	2.0
DELTA/VISTA	2881	3.0
DELTA/VISTA Automated Escaping Test	2880	3.01
INGENIO Platform	2869	2.06
LIVIAN	2943	1.02
MINI	2840	6.0
PDF Report Generator	2905	1.03
PRIZM AVT	2849	3.1
PRIZM/PRIZM 2/CONTAK CD 2/VITALITY DS,FL	2844	4.0
Pulsar/Discover/Meridian/CONTAK TR	2890	5.1

On the right side of the window, there are controls for printing: 'Printer' set to 'Internal', 'Number of Copies' set to '1', and a 'Print' button. Below the table, there is an information icon and a note: 'To save the About report, insert a USB drive then press Print.'

At the bottom of the window, there are three buttons: 'Utilities', 'About', and 'Select PG'. The 'About' button is highlighted with a blue border.

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