

**MANDATORY SOFTWARE UPGRADE**  
**ACL ACUSTAR SYSTEM SOFTWARE VERSION 2.2.2**

July 2013

Dear Valued ACL AcuStar Customer:

Instrumentation Laboratory is releasing mandatory System Software V2.2.2 for implementation on the ACL AcuStar. The refined system software is mandatory because it includes an important modification that corrects a software defect present only in System Software V2.2.1.

We internally identified that, in remote instances, System Software V2.2.1 may not cancel or flag certain tests that are affected by a hardware error, creating the potential for an erroneously low test result. System Software V2.2.2 corrects this issue.

Please contact your local ACL AcuStar distributor if you have additional questions regarding this mandatory upgrade.

Sincerely,



Carol Marble  
Regulatory Affairs Director