

Ihre Zeichen

Unsere Zeichen
QMV-kroTel. (+49 74 61) 95-
2182Fax (+49 74 61) 95-
1555Datum
2013-11-22

Safety Information

**Product: TARGON FN TELESCREWS
K0820T, K0822T, K0823T, K0824T, K0825T, K0826T, K0828T**

Internal tests have shown that Targon FN screws can occasionally be difficult to screw into the Targon FN femoral plates. In order to avoid any uncertainty in the use of the product, we would like to inform you of the following facts.

It is an inherent feature of the system that the screw heads of the TeleScrews may protrude when they are screwed into the Targon FN plate. This does not have any functional effects, as the screw depth is only conditionally a criterion for the stability of the screw connection.

Due to the design, the required tightening torque of 8 Nm can be achieved after just half a turn of the screw, depending on the tolerance of the TeleScrew and plate. To ensure a secure connection between the Targon FN plate and the TeleScrews, it is generally recommended to check that the threads of both the Targon FN plate and the TeleScrew are free from any bone particles, blood or tissue before screwing them in. In case of doubt and if you encounter any problems screwing the TeleScrew into the plate, take the TeleScrew out again and clean the thread of the plate and the TeleScrew.

Vorsitzender des Aufsichtsrates:
Prof. Dr. h.c. Ludwig Georg Braun

Vorstand:
Prof. Dr. Hanns-Peter Knaebel
(Vorsitzender)
Dr. Harald Stallforth
(stellv. Vorsitzender)
Dr. Joachim Schulz

Sitz der Gesellschaft: Tuttlingen
Reg. Gericht: Stuttgart HRB 726261
USt. Id.-Nr. DE812160059
WEEE-Reg.-Nr. DE 65109852

Bankverbindungen:
Deutsche Bank AG Tuttlingen
BLZ 653 700 75 Konto 21 22 000 00
IBAN DE44 6537 0075 0212 2000 00
SWIFT / BIC DEUTDE33
Baden-Württembergische Bank
BLZ 600 501 01 Konto 487 1905
IBAN DE31 6005 0101 0004 8719 05
SWIFT / BIC SOLADE33

Hausanschrift:
Aesculap AG
Am Aesculap-Platz
78532 Tuttlingen
Deutschland

According instructions have also been added to the operating manual. Your Aesculap Consultant for Medical Devices will provide you with the updated version as soon as it is available. If you have any questions, please contact one of the following people:

Dr. Andreas Doser

☎ + 49 7461-95 2236

Mr. Robert Weber

☎ + 49 07461-95 2605.

Please make sure that each user of the above mentioned product and any other person that has to be informed within your organization, will take notice of this safety information.

If you have forwarded the involved product to a third party, please provide a copy of this information to this party and also note Aesculap AG about this circumstance.

Please acknowledge the receipt of the safety information to:

**Aesculap AG
Dept. QMV
Ms Kerstin Rothweiler
Am Aesculap-Platz
DE-78532 Tuttlingen
Fax: + 49 7461 95-1555
kerstin.rothweiler@aesculap.de**

Aesculap AG

Quality Management Vigilance

Please complete this form and return to Ms Kerstin Rothweiler, Dept. QMV
Fax: +49 7461-95 1555 or E-mail: kerstin.rothweiler@aesculap.de

Acknowledgement
Safety Information
TARGON FN TELESCREWS
KO820T, KO822T, KO823T, KO824T, KO825T, KO826T, KO828T

1. We acknowledge receipt of the safety information.

2. Please mark accordingly:

User:

We have taken note of the safety information.

Retailer:

We have forwarded the safety information to our customers and attach an according customer list.

Facility _____ Location _____

Name _____

Department _____ Tel. _____

Signature _____ Date _____