

**Cetro**

FIELD SAFETY NOTICE

Medical Device Safety Advisory Notice and Product Recall

Product concerned:	Cetro Scalpel, art. No. 9602
Date:	12 Feb. 2014
Action:	Product recall

Product Advisory Notice:

We have been noticed that our Scalpel with art. no. 9602, used for scalp blood testing, may have a faulty knife blade.

There is a risk the blade may loosen. It concerns all manufactured products with lot no:s starting with the letters CSC. (For example CSC-240).

Product Recall:

To prevent occasional accidents all Cetro Scalpels, mentioned above, is immediately recalled.

You need to take the following actions:

Please return the products, as soon as possible, to:

Cetro AB
Fräsarvägen 30
SE-142 50 Stockholm
Sweden

Your returned products will be replaced by products from another manufacturer, without any costs for you.

We ask you to disperse this safety notice to all people concerned within your organisation, and to send us a mail when you have received and handled this advisory notice.

The undersigned person declares that this Field Safety Notice has been reported to concerned authority.


CETRO AB

Staffan Lavér
Managing Director

Cetro AB
Box 12324
SE-102 28 Stockholm

tele. +46 08-44 82 103
fax -44 82 104

e-mail: cetro@cetro.se



Cetro AB
Box 12324
SE-102 28 Stockholm

tele. +46 08-44 82 103
fax -44 82 104

e-mail: cetro@cetro.se