

30th January 2014

URGENT FIELD SAFETY NOTICE

Oxoid Legionella pneumophila 2-14 Reagent, Lot: 1265366

This reagent is contained within the following products:

Oxoid Legionella Latex Test (DR0800M). Lots: 1314413 and 1301907

Oxoid Legionella pneumophila serogroup 2-14 Test Reagent (DR0802M). Lot: 1265362

Customers are to be advised of the following:

An internal technical investigation has confirmed that Legionella pneumophila 2-14 reagent Lot:1265366, may return false negative results.

Continued use of these lots may result in a delay to the correct serogroup characterisation of Legionella isolates.

Customers who have used these lots should review results reported and consider retesting and/or seek appropriate expert advice.

Accordingly, in keeping with our Quality Policy, we request that you destroy any remaining inventory of the lots listed above and contact Customer Services or your local distributor who will issue a replacement.

This notice needs to be passed on to all who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred.

If you experience any problems complying with this notice or you require any further information, please contact our Technical Support Department on +44 (0)1256 694238, or at oxid.techsupport@thermofisher.com

Yours faithfully,



Colin Booth
Global Director, Quality Assurance and Regulatory Affairs, Microbiology Products