

To the attention of Quality Assurance Dpt or
Regulatory Affairs Dpt or Management

Saint Priest, February 24th 2017

Sujet: **URGENT – FIELD SAFETY NOTICE – RECALL NOTIFICATION LETTER**

Medical device:

PANTA® and PANTA® XL nails

Reference:

**500050 - 500080 - 500150 - 500180 - 500250 - 500280 - 500350 - 500380 - 510111 - 510141 -
510211 - 510241 - 510311 - 510341**

Legal Manufacturer :

*NEWDEAL SAS, Immeuble Séquoia 2 - 97 allée Alexandre Borodine -
Parc Technologique de la Porte des Alpes - 69800 Saint Priest – France.*

Batch involved:

All non-expired and unused products listed on Appendix 1

Madam, Sir,

Newdeal SAS, a company within Integra LifeSciences Group, has recently identified, through an internal evaluation, the possibility of sealing defect for the Panta® or Panta® XL packaging. The defect is a non-homogeneous seal and if it were not completely sealed, the sterility of the packaging or the nail itself could be compromised.

Loss of sterility may result in a wound infection that is significant but reversible, requiring intervention beyond standard-of-care. The package defect might not be easily detectable upon visual inspection prior to use but an adverse health consequence is unlikely to occur based on our health hazard evaluation.

The review of the available clinical data on the Panta® or Panta® XL nails does not raise an abnormal infection rate, consequently no specific follow up for patient implanted is required.

While no adverse event or patient injury has been reported due any package defect, Newdeal SAS has made the decision to conduct a voluntary recall of any unused and unexpired products listed on appendix 1.

We are notifying you of the recall as our records indicate that you have been supplied with some Panta® or Panta® XL nails listed on appendix 1.

We kindly ask you to examine your inventory to determine if you have Panta® or Panta® XL nails listed on appendix 1, please quarantine them.

We also kindly ask you to contact the final customers who may have the affected products and provide them with this letter. If they have affected product, they have to stop using them immediately and remove them from service.

Once the audit of your inventory and your final customers' inventory achieved, please sign and return the "Recall acknowledgment and Return Form" enclosed, by which you confirm that you have received this Recall notification and you intend to fully comply with this Recall notification.

Page 1 of 2

Newdeal

Siège Social : Immeuble Séquoia 2 • 97 allée Alexandre Borodine • Parc Technologique de la Porte des Alpes •
69800 Saint Priest • France

33 (0)4 37 47 51 51 office • 33 (0)4 37 47 51 52 fax • integralife.com

Société par Actions Simplifiée au capital de 1.000.000 € • NAF 4646Z • 412 111 510 RCS Lyon

With this form, you will ensure that all the devices Panta® or Panta® XL nails affected, will be sent back including those already shipped to your customers. You also confirm that this notification has been forwarded to every concerned customer.

Integra Customer Service will contact you upon receipt of this information to organize the return of the concerned products (Return Merchandise Authorization number assignment).

The receipt of this form ensures that Newdeal has achieved a level of effectiveness in communicating this information.

We also recommend that you keep a copy of this notification and a signed copy of the acknowledgement form for your records.

National Competent Authorities may perform audits of field actions of this nature to verify that our customers have been notified and understand the nature of the field action being taken.

The National Competent Authority of your country has been alerted of this Field Safety Corrective Action.

Please feel free to contact me for any additional questions. Your cooperation is appreciated, and we thank you for your continued collaboration.

Yours Sincerely,

Marilyse Latour
NEWDEAL SAS
Quality Assurance and Regulatory Affairs
Manager

In attached file: Recall acknowledgment and Return Form (2 pages)
Appendix 1 – List of references and batches involved (2 pages)

RECALL ACKNOWLEDGMENT AND RETURN FORM

Medical device:

PANTA® and PANTA® XL nail

Référence :

500050 - 500080 - 500150 - 500180 - 500250 - 500280 - 500350 - 500380 - 510111 - 510141 - 510211 - 510241 - 510311 - 510341

Legal Manufacturer :

NEWDEAL SAS, Immeuble Séquoia 2 - 97 allée Alexandre Borodine - Parc Technologique de la Porte des Alpes - 69800 Saint Priest - France.

Batch involved:

**All non-expired and unused products listed on Appendix 1
February 2017**

Please send the form back to:

By fax/telecopy: +33 (0)4 37 47 51 52

Or by e-mail: marilyse.latour@integralife.com

With this form, I confirm that:

I have received, read and understood the information provided in the Integra Recall notification regarding PANTA® and PANTA® XL nails.

I have transferred this recall letter to the persons to whom I have sold and/or place on consignment the concerned products. I ensure that the form is duly returned to me signed by these persons.

My inventory and my final customers' inventory have been reviewed and the results are as follow (*please tick the appropriate answer*):

Yes, I do have affected product(s) in my inventory or my final customers' inventory. These affected product(s) have been isolated and will be sent back.

Please indicate quantity, lot numbers and circle the reference involved in the table below:

No, I do not have the affected product in my inventory.

I ensure that all the affected products, including those I had already sent to my customers are being quarantined and will be shipped back to Integra.

Distributor / Healthcare facility name

Contact Name

Street Address

City, Country, Postal Code

Telephone

Email

Signature

Page 1 of 2

Newdeal

Siège Social : Immeuble Séquoia 2 • 97 allée Alexandre Borodine • Parc Technologique de la Porte des Alpes • 69800 Saint Priest • France

33 (0)4 37 47 51 51 office • 33 (0)4 37 47 51 52 fax • integralife.com

Société par Actions Simplifiée au capital de 1.000.000 € • NAF 4646Z • 412 111 510 RCS Lyon

RECALL ACKNOWLEDGMENT AND RETURN FORM

Medical device:

PANTA® and PANTA® XL nail

Référence :

500050 - 500080 - 500150 - 500180 - 500250 - 500280 - 500350 - 500380 - 510111 - 510141 - 510211 - 510241 - 510311 - 510341

Legal Manufacturer :

NEWDEAL SAS, Immeuble Séquoïa 2 - 97 allée Alexandre Borodine - Parc Technologique de la Porte des Alpes - 69800 Saint Priest - France.

Batch involved:

***All non-expired and unused products listed on Appendix 1
February 2017***

| Reference | Affected batch | Quantity |
|-----------|----------------|----------|
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |

APPENDIX 1 – List of references and batches involved by the recall

| Model # | Batch | Model # | Batch | Model # | Batch |
|---------|--------|---------|-------|---------|--------|
| 500050 | F33A | 500150 | F64Y | 500250 | F500 |
| 500050 | F33B | 500150 | F64Z | 500250 | F5M1 |
| 500050 | F3BN | 500150 | F6L3 | 500250 | F650 |
| 500050 | F3Y5 | 500150 | F7CA | 500250 | F825 |
| 500050 | F3Y6 | 500150 | F7CB | 500250 | F826 |
| 500050 | F4S4 | 500150 | F7CC | 500250 | F827 |
| 500050 | F56Y | 500150 | F7CD | 500250 | F828 |
| 500050 | F5LZ | 500150 | F821 | 500250 | F8T4 |
| 500050 | F64W | 500150 | F822 | 500250 | FGLQ |
| 500050 | F7C8 | 500150 | F8T1 | 500250 | FGSM |
| 500050 | F7C9 | 500150 | F8T2 | 500280 | F27T |
| 500050 | F8SZ | 500150 | F9D6 | 500280 | F33D |
| 500050 | F9PC | 500150 | F9D7 | 500280 | F3BP |
| 500050 | FCUG | 500150 | FEDX | 500280 | F3MA |
| 500050 | FCUH | 500150 | FFB9 | 500280 | F4H6 |
| 500050 | FEDW | 500150 | FG4T | 500280 | F5M2 |
| 500050 | FEVR | 500150 | FGSL | 500280 | F6L4 |
| 500050 | FGAX | 500150 | FHLR | 500280 | F829 |
| 500050 | FGLN | 500180 | F1ZS | 500280 | F82A |
| 500050 | FHLP | 500180 | F33C | 500280 | F82B |
| 500080 | F0F1 | 500180 | F3M9 | 500280 | F82C |
| 500080 | F64X | 500180 | F4H4 | 500280 | F8T5 |
| 500080 | F81Y | 500180 | F4ZZ | 500280 | FFBB |
| 500080 | F81Z | 500180 | F68B | 500350 | E2W1/1 |
| 500080 | FEVS | 500180 | F823 | 500350 | F3BQ |
| 500080 | FGLP | 500180 | F824 | 500350 | F3MB |
| 500080 | FGSK | 500180 | F8T3 | 500350 | F3MD |
| 500150 | F15L | 500180 | FDEC | 500350 | F3Y8 |
| 500150 | F15L/1 | 500180 | FEDY | 500350 | F4H7 |
| 500150 | F1ZQ | 500180 | FFBA | 500350 | F5M3 |
| 500150 | F27U | 500180 | FFQL | 500350 | F651 |
| 500150 | F27V | 500250 | EP14 | 500350 | F6TT |
| 500150 | F3Y7 | 500250 | F0WQ | 500350 | FF3B |
| 500150 | F5M0 | 500250 | F15M | 500380 | F3BR |
| 500150 | F60W | 500250 | F4H5 | 500380 | F3MC |

| Model # | Batch |
|---------|-------|
| 500380 | F4S5 |
| 500380 | F5M4 |
| 500380 | F60X |
| 500380 | F652 |
| 500380 | F82D |
| 500380 | F82E |
| 500380 | F8T6 |
| 500380 | FAF2 |
| 510111 | F3H8 |
| 510111 | F4XA |
| 510111 | F507 |
| 510111 | F5M5 |
| 510111 | F68C |
| 510111 | F82F |
| 510111 | F82G |
| 510111 | F8T7 |
| 510111 | F9DC |
| 510111 | F9DD |

| Model # | Batch |
|---------|-------|
| 510141 | ELSS |
| 510141 | EPFZ |
| 510141 | F3H9 |
| 510141 | F82H |
| 510141 | F82J |
| 510141 | F8LG |
| 510141 | FEVT |
| 510141 | FFBC |
| 510211 | F6JA |
| 510211 | F82K |
| 510211 | F82L |
| 510211 | F8C8 |
| 510211 | F8T8 |
| 510211 | FFQN |
| 510211 | FG4V |
| 510241 | F5M6 |
| 510241 | F82M |
| 510241 | F82N |

| Model # | Batch |
|---------|--------|
| 510241 | F8C9 |
| 510241 | F8T9 |
| 510311 | EPG0/1 |
| 510311 | EPG0/2 |
| 510311 | EPG0/3 |
| 510311 | EPG0/4 |
| 510311 | EPG0/5 |
| 510311 | F82P |
| 510311 | F82Q |
| 510311 | F8TA |
| 510341 | F4H8 |
| 510341 | F508 |
| 510341 | F657 |
| 510341 | F8TB |