

## Update - Urgent Field Safety Notice

<b>Affected Product</b>	<b>ONLY CERTAIN CODES - Note: see Addendum "A" for full details.</b>  <b>Feeding Tubes, Nelaton Catheters, Levins Tubes, Rectal Catheter, Ryles Tubes, Suction Catheters, Stomach Tubes, Medical Packs.</b>
<b>Type of action</b>	<b>Voluntary Recall / Product Disposal</b>
<b>Pennine Healthcare Ref:</b>	<b>PHFSN0219</b>
<b>Product Code(s)</b>	<b>Various – As per Addendum "A"</b>
<b>Lot Number</b>	<b>Various – As per Addendum "A"</b>  <b>An example in this range of LOT numbers is: "06D15"</b>

4th June 2019

Dear Customer,

### **Description of the problem:**


This updated communication is to inform you that alternative actions can be taken in cases where a lack of the devices impacted by this recall may impact patient safety. The alternative actions are as detailed below.

**Note: Where patient safety is not compromised by the lack of availability of this device, the actions detailed in the original FSN dated 8<sup>th</sup> May 2019 must be followed (repeated below).**


The risks posed by this issue remain the same as those detailed in the original FSN.

### **Actions to be taken by the user if devices must be used due to an increase in patient safety risk because alternative non-defective devices are unavailable when needed:**

1. All impacted devices must be inspected immediately prior to the removal from their primary packaging.
2. Using the pictures below for reference, manipulate the distal tip of the device in order to clearly visualize the tip through the film packaging.
3. Identify if side eyes are present.
4. Where side eyes are present, continue to use the devices as normal.
5. We recommend that a second visual check is performed immediately prior to the insertion of the tube.
6. Where no side eyes are present, do not use the device, select an alternative unit and repeat the checks as detailed above (for any rejected units, follow the steps detailed in the original FSN (repeated below).



**Example of correct product with side eyes**



**Example of incorrect product without side eyes**

**Original actions to be taken by the distributor / user:**

1. Identify and quarantine all stock of the affected batch / LOT numbers immediately.
2. If you are a distributor you are required to confirm to Pennine that you have provided this FSN to all your customers and contained any returned stock.
3. Complete the attached Reply Form (Addendum B) to confirm that you have read and understood the contents of this Field Safety Notice.
4. Upon receipt of the completed reply form, Pennine Healthcare will arrange for the return or confirm how destruction of the affected product should be organised. (Proof of destruction and records of type and quantities involved will be needed).
5. For further clarification of quarantined or affected stock replacements or credits, please contact your local distributor or if supplied directly from Pennine our customer services using:

- UK markets – [customer.services@penninehealthcare.co.uk](mailto:customer.services@penninehealthcare.co.uk) tel. +44 (0)1332 794 880
  - International Markets - [exports@penninehealthcare.co.uk](mailto:exports@penninehealthcare.co.uk) tel. +44 (0)1332 794 880
6. Maintain awareness of this Field Safety Notice until all stock of the affected LOT batch has been destroyed or returned to Pennine Healthcare.

**Identification of affected products:**

Use the Product Codes and LOT numbers as detailed in Addendum A to identify the affected devices.

**Transmission of this Field Safety Notice:**

This notice should be passed on to all those who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred.

Please maintain awareness on this notice until all required actions have been performed within your organisation.

We confirm that this field safety corrective action has been communicated to the relevant competent authorities and distributors.

Pennine is committed to providing quality products to our customers and we sincerely apologise for any inconvenience this notice may cause.

For and on behalf of Pennine Healthcare:



Simon Motlib  
Quality Assurance Manager

**Addendum A - Affected Medical Devices – (NB - Only these product codes are affected)**

Product Code	Batch Numbers
Medical Pack FRP-04400	02C16,09C16,10M15,15D16,17B16,
Medical Pack FRP-04700	05K18,06L18,12B18,19M18,29K18,
Medical Pack FRP-04701	04B19,14C19,17D19,23A19,
<b>Feeding Tubes (FT)</b>	
FT-1608/105	01C17,01M17,06C18,07L16,09E18,09L17,11K18,12J17,18J14,
FT-1610/105	09F14,11H17,15J15,17G15,17L14,28A19,30K17,
FT-1908/105	03J14,04H16,04J15,06H15,08A19,09B15,09F14,09G15,10G17,11B19,11G14,11L15,13E15,15D16,16K17,16L18,17D18,18C19,18K18,21L17,22M14,25G17,27E14,27L18,29H17,30G14,
FT-1908/60	01C16,02G15,04A18,04D19,04M18,06K14,12C18,13L18,20F17,20J18,25K18,26D16,
FT-1910/105	03G14,18J14,
FT-1910/60	04M18,10A19,20J18,25K18,
<b>Levins Tubes (LT)</b>	
LT-2208	01C17,01E15,01M17,02B18,04G18,05F15,07H14,07L16,09F17,10G17,11H17,13F16,16B15,17F15,21A15,22H16,22K18,25E18,28C18,28K14,30A15,30D18,31C15,
LT-2210	01B17,01E15,02M16,04G17,05F15,09F17,11H17,12K16,13F18,14D15,15A15,17D18,17F15,17H18,18B16,19C15,23E16,24L17,26E15,26F14,27B18,27C17,28A19,29A16,29G14,30D18,
LT-2212	01B17,01E15,01M17,03E17,03K17,04B16,05F15,06D18,09F14,12M14,15A15,15A18,16B15,16F15,16M16,17B16,17K14,18L14,18M18,19H15,19M16,22A16,22D16,22E15,22F17,22H16,23D18,23E16,25E18,26K15,27G18,29B16,30D18,31C17,18J14,
LT-2214	01B17,01C17,02D19,02M16,03G14,03K16,04B16,04E18,06H15,06L14,06M17,07F18,08F15,09E14,10D17,10L16,11B16,12E15,12M14,13A16,13K17,19E14,19F18,20K14,20K15,21A15,21D16,22H14,22M16,24A18,25F14,25J18,26E15,26J17,26K15,27K16,27L14,29F17,29K18,31E16,
LT-2216	01J16,01L18,02D19,02E14,03E17,03G17,03H17,03K17,04J15,05A16,05B19,05M16,06C17,06C19,06D18,06K15,06M17,07L16,11B16,11D18,11G14,11J14,11K18,12A18,12E14,13F16,14A19,14C18,14D15,15A15,15D16,16F17,16L16,16L18,17B16,17E16,17E18,17L14,19C15,19E14,19J14,19J16,19J17,22C19,22F15,22G14,22G16,22M14,23G15,24B16,26E15,27E14,27F18,27G17,27G18,28B17,28C17,28K15,29B16,30A18,30D18,30K17,
LT-2218	01K18,01M17,02E16,02M16,03K17,04E18,05A16,06C17,06D17,09G15,10B16,10C15,10D19,10G17,11D16,12E14,13F18,13K15,14C18,15A15,15B17,15J17,16F15,16G18,16L18,17E16,17L16,18H14,18L15,19J14,20D15,21A19,22F17,22G14,22M14,22M16,24A18,26E15,26H16,27L14,28H15,28K15,29F16,29K18,30G18,

LT-2220	01J16,02G14,03E18,06A17,08F18,09K17,13G16,17A17,17G17,19B15,27F15,29E15,31C17,
LT-2408	04B16,04H16,04J15,05D16,08A19,11G14,11H16,12J17,18H16,18J14,18L14,22G15,25F14,25J14,26J17,30J16,
LT-2410	03G14,06C18,08A19,08D16,09F14,15K18,16J16,17G15,17L14,18J14,23G15,29A16,30K17,
LT-2412	01M17,03G14,04J15,06C18,06E14,08D16,11H16,14D15,15K18,16G15,16M16,17L14,18G18,19J14,20C18,24J15,26J16,27C17,29B16,30D18,31C17,
LT-2414	02E14,03D17,03G14,03M15,04B16,04H16,05D16,05K15,08J16,09F17,11H16,14C18,14D15,14M17,15C17,15D16,15J15,15K18,16F14,16G15,17H17,17L14,19G18,19J14,22M16,25G17,25J14,27K16,30A15,30E17,
LT-2416	02E14,03G14,04H16,05D16,06C18,06D17,06K15,07L16,08M15,10B15,11A17,11H16,11K18,14A19,14D15,15D16,16F14,16F17,16J15,17H17,18L14,19J14,19J16,20M17,22E17,22G15,22M14,23G18,24B16,26J16,27C17,27G17,29F16,30D18,30K17,
LT-2418	03A18,03H17,03M15,04B16,05F14,08A19,09B15,10C15,10G17,11D16,11H17,12A18,14D15,15D16,15G14,15K18,16J15,17D18,17L14,19H15,19J14,20K16,22G15,22M16,25G17,26C18,28G16,29B16,30E17,30J16,31C17,
LT-2420	02G14,02L17,03B15,10L14,16J16,18J14,
<b>Nelaton Catheters (NC)</b>	
NC-1220/SW	02J16,03F16,03G14,03K14,04A19,04D17,04G18,06B17,07C18,07D15,07L17,08E18,08F18,08L16,09A15,09L15,10D18,10G15,11K17,16D15,16G18,16L15,17E17,18B19,18K18,20J17,20L14,21J16,22H16,22H18,22J14,24B15,24G14,25E16,26D18,28J18,29E15,30K15,
NC-1220/SW/1	04A17,06B17,07D15,10C17,11J18,20L14,24B16,25J15,26K15,
<b>Rectal Catheters (RC)</b>	
RC-1820	02C16,02J16,03F16,04A17,04A19,04D17,05F14,06B17,07C18,07D15,08F18,10D18,10M18,12D16,13M17,15D19,16G18,16M14,17E17,18B19,19K17,20G17,22H18,23B17,24B15,24F14,25H15,26K15,26K16,28J18,29F15,
RC-1820/1	01H14,06E15,07C18,07K15,09A15,09H17,10D18,11A18,11J18,12D16,16L15,17E17,18B19,20G17,22H16,22H18,23B17,26K16,28H14,28J18,29F15,29F16,30L17,
<b>Ryles Tubes (RT)</b>	
RT-2008	01M17,02B18,03M14,04B16,04H16,04J15,05C15,05E15,08A15,08B17,08J16,09E18,09L17,11B19,11H17,12E14,13E15,14A19,15B17,16B15,16L18,17H18,18C19,18G18,18H16,19G17,20D15,21L17,22B18,22K18,23E16,26H16,26J17,27E14,28K14,28K15,31C15,31E16,
RT-2008/BULK	06L15,16C17,21C16,
RT-2008/L	04G18,09F14,11K18,17K14,22H16,22K18,28C18,
RT-2010	01E15,02B18,02M16,03J14,04G17,04J15,05A16,05F15,06C18,06K15,06L14,07B17,07L16,08A15,09C15,09E16,09F14,09F17,10D17,11A17,11H17,13B18,13F18,13M18,14D15,14E15,14J15,15K18,16J16,16L16,16L18,17D18,17G15,17H18,18C16,18C19,18L15,19E14,19G17,20M17,21F16,23B15,24J18,24L17,25C15,25E18,25G17,26H16,27C17,27K16,27L14,27L18,28A19,28K14,29G14,30A15,30G18,31E16,
RT-2010/BULK	10A17,15C16,25J17,27C19,14M15,

RT-2010/L	03L15,04G17,04H15,05F15,06K15,08D16,08K18,09B15,09E16,09F17,13B18,13E15,13F16,13F18,13M18,15D15,15J17,16K14,16L16,17L14,18H16,19G17,19H15,20B17,22F15,22K18,24A18,25E18,26C18,26H16,29A16,31G14,
RT-2012	01L18,01M17,03E17,03F14,04B16,04H15,05C15,06C17,06C18,06D17,06E14,06K15,06L14,07C16,07H14,08A15,08J16,08M15,09E16,09E17,09F14,09F17,10D17,10L16,11D18,11G18,11K17,11K18,11L15,12K16,13E15,13F16,13F18,13G16,13J18,14A19,14K15,15A18,15J17,16G15,16K14,16L16,16M16,17A17,17H17,18D19,18G18,18H16,18L14,18M18,19G17,19H15,20B17,20C18,20D15,20M17,21L17,22B18,22C19,22D16,22F17,22H14,22J14,22K18,22L18,23D18,23K17,24A18,24H18,24J18,24L16,24L17,25E18,26C19,26J16,26K15,27C17,27E15,27G17,27H15,27L14,28A19,29C16,29F15,29G14,30G18,31C15,02K14,
RT-2012/BULK	02C17,06L15,12A18,21C16,
RT-2012/L	05J17,07C16,08A15,08K18,09E16,09F14,11D16,11G18,11H16,13E15,13F16,13F18,13J18,14A19,15D15,15J17,16F15,16G15,16L16,18B16,18D19,18M18,19H15,20B17,21F16,21L17,22L18,23E16,24H18,24J15,24J18,27E14,27G17,30D18,
RT-2014	01D15,01J16,01K18,02B18,02E14,02F14,03A18,03H17,03J18,03K16,04D19,04E18,04G18,04J14,05H15,05M16,06C15,06M17,07F18,07H14,07K15,08A15,08A19,08F15,09G15,09L18,11B19,11D17,12E14,12E17,12F17,12L15,12M14,13B18,13E15,13K17,13L17,14C18,14F16,14K16,14M17,15A15,15B17,15C17,15J15,15J17,15M15,16F14,16K17,16L16,16L18,16M16,17A17,17E18,17H17,18C16,18H14,19G18,19H15,19J14,20D15,20F18,20K14,20K15,21A19,21D16,22K18,23K17,24A18,25F14,25G17,25J18,26H16,26J17,26L15,27B18,27K16,27L18,28K15,29A16,29B16,29F15,29F17,29G14,29K18,29L18,30A15,31C17,31K17,
RT-2014/BULK	03J15,08F17,10D17,11C16,14J17,15F15,17H16,23D15,24B16,27E16,
RT-2014/L	01K18,02D19,02F14,03G14,03H17,03K16,05A16,05D16,05K15,07L18,09C15,09E18,11H16,11J17,13L17,14M17,15A15,15B17,15J17,16M16,17A17,17H18,17L14,19J14,20D15,22K18,25H15,26K15,27B18,27K16,27L14,29A16,29B16,29F17,29L18,
RT-2016	01E14,01E15,01J16,01L18,02F14,03E17,03G14,03H15,03H17,03J18,03K17,03L15,03M14,04G17,05B19,05C15,05M16,06C17,06C18,06C19,06D17,06D18,07H14,07L16,08A15,08A19,08F15,08M15,09G15,09M14,10B16,10C15,11D16,11D18,11G18,11H16,11J17,12A18,12E14,13M18,14B18,16F17,16J15,16K17,16M16,17E18,17F15,18C19,18H14,18H16,18K18,18L14,19H15,19J17,20C18,20D15,20K14,20L18,20M17,21L17,22E15,22G14,22M14,23E17,23G18,24J15,24L16,24L17,25C15,25F14,25H15,25J18,26A17,27F18,27G17,27G18,28B17,28C17,28K14,29F16,30A15,30A18,30K17,31E16,03J14,
RT-2016/L	01E15,02D19,02J14,03J18,04G17,04H16,08A19,09E16,10B15,11G18,11K18,11L15,12A18,12K16,13A16,13F16,13M18,14A19,14B18,14H14,15A15,16K17,17E18,17K14,17L14,19H15,19J17,20K15,21L17,22F17,23G15,24L16,25J18,26A17,27F18,28B17,29C16,15M18,
RT-2018	01C17,01D15,01E15,01F18,01M17,02B18,02M16,03H17,03J18,03L15,03M14,04B16,07B17,07C16,08A19,08E14,08F15,09L17,10B16,10G17,10L16,11B19,11D17,11H17,12A18,12E17,12F14,13F16,13F18,13G18,13K17,13M18,14C18,15M15,16L18,17C16,17D18,17E16,17H18,17L16,18H14,18K18,18L15,21A15,21A19,21D16,21G15,22E15,22H16,22M14,24A18,25C15,25F14,25G17,26C19,26H16,27B18,27H15,27L14,28C18,28K14,29K18,30G18,31C17,02K14,11J14,



RT-2018/L	01E14,03H15,03M14,07B17,08K18,09L17,13E15,17H18,17L16,19H15,20K15,20K16,24A18,26C19,26J17,26K15,29B16,
RT-2020	02B17,03B15,03C16,03F14,05K15,06M18,07C19,09K17,10D18,10J18,12B16,13G18,13K14,14D15,15A16,16A17,16K17,17H15,19J17,20F14,26K16,26L15,27F15,30A18,30D15,30H16,31C15,
RT-2108	01M17,02B18,02D19,02E16,03H15,04G18,04J15,06E14,07C16,08A19,09B15,09C15,09E18,09G15,09L17,10L16,11A17,11B19,11K18,11L15,12F17,12J17,13F16,14A19,15D16,17D18,18G18,18L14,19E14,19G17,20K15,21A15,25E18,25G17,26A17,26J17,27K16,27L18,28C18,29F16,30D18,31E16,
RT-2108/EF	02D19,
RT-2110	02B18,04G17,05A16,05F15,06C18,06L14,07C16,08A15,09B15,09E17,09F14,09F17,09L17,11A17,12A18,12E15,13B18,13F18,13G18,15D15,15D16,15J17,15M15,16F14,16L18,17D18,17F15,18B16,18C19,19E14,19G17,19H15,20K15,20K16,20M17,21L17,22G15,22M14,24L17,25J14,26C15,26C18,26E15,26J17,26L15,27C17,27G17,27K16,27L14,27L18,28A19,28B17,30D18,31E16,
RT-2110/EF	28C19,
RT-2112	03H15,03M15,04J15,05J17,06C15,06D17,06D18,06E14,06K15,06L14,06M17,07F16,07H14,07L18,08A15,09B15,09E16,09E17,09F14,09F17,10L16,11A17,11D16,11G18,11J14,11K17,11L15,12E15,13A16,13G16,13J18,13M18,14K15,15A18,15C17,15K18,16F14,16G15,16J16,16L16,16M15,17A17,18G18,18H14,19F18,19G17,19H15,20B17,20C18,20K16,20M17,21A19,21L17,22A16,22B18,22C19,22D16,22F17,22H16,22L18,22M14,23D18,24A18,24H18,24J15,24J18,25E18,26J16,27C17,27E14,27E15,27H15,28A19,28B17,28L14,29B16,29F15,29G14,29H17,30D15,30D18,31K17,
RT-2114	01B17,02B18,02C17,02E16,03A18,03D17,03J14,03J18,03K17,03M14,03M15,04D19,04E18,04H16,04J15,05B19,05H15,05J17,05K15,06E14,06M17,07C16,07F18,07L16,08F15,09C15,09F17,09L17,09M15,10B16,11D17,11G18,11J17,11L15,12E15,12E17,12M14,13A16,13B18,13F16,13K17,14A19,14D15,14J15,14K16,14M17,15C17,15M14,16F15,16F17,16G15,16J16,16K14,16L18,17E18,17K14,18C16,18C19,18H14,19G18,19J17,20F18,20K15,21F16,22F15,22K18,22M14,23D18,23K17,24A18,24L16,25F14,25G17,25H15,25J14,25J18,26C15,26J17,27B18,27E14,27G18,27L18,28C18,28L14,29F15,29F16,29F17,29G14,29H17,29L18,30A15,31K17,
RT-2116	01C16,01F18,02C17,02F14,02K14,03E17,03J18,03K17,03L15,03M14,04G17,05A16,05B19,06C17,06C19,06D17,06D18,07F16,07H14,07M17,09C15,09M14,10B15,10G15,10L16,11A17,11G18,11J17,11K18,12M14,14A16,14A19,14B18,14J15,15B17,15D15,16F17,16K17,16L18,17D18,17E18,17F15,18H16,19J17,20C14,20C18,20K16,20L18,21F16,21L17,22C19,22F15,22H16,23E17,24A18,25H15,25J18,26C15,26K15,27F18,27G18,27H15,28B17,28L14,29H17,29L18,30A18,30D18,
RT-2118	01B17,01C17,01E14,01K18,01M17,02B18,02E16,03A18,03F14,03G14,03J18,03M15,05F15,05H15,07F16,08K18,09L17,11B19,11K17,12A18,12M14,13A16,13E14,13K15,13M18,14A19,14C18,14J15,15G14,15J17,16B15,16G18,16K14,16L18,17C16,17K14,19G17,21A15,21A19,22B18,22F17,22H16,23K17,24A18,26C15,27K16,27L18,28B17,28C18,28C19,28L14,29F16,
RT-2308	01D15,02E16,06K15,08J16,09G15,09L17,11K18,18H16,18L15,19G17,23E16,26H16,29B16,29H17,30A18,31E16,

RT-2308/EF	08A19,
RT-2310	02E16,03J14,04G17,04H15,06C18,07H14,08A19,08M15,09E18,09L17,09L18,11A17,11H17,13F16,13F18,15D15,16B15,17D18,17H18,19G17,19J16,22E17,22G16,22H14,23E16,24A18,27C17,27G17,28C19,30D15,
RT-2310/EF	28C19,
RT-2312	01B17,03F14,04H15,05J17,07L18,09F17,10G18,11D18,11H16,11K17,13F16,13F18,15A18,15J17,15K18,16L16,18D19,18G18,20M17,22B18,22D16,22F17,22H14,22K18,22L18,23E16,24A18,24H18,25E18,25F14,26C19,27C17,28A19,29H17,
RT-2314	01C17,02B18,04D19,04H16,04J14,05A16,05H15,05M16,07F18,08F15,09E16,09L17,09L18,10B16,11G18,14K16,14M17,15B17,15C17,16F17,16G15,16L16,16M16,17A17,17H18,18H16,19G18,19J17,20C18,20D15,20F18,20K14,20K15,21A15,21A19,21F16,22H14,23B15,23D18,24A18,25G17,26L15,27L14,28C18,28G16,29C16,29F17,29H17,30E17,31C15,31C17,
RT-2314/EF	20F18,27G18,28C18,
RT-2316	01B17,01J16,03H15,03J14,03J18,03K17,04G17,05A16,08A15,09L18,10B15,11D16,11G18,11H16,11J17,13M18,14H18,17B16,17E16,17E18,18K18,19C15,22B18,22E17,22G16,24L16,27C17,27G17,28A19,28B17,30D18,
RT-2316/EF	02D19,06D18,12K17,23G18,
RT-2318	03H15,03K17,10G17,11D17,11L15,12E17,12K16,13K17,14H18,15B17,15D16,16G18,16L18,17D18,17E16,21A19,22F17,28G16,29H17,30E17,31C17,
RT-2318/EF	04G18,26C18,27G18,
RT-2320	19J17,22B17,29H17,29K18,
<b>Suction Catheters (SC)</b>	
SC-1020	01E14,01H14,02C16,03B17,04G18,05F17,06C15,06E16,07C18,07K15,07L14,08E18,09A17,09L15,10D18,10M18,11A18,11H17,11K17,12B15,12D16,12D19,12K16,13D17,14C17,15J16,16L18,17H17,18K18,19G17,20F16,24B16,24G14,25A16,25F14,25H15,26D18,26F17,27L14,28J18,29F15,29L16,30K18,30L15,31G15,31H17,14K14,
<b>Stomach Tubes (ST)</b>	
ST-2408	01E15,02E16,04G18,09M15,10D17,11H17,12F17,16B15,17D18,22K18,
ST-2410	02B18,02K14,02M16,04F14,04G17,05C15,07H14,08D16,09E18,09F17,11H17,12J17,15M15,16G15,17C16,17H18,18B16,19G17,20K15,20M17,22E17,25G17,28A19,28C19,30G18,30K17,
ST-2412	03H15,06C17,06D17,10C15,12M14,14L14,20B17,22A16,23D18,24A18,24L16,25F14,26K15,28A19,31K17,
ST-2414	04J14,05M16,07K15,09E14,10D17,11D16,14C18,15C17,15J15,16M16,19F14,19F18,20D15,21A15,23D18,29A19,31E16,31K17,
ST-2416	06C18,06D17,10B15,11A17,11D16,11G18,11K18,14A19,14D15,16L16,20K14,22C19,22J14,22M14,26C15,29A19,29F16,30K17,
ST-2418	03H17,03M15,04G18,06K15,07F18,09A19,09B15,09F14,10B15,13M18,14C18,18K18,20D15,21A15,21A19,22F17,22H16,22M16,26C19,26J17,27B18,27E14,27G18,29F16,30K17,31C17,



ST-2420	02M16,03A19,03D17,03G18,05C15,05K15,06C18,09K17,10L14,11B15,14D15,14F16,18M15,23F17,25A18,31C15,
<b>X-ray Levins Tubes (XLT)</b>	
XLT-2208	01C17,01E15,02D19,03H15,04G18,05A16,05D16,06L14,07C16,09E18,09M15,10D17,10G17,10L16,11A17,11B19,11G14,11L15,12E14,12J17,14F16,16B15,16K14,18G18,19E14,21L17,22B18,23E16,24L16,25E18,25F14,25G17,25H15,26A17,28C18,28K14,29H17,30G14,31E16,
XLT-2210	02B18,03G14,04G17,05A16,06C18,06D17,07C16,08D16,08K18,09A19,09B15,09C15,09F14,09L18,11A17,12E14,12E17,13B18,13G16,13G18,13M18,14E15,14L14,15K18,16B15,16K14,16L16,17F15,18C19,18J14,18L15,19E14,19F18,19G17,20B17,20K15,22D16,22G16,22K18,23E16,23K17,24A18,24L17,25E18,25H15,26C18,26E15,26H16,26J17,27G17,27K16,27L18,28B17,28K14,29A16,29G14,29H17,30D14,
XLT-2212	01B17,03K17,04J15,05F15,05J17,06D18,06E14,06L14,07C16,07H14,07L18,08K18,09B15,09M15,10D17,11D16,12J18,12M14,13A16,13F18,13G16,14D15,14E15,15A18,16B15,16G15,16L16,16M15,16M16,17K14,18G18,18J14,18M18,19G17,19H15,20B17,20C18,20K16,20M17,22B18,22C19,22F17,22H16,22K18,23E16,23K17,24H18,24J18,25E18,26K15,27E14,27H15,29G14,31C17,31E16,
XLT-2212/BULK	11L14,11M14,19E14,
XLT-2214	01C17,01F18,02B18,02D19,02E14,02M16,03A18,03H15,03J18,03K16,04E18,05B19,05J17,06H15,06M17,07F18,07K15,08A15,08A19,09E17,09E18,09F14,10B16,10D17,11G18,11H16,12F17,12L15,13A16,13K17,14K16,14M17,15B17,15C17,15D16,15J15,15K18,15M15,16G15,16L18,17A17,17E16,17L14,18C19,18H14,18L15,18M18,19E14,19F18,19J17,20C18,20K14,21D16,21F16,22G16,23B15,23K17,24A18,24J15,25G17,25J14,26E15,26H16,26J17,26K15,26L15,27G18,27L14,27L18,29A16,29B16,29F16,29F17,29G14,29L18,30A15,30D18,31C15,
XLT-2216	01B17,01F18,01J16,01L18,02F14,03E17,03G14,03J14,03K17,04G17,04H15,04J15,05B19,05M16,06C17,06C19,07C16,07F16,07H14,07M17,08A19,08F15,09C15,09G15,10B16,11D16,11G14,11G18,11H16,11J14,11J17,11K18,11L15,12A18,12K16,12M14,13A16,13G16,13K15,13M18,14B18,14D15,14E15,14H18,15A15,15M15,16F14,16F17,16J15,16K17,16L16,16L18,17E16,17E18,18K18,18L14,19E14,19J17,20C18,20D15,20K14,20M17,21F16,22B18,22C19,22F15,22G14,22M14,23B15,23G18,24L16,24L17,25J18,26A17,27F18,27G17,27L14,28C17,28K15,29H17,30A15,30A18,30D18,30E17,30K17,31E16,
XLT-2218	01B17,01K18,01M17,02K14,02M16,03A18,03J18,04E18,05A16,05B19,06H15,06L14,08A15,08A19,08J16,08K18,09G15,10G17,11D16,11D17,11H16,11L15,12E15,12K16,13B18,13F16,13F18,13K15,13K17,14C18,15J17,15M15,16F14,16G18,16J15,17E16,17H18,18K18,18L15,20K14,20K16,21A19,21D16,22A16,22F17,22G14,22L18,22M16,23B15,23E16,23K17,25E18,25G17,25J18,26C18,26E15,26H16,26J17,26L15,28C17,28G16,28K15,29B16,29H17,29K18,30J16,31C15,
XLT-2220	01B17,03B15,03E18,05C15,05K15,06L14,09K17,16L15,19A15,20F14,23B16,23F17,30D15,31C15,
XLT-2408	03A18,03J14,04G18,07L16,09B15,11G14,14C18,18L14,26A17,27L18,28K14,28K15,
XLT-2410	09F17,13G16,16B15,24J18,25H15,29F16,29G14,
XLT-2412	01L18,08J16,09F17,15D16,18H16,18J14,24J18,26C19,27C17,27H15,03G14,
XLT-2414	01C17,03D17,04E18,06H15,07L18,08A19,09E17,09M15,11B16,11G18,11H16,15K18,16J16,21A19,21D16,23G14,27G18,27L18,29C16,29H17,29K18,30A15,31K17,

XLT-2416	01B17,01F18,02D19,02E14,02E16,03E17,06C15,06C19,06D18,08A19,09B15,09M15,10B15,10B16,11G14,11L15,12E14,13K15,13M18,14A19,16F17,17E18,18K18,19C15,20C18,20D15,20L18,21A15,22G14,23E17,23G15,23G18,24B16,24L17,25F14,27F18,27L14,28C17,28C19,28K14,29C16,30A18,
XLT-2418	01D15,01E14,02C17,03M15,07F16,09A19,09B15,09C15,09E17,10B16,12E17,18K18,19C15,19H15,19J14,22F15,22F17,22G14,22G15,22H14,23B16,26C18,26C19,27G18,29C16,29H17,30K17,07H14,
XLT-2420	03A19,03E16,04G18,07C19,08G15,10L14,12B16,12E17,13K14,17K14,18K18,19H16,19J17,19L14,23G14,25A18,26L14,26L15,29L17,30A18,31C15,
<b>X-Ray Stomach Tubes (XST)</b>	
XST-2408	01D15,01E15,06K15,08M15,09B15,12F17,18B16,18H14,19G17,22K18,25J14,
XST-2410	19G17,26L15,29A16,
XST-2412	26L15,27L14,
XST-2414	01C17,07K15,14K16,29F17,
XST-2416	04G17,06D17,13K15,18H16,
XST-2418	02C17,06K15,19G17,22H16,
XST-2418/BULK	01B18,10G17,16H18,19K15,20M16,21G15,25H17,26K18,29G16,29H18,
XST-2420	05C15,

## Addendum B

### Updated Field Safety Notice Acknowledgement Form

Name	
Position	
Company/Institute	
Address	

--	--

We do not have any stock of the affected product. (✓ tick if this is your answer)	
---	--

Product Code and Quantity of Affected Product in Stock		
Product Code	LOT Number(s)	Qty.

NB: Please contact your local distributor representative or Pennine directly to confirm quantity and type of quarantined stock before completing returns or plans for destruction. **Evidence will be required prior to and after destruction.**

I hereby confirm that I have read and understood the contents of this Updated Field Safety Notice.	
Sign	
Date	

Please return the completed reply form to: [recalls@penninehealthcare.co.uk](mailto:recalls@penninehealthcare.co.uk)