

August 2022

**URGENT – FIELD SAFETY NOTICE**

Type of Action	Recall Notice
Teleflex Reference	EIF-000513
Product Codes	Commercial Name
19012T	Iso-Gard® Filter S with Expandi-Flex/Elbow, Tethered Cap
191667-000100	Filter + Catheter Mount
19211; 19212	Iso-Gard® Filter S
19261; 19272	Iso-Gard® Filter S with Expandi-Flex
19261T; 19262T; 19272T	Iso-Gard® Filter S with Expandi-Flex, Tethered Cap
19211T; 19212T	Iso-Gard® Filter S, Tethered Cap

Dear Customer,

**Details of impacted devices**

Teleflex Medical Europe Limited, a subsidiary of Teleflex Incorporated, has initiated a voluntary Field Safety Corrective Action (“FSCA”) for the above listed products; refer to Appendix 2 for Unique Device Identifier (“UDI”) and lot number information

**Description of the issue & immediate actions required**

Teleflex is initiating this voluntary recall for the above-mentioned products due to reports of the device splitting or detaching during use. The devices are intended for connection to respiratory equipment in intensive care units and operating theatres to protect the equipment from potential airborne contaminants. If the filter splits or detaches during use, the functionality and efficiency of the device may be compromised, potentially resulting in a leakage and insufficient air supply to the patient with potential desaturation.

As of August 2022, a total of 36 complaints reporting the filter splitting or detaching have been received for the products in scope of this Field Safety Notice. Of these 36 complaints, 4 complaints reported injuries including patient desaturation. No deaths or long-term patient injuries have been reported at this time.

Our records indicate you have received products that are in scope of this Field Safety Notification.

Depending on your location please adhere to the following Action list:

Device location	Action List Number
Medical facilities (hospitals, medical staff, etc.)	<b>1</b>
Distributors	<b>2</b>

### Action list number 1 – Medical facilities

1. We request that you immediately check your inventory for product within the scope of this FSCA. Users should cease use and distribution of affected product and immediately quarantine the affected product.
2. If you have impacted product, mark the applicable checkbox on the Acknowledgement Form (Appendix 1) and contact Teleflex Customer Service by calling the phone number provided below. Teleflex Customer Service will issue a return goods authorisation (RGA) number to you. Write the (RGA) number into the respective field in the Acknowledgement Form and promptly return this form by email to Teleflex Customer Service.
3. If you do not have impacted product, mark the applicable checkbox on the Acknowledgement Form (Appendix 1) and return the form to Teleflex at the contact details provided.
4. Teleflex (or your local dealer) will issue a credit note upon receipt of the returned affected product.

### Action list number 2 – Distributors

1. Provide this field safety notice to all customers who have received impacted product. Each of your customers is then required to complete the Acknowledgement Form and return it to you.
2. We request that you immediately check your inventory for impacted product. Cease use and distribution of impacted product and immediately quarantine the affected product. You may then return all product in scope.
3. As a distributor, you are then required to confirm to Teleflex that you have completed the field activity outlined in actions 1 and 2 of this Action List Number 2. Upon completion of your actions, please forward the completed Acknowledgement Form to Teleflex Customer Service.  
**Important** - Please ensure you only list batch numbers in scope of this Field Safety Notice when completing this form.
4. Please be aware that all European Economic Area/Switzerland, United Kingdom (EEA/CH/UK) and Turkey (TR) Competent Authorities in which Teleflex distribute directly will be notified by Teleflex.
5. If you have further distributed product outside of your country, please notify Teleflex Customer Service by return email to the email address below.
6. If you are a distributor and/or have a reporting responsibility within or outside the EEA/CH/UK/TR area, please notify your local Competent Authority of this action. Please forward the notification and all communication with your local competent authority to Teleflex.

Adverse reactions or quality problems experienced with the use of this product should be reported to Teleflex Customer Service at the contact information below.

### Teleflex

Teleflex informs impacted customers, employees of Teleflex and distributors of this Field Safety Corrective Action.

### Transmission of this Field Safety Notice

This notice should be passed on to all persons who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred.

Please consider end users, clinicians, risk managers, supply chain/distribution centres, etc., in the circulation of this notice. Please maintain awareness of this notice until all required actions have been completed in your organisation.

**Contact reference person**

Should you require any further information or support concerning this issue, please contact:

**Customer Service:**

**Contact:** Shane Kenny

**Telephone:** +353 (0)86 3479154

**Email:** [Recalls.Intl@teleflex.com](mailto:Recalls.Intl@teleflex.com)

Teleflex is committed to providing high quality, safe and effective products. We regret any inconvenience this action may cause your operations. If you have any other questions, please contact your local Teleflex sales representative or Teleflex Customer Service.

*For and on behalf of Teleflex,*

*Padraig Hegarty*

**Padraig Hegarty VP, Global QA (Manufacturing)**

Appendix 1

Customer No  
\_\_\_\_\_

**FIELD SAFETY CORRECTIVE ACTION**  
**ACKNOWLEDGEMENT FORM**

**PRODUCT FIELD ACTION BY TELEFLEX – IMMEDIATE ATTENTION REQUIRED**

Ref. EIF-000513

**RETURN COMPLETED FORM IMMEDIATELY TO:**

Email: [Recalls.Intl@teleflex.com](mailto:Recalls.Intl@teleflex.com)

<input type="checkbox"/> We confirm receipt of this FSN and completed the required actions contained therein. We confirm that our inventory does <b>NOT</b> include products impacted by this Field Action.	<input type="checkbox"/> We confirm receipt of this FSN and completed the required actions contained therein. We confirm our inventory <b>DOES</b> include products impacted by this Field Action. The use and further distribution of the affected products is stopped. All products are put on hold and the amount below will be returned. <b>Return Goods Authorisation No</b> _____
---	--

**Complete this Acknowledgement form and return immediately by using the contact information above.**

Product code	Lot/batch number	Quantity returning
<b>Important</b> - Please ensure you only list batch numbers in scope of this Field Safety Notice when completing this form.		
<ul style="list-style-type: none"> <li>Include a copy of the <b>completed Acknowledgement Form</b> in the returns package with the returned units</li> <li>Ensure the <b>RGA number is clearly visible</b> on the returns package</li> <li>Please label returns as <b>"Field Safety Returns"</b></li> </ul> <p><b>Note:</b> Non FSCA product returns should be processed per standard product return processes.</p>		

INSTITUTION NAME (EG NAME OF HOSPITAL, HEALTH CARE ORGANISATION)	
INSTITUTION ADDRESS	Phone/FAX
FORM COMPLETED BY	Stamp
PRINT NAME: _____  SIGNATURE: _____	
DATE	

**Appendix 2 - Product Code / Lot Number / Unique Device Identifier (UDI)**

Product Code	Batch	UDI Number
19211	19LT34	(01)04026704347995(17)241128(10)19LT34
	KMH20K0207	(01)04026704347995(17)250928(10)KMH20K0207
	KMH20L0238	(01)04026704347995(17)251028(10)KMH20L0238
	KMH20L0243	(01)04026704347995(17)251028(10)KMH20L0243
	KMH20M0093	(01)04026704347995(17)251128(10)KMH20M0093
	KMH21A0321	(01)04026704347995(17)251228(10)KMH21A0321
	KMH21C0100	(01)04026704347995(17)260228(10)KMH21C0100
	KMH21C0109	(01)04026704347995(17)260228(10)KMH21C0109
	KMH21C0118	(01)04026704347995(17)260228(10)KMH21C0118
	KMH21C0292	(01)04026704347995(17)260228(10)KMH21C0292
	KMH21C0357	(01)04026704347995(17)260228(10)KMH21C0357
	KMH21D0323	(01)04026704347995(17)260328(10)KMH21D0323
	KMH21D0354	(01)04026704347995(17)260328(10)KMH21D0354
	KMH21E0118	(01)04026704347995(17)260428(10)KMH21E0118
	KMH21F0369	(01)04026704347995(17)260528(10)KMH21F0369
	KMH21G0338	(01)04026704347995(17)260628(10)KMH21G0338
	20FT22	(01)04026704347995(17)250528(10)20FT22
	20FT20	(01)04026704347995(17)250528(10)20FT20
	KMH21G0388	(01)04026704347995(17)260628(10)KMH21G0388
	KMH21H0027	(01)04026704347995(17)260728(10)KMH21H0027
	KMH21H0077	(01)04026704347995(17)260728(10)KMH21H0077
	20GT15	(01)04026704347995(17)250628(10)20GT15
	KMH21H0188	(01)04026704347995(17)260728(10)KMH21H0188
	KMH21H0214	(01)04026704347995(17)260728(10)KMH21H0214
	KMH21J0328	(01)04026704347995(17)260828(10)KMH21J0328
	KMH21K0029	(01)04026704347995(17)260928(10)KMH21K0029
	KMH20K0194	(01)04026704347995(17)250928(10)KMH20K0194
	KMH20K0238	(01)04026704347995(17)250928(10)KMH20K0238
	KMH21K0083	(01)04026704347995(17)260928(10)KMH21K0083
	KMH21K0401	(01)04026704347995(17)260928(10)KMH21K0401
	KMH20K0253	(01)04026704347995(17)250928(10)KMH20K0253
	KMH20K0259	(01)04026704347995(17)250928(10)KMH20K0259
	KMH20K0281	(01)04026704347995(17)250928(10)KMH20K0281
	KMH20K0273	(01)04026704347995(17)250928(10)KMH20K0273
	KMH20K0283	(01)04026704347995(17)250928(10)KMH20K0283
	KMH20K0409	(01)04026704347995(17)250928(10)KMH20K0409
	KMH20K0337	(01)04026704347995(17)250928(10)KMH20K0337
	KMH20K0414	(01)04026704347995(17)250928(10)KMH20K0414
	KMH20K0415	(01)04026704347995(17)250928(10)KMH20K0415
	KMH20K0423	(01)04026704347995(17)250928(10)KMH20K0423
	KMH20K0430	(01)04026704347995(17)250928(10)KMH20K0430
	KMH20K0431	(01)04026704347995(17)250928(10)KMH20K0431
	KMH20K0434	(01)04026704347995(17)250928(10)KMH20K0434
	KMH20K0438	(01)04026704347995(17)250928(10)KMH20K0438
	KMH20L0200	(01)04026704347995(17)251028(10)KMH20L0200
	KMH20L0204	(01)04026704347995(17)251028(10)KMH20L0204
	KMH20L0239	(01)04026704347995(17)251028(10)KMH20L0239
	KMH20L0240	(01)04026704347995(17)251028(10)KMH20L0240
	KMH20L0241	(01)04026704347995(17)251028(10)KMH20L0241
	KMH20L0242	(01)04026704347995(17)251028(10)KMH20L0242
KMH20L0252	(01)04026704347995(17)251028(10)KMH20L0252	
KMH20L0302	(01)04026704347995(17)251028(10)KMH20L0302	
KMH20L0296	(01)04026704347995(17)251028(10)KMH20L0296	
KMH20L0322	(01)04026704347995(17)251028(10)KMH20L0322	
KMH20L0324	(01)04026704347995(17)251028(10)KMH20L0324	
KMH20L0421	(01)04026704347995(17)251028(10)KMH20L0421	
KMH20M0098	(01)04026704347995(17)251128(10)KMH20M0098	
KMH20M0103	(01)04026704347995(17)251128(10)KMH20M0103	

**Appendix 2 - Product Code / Lot Number / Unique Device Identifier (UDI)**

Product Code	Batch	UDI Number
19211	KMH20M0244	(01)04026704347995(17)251128(10)KMH20M0244
	KMH20M0247	(01)04026704347995(17)251128(10)KMH20M0247
	KMH20M0267	(01)04026704347995(17)251128(10)KMH20M0267
	KMH20M0269	(01)04026704347995(17)251128(10)KMH20M0269
	KMH20M0292	(01)04026704347995(17)251128(10)KMH20M0292
	KMH20M0301	(01)04026704347995(17)251128(10)KMH20M0301
	KMH20M0316	(01)04026704347995(17)251128(10)KMH20M0316
	KMH20M0461	(01)04026704347995(17)251128(10)KMH20M0461
	KMH21A0330	(01)04026704347995(17)251228(10)KMH21A0330
	KMH21A0349	(01)04026704347995(17)251228(10)KMH21A0349
	KMH21A0365	(01)04026704347995(17)251228(10)KMH21A0365
	KMH21A0379	(01)04026704347995(17)251228(10)KMH21A0379
	KMH21A0407	(01)04026704347995(17)251228(10)KMH21A0407
	KMH21A0430	(01)04026704347995(17)251228(10)KMH21A0430
	KMH21A0434	(01)04026704347995(17)251228(10)KMH21A0434
	KMH21C0090	(01)04026704347995(17)260228(10)KMH21C0090
	KMH21C0101	(01)04026704347995(17)260228(10)KMH21C0101
	KMH21C0108	(01)04026704347995(17)260228(10)KMH21C0108
	KMH21C0105	(01)04026704347995(17)260228(10)KMH21C0105
	KMH21C0114	(01)04026704347995(17)260228(10)KMH21C0114
	KMH21C0112	(01)04026704347995(17)260228(10)KMH21C0112
	KMH21C0122	(01)04026704347995(17)260228(10)KMH21C0122
	KMH21C0124	(01)04026704347995(17)260228(10)KMH21C0124
	KMH21C0125	(01)04026704347995(17)260228(10)KMH21C0125
	KMH21C0303	(01)04026704347995(17)260228(10)KMH21C0303
	KMH21C0282	(01)04026704347995(17)260228(10)KMH21C0282
	KMH21C0286	(01)04026704347995(17)260228(10)KMH21C0286
	KMH21C0306	(01)04026704347995(17)260228(10)KMH21C0306
	KMH21C0296	(01)04026704347995(17)260228(10)KMH21C0296
	KMH21C0315	(01)04026704347995(17)260228(10)KMH21C0315
	KMH21C0318	(01)04026704347995(17)260228(10)KMH21C0318
	KMH21C0337	(01)04026704347995(17)260228(10)KMH21C0337
	KMH21C0320	(01)04026704347995(17)260228(10)KMH21C0320
	KMH21C0341	(01)04026704347995(17)260228(10)KMH21C0341
	KMH21C0348	(01)04026704347995(17)260228(10)KMH21C0348
	KMH21C0349	(01)04026704347995(17)260228(10)KMH21C0349
	KMH21C0373	(01)04026704347995(17)260228(10)KMH21C0373
	KMH21C0372	(01)04026704347995(17)260228(10)KMH21C0372
	KMH21C0375	(01)04026704347995(17)260228(10)KMH21C0375
	KMH21C0384	(01)04026704347995(17)260228(10)KMH21C0384
	KMH21C0395	(01)04026704347995(17)260228(10)KMH21C0395
	KMH21D0170	(01)04026704347995(17)260328(10)KMH21D0170
	KMH21D0305	(01)04026704347995(17)260328(10)KMH21D0305
	KMH21D0325	(01)04026704347995(17)260328(10)KMH21D0325
	KMH21D0309	(01)04026704347995(17)260328(10)KMH21D0309
	KMH21D0326	(01)04026704347995(17)260328(10)KMH21D0326
	KMH21D0333	(01)04026704347995(17)260328(10)KMH21D0333
	KMH21D0357	(01)04026704347995(17)260328(10)KMH21D0357
	KMH21D0361	(01)04026704347995(17)260328(10)KMH21D0361
	KMH21D0363	(01)04026704347995(17)260328(10)KMH21D0363
KMH21E0001	(01)04026704347995(17)260428(10)KMH21E0001	
KMH21E0006	(01)04026704347995(17)260428(10)KMH21E0006	
KMH21E0162	(01)04026704347995(17)260428(10)KMH21E0162	
KMH21E0193	(01)04026704347995(17)260428(10)KMH21E0193	
KMH21E0268	(01)04026704347995(17)260428(10)KMH21E0268	
KMH21F0352	(01)04026704347995(17)260528(10)KMH21F0352	
KMH21F0376	(01)04026704347995(17)260528(10)KMH21F0376	
KMH21F0371	(01)04026704347995(17)260528(10)KMH21F0371	

**Appendix 2 - Product Code / Lot Number / Unique Device Identifier (UDI)**

Product Code	Batch	UDI Number
19211	KMH21F0378	(01)04026704347995(17)260528(10)KMH21F0378
	KMH21F0381	(01)04026704347995(17)260528(10)KMH21F0381
	KMH21F0396	(01)04026704347995(17)260528(10)KMH21F0396
	KMH21F0387	(01)04026704347995(17)260528(10)KMH21F0387
	KMH21F0389	(01)04026704347995(17)260528(10)KMH21F0389
	KMH21F0401	(01)04026704347995(17)260528(10)KMH21F0401
	KMH21F0404	(01)04026704347995(17)260528(10)KMH21F0404
	KMH21G0291	(01)04026704347995(17)260628(10)KMH21G0291
	KMH21F0407	(01)04026704347995(17)260528(10)KMH21F0407
	KMH21F0410	(01)04026704347995(17)260528(10)KMH21F0410
	KMH21G0016	(01)04026704347995(17)260628(10)KMH21G0016
	KMH21G0009	(01)04026704347995(17)260628(10)KMH21G0009
	KMH21G0012	(01)04026704347995(17)260628(10)KMH21G0012
	KMH21G0021	(01)04026704347995(17)260628(10)KMH21G0021
	KMH21G0026	(01)04026704347995(17)260628(10)KMH21G0026
	KMH21G0025	(01)04026704347995(17)260628(10)KMH21G0025
	KMH21G0035	(01)04026704347995(17)260628(10)KMH21G0035
	KMH21G0030	(01)04026704347995(17)260628(10)KMH21G0030
	KMH21G0036	(01)04026704347995(17)260628(10)KMH21G0036
	KMH21G0051	(01)04026704347995(17)260628(10)KMH21G0051
	KMH21G0042	(01)04026704347995(17)260628(10)KMH21G0042
	KMH21G0046	(01)04026704347995(17)260628(10)KMH21G0046
	KMH21G0050	(01)04026704347995(17)260628(10)KMH21G0050
	KMH21G0052	(01)04026704347995(17)260628(10)KMH21G0052
	KMH21G0054	(01)04026704347995(17)260628(10)KMH21G0054
	KMH21G0060	(01)04026704347995(17)260628(10)KMH21G0060
	KMH21G0055	(01)04026704347995(17)260628(10)KMH21G0055
	KMH21G0056	(01)04026704347995(17)260628(10)KMH21G0056
	KMH21G0082	(01)04026704347995(17)260628(10)KMH21G0082
	KMH21G0064	(01)04026704347995(17)260628(10)KMH21G0064
	KMH21G0070	(01)04026704347995(17)260628(10)KMH21G0070
	KMH21G0074	(01)04026704347995(17)260628(10)KMH21G0074
	KMH21G0086	(01)04026704347995(17)260628(10)KMH21G0086
	KMH21G0088	(01)04026704347995(17)260628(10)KMH21G0088
	KMH21G0091	(01)04026704347995(17)260628(10)KMH21G0091
	KMH21G0098	(01)04026704347995(17)260628(10)KMH21G0098
	KMH21G0093	(01)04026704347995(17)260628(10)KMH21G0093
	KMH21G0104	(01)04026704347995(17)260628(10)KMH21G0104
	KMH21G0114	(01)04026704347995(17)260628(10)KMH21G0114
	KMH21G0119	(01)04026704347995(17)260628(10)KMH21G0119
	KMH21G0172	(01)04026704347995(17)260628(10)KMH21G0172
	KMH21G0158	(01)04026704347995(17)260628(10)KMH21G0158
	KMH21G0159	(01)04026704347995(17)260628(10)KMH21G0159
	KMH21G0184	(01)04026704347995(17)260628(10)KMH21G0184
	KMH21G0186	(01)04026704347995(17)260628(10)KMH21G0186
	KMH21G0279	(01)04026704347995(17)260628(10)KMH21G0279
	KMH21G0271	(01)04026704347995(17)260628(10)KMH21G0271
	KMH21G0282	(01)04026704347995(17)260628(10)KMH21G0282
	KMH21G0290	(01)04026704347995(17)260628(10)KMH21G0290
	KMH21G0300	(01)04026704347995(17)260628(10)KMH21G0300
KMH21G0296	(01)04026704347995(17)260628(10)KMH21G0296	
KMH21G0307	(01)04026704347995(17)260628(10)KMH21G0307	
KMH21G0319	(01)04026704347995(17)260628(10)KMH21G0319	
KMH21G0314	(01)04026704347995(17)260628(10)KMH21G0314	
KMH21G0320	(01)04026704347995(17)260628(10)KMH21G0320	
KMH21G0335	(01)04026704347995(17)260628(10)KMH21G0335	
KMH21G0323	(01)04026704347995(17)260628(10)KMH21G0323	
KMH21G0332	(01)04026704347995(17)260628(10)KMH21G0332	

**Appendix 2 - Product Code / Lot Number / Unique Device Identifier (UDI)**

Product Code	Batch	UDI Number
19211	KMH21G0345	(01)04026704347995(17)260628(10)KMH21G0345
	KMH21G0356	(01)04026704347995(17)260628(10)KMH21G0356
	KMH21G0358	(01)04026704347995(17)260628(10)KMH21G0358
	KMH21G0361	(01)04026704347995(17)260628(10)KMH21G0361
	KMH21G0359	(01)04026704347995(17)260628(10)KMH21G0359
	KMH21G0365	(01)04026704347995(17)260628(10)KMH21G0365
	KMH21G0369	(01)04026704347995(17)260628(10)KMH21G0369
	KMH21G0368	(01)04026704347995(17)260628(10)KMH21G0368
	KMH21G0371	(01)04026704347995(17)260628(10)KMH21G0371
	KMH21G0373	(01)04026704347995(17)260628(10)KMH21G0373
	KMH21G0376	(01)04026704347995(17)260628(10)KMH21G0376
	KMH21G0378	(01)04026704347995(17)260628(10)KMH21G0378
	KMH21G0381	(01)04026704347995(17)260628(10)KMH21G0381
	KMH21G0390	(01)04026704347995(17)260628(10)KMH21G0390
	KMH21G0391	(01)04026704347995(17)260628(10)KMH21G0391
	KMH21G0425	(01)04026704347995(17)260628(10)KMH21G0425
	KMH21H0002	(01)04026704347995(17)260728(10)KMH21H0002
	KMH21G0438	(01)04026704347995(17)260628(10)KMH21G0438
	KMH21H0014	(01)04026704347995(17)260728(10)KMH21H0014
	KMH21H0008	(01)04026704347995(17)260728(10)KMH21H0008
	KMH21H0019	(01)04026704347995(17)260728(10)KMH21H0019
	KMH21H0016	(01)04026704347995(17)260728(10)KMH21H0016
	KMH21H0020	(01)04026704347995(17)260728(10)KMH21H0020
	KMH21H0021	(01)04026704347995(17)260728(10)KMH21H0021
	KMH21H0029	(01)04026704347995(17)260728(10)KMH21H0029
	KMH21H0031	(01)04026704347995(17)260728(10)KMH21H0031
	KMH21H0033	(01)04026704347995(17)260728(10)KMH21H0033
	KMH21H0038	(01)04026704347995(17)260728(10)KMH21H0038
	KMH21H0037	(01)04026704347995(17)260728(10)KMH21H0037
	KMH21H0055	(01)04026704347995(17)260728(10)KMH21H0055
	KMH21H0045	(01)04026704347995(17)260728(10)KMH21H0045
	KMH21H0046	(01)04026704347995(17)260728(10)KMH21H0046
	KMH21H0044	(01)04026704347995(17)260728(10)KMH21H0044
	KMH21H0053	(01)04026704347995(17)260728(10)KMH21H0053
	KMH21H0052	(01)04026704347995(17)260728(10)KMH21H0052
	KMH21H0161	(01)04026704347995(17)260728(10)KMH21H0161
	KMH21H0159	(01)04026704347995(17)260728(10)KMH21H0159
	KMH21H0167	(01)04026704347995(17)260728(10)KMH21H0167
	KMH21H0170	(01)04026704347995(17)260728(10)KMH21H0170
	KMH21H0200	(01)04026704347995(17)260728(10)KMH21H0200
	KMH21H0208	(01)04026704347995(17)260728(10)KMH21H0208
	KMH21H0210	(01)04026704347995(17)260728(10)KMH21H0210
	KMH21H0221	(01)04026704347995(17)260728(10)KMH21H0221
	KMH21H0223	(01)04026704347995(17)260728(10)KMH21H0223
	KMH21J0157	(01)04026704347995(17)260828(10)KMH21J0157
	KMH21J0164	(01)04026704347995(17)260828(10)KMH21J0164
	KMH21J0165	(01)04026704347995(17)260828(10)KMH21J0165
	KMH21J0167	(01)04026704347995(17)260828(10)KMH21J0167
	KMH21J0170	(01)04026704347995(17)260828(10)KMH21J0170
	KMH21J0168	(01)04026704347995(17)260828(10)KMH21J0168
KMH21J0171	(01)04026704347995(17)260828(10)KMH21J0171	
KMH21J0177	(01)04026704347995(17)260828(10)KMH21J0177	
KMH21J0175	(01)04026704347995(17)260828(10)KMH21J0175	
KMH21J0186	(01)04026704347995(17)260828(10)KMH21J0186	
KMH21J0181	(01)04026704347995(17)260828(10)KMH21J0181	
KMH21J0188	(01)04026704347995(17)260828(10)KMH21J0188	
KMH21J0191	(01)04026704347995(17)260828(10)KMH21J0191	
KMH21J0192	(01)04026704347995(17)260828(10)KMH21J0192	



**Appendix 2 - Product Code / Lot Number / Unique Device Identifier (UDI)**

Product Code	Batch	UDI Number
19211	KMH21J0195	(01)04026704347995(17)260828(10)KMH21J0195
	KMH21J0199	(01)04026704347995(17)260828(10)KMH21J0199
	KMH21J0207	(01)04026704347995(17)260828(10)KMH21J0207
	KMH21J0201	(01)04026704347995(17)260828(10)KMH21J0201
	KMH21J0204	(01)04026704347995(17)260828(10)KMH21J0204
	KMH21J0286	(01)04026704347995(17)260828(10)KMH21J0286
	KMH21J0239	(01)04026704347995(17)260828(10)KMH21J0239
	KMH21J0240	(01)04026704347995(17)260828(10)KMH21J0240
	KMH21J0255	(01)04026704347995(17)260828(10)KMH21J0255
	KMH21J0256	(01)04026704347995(17)260828(10)KMH21J0256
	KMH21J0262	(01)04026704347995(17)260828(10)KMH21J0262
	KMH21J0265	(01)04026704347995(17)260828(10)KMH21J0265
	KMH21J0268	(01)04026704347995(17)260828(10)KMH21J0268
	KMH21J0287	(01)04026704347995(17)260828(10)KMH21J0287
	KMH21J0289	(01)04026704347995(17)260828(10)KMH21J0289
	KMH21J0305	(01)04026704347995(17)260828(10)KMH21J0305
	KMH21J0291	(01)04026704347995(17)260828(10)KMH21J0291
	KMH21J0294	(01)04026704347995(17)260828(10)KMH21J0294
	KMH21J0306	(01)04026704347995(17)260828(10)KMH21J0306
	KMH21J0312	(01)04026704347995(17)260828(10)KMH21J0312
	KMH21J0317	(01)04026704347995(17)260828(10)KMH21J0317
	KMH21J0314	(01)04026704347995(17)260828(10)KMH21J0314
	KMH21J0333	(01)04026704347995(17)260828(10)KMH21J0333
	KMH21J0334	(01)04026704347995(17)260828(10)KMH21J0334
	KMH21K0001	(01)04026704347995(17)260928(10)KMH21K0001
	KMH21K0007	(01)04026704347995(17)260928(10)KMH21K0007
	KMH21K0003	(01)04026704347995(17)260928(10)KMH21K0003
	KMH21K0011	(01)04026704347995(17)260928(10)KMH21K0011
	KMH21K0015	(01)04026704347995(17)260928(10)KMH21K0015
	KMH21K0019	(01)04026704347995(17)260928(10)KMH21K0019
	KMH21K0018	(01)04026704347995(17)260928(10)KMH21K0018
	KMH21K0030	(01)04026704347995(17)260928(10)KMH21K0030
	KMH21K0037	(01)04026704347995(17)260928(10)KMH21K0037
	KMH21K0078	(01)04026704347995(17)260928(10)KMH21K0078
	KMH21K0080	(01)04026704347995(17)260928(10)KMH21K0080
	KMH21K0092	(01)04026704347995(17)260928(10)KMH21K0092
	KMH21K0100	(01)04026704347995(17)260928(10)KMH21K0100
	KMH21K0112	(01)04026704347995(17)260928(10)KMH21K0112
	KMH21K0182	(01)04026704347995(17)260928(10)KMH21K0182
	KMH21K0237	(01)04026704347995(17)260928(10)KMH21K0237
	KMH21K0403	(01)04026704347995(17)260928(10)KMH21K0403
	KMH21K0407	(01)04026704347995(17)260928(10)KMH21K0407
	KMH21K0408	(01)04026704347995(17)260928(10)KMH21K0408
	KMH21K0409	(01)04026704347995(17)260928(10)KMH21K0409
	KMH21L0012	(01)04026704347995(17)261028(10)KMH21L0012
	KMH21L0370	(01)04026704347995(17)261028(10)KMH21L0370
	KMH21M0054	(01)04026704347995(17)261128(10)KMH21M0054
	KMH21M0103	(01)04026704347995(17)261128(10)KMH21M0103
	KMH21M0206	(01)04026704347995(17)261128(10)KMH21M0206
	KMH21M0135	(01)04026704347995(17)261128(10)KMH21M0135
KMH21M0228	(01)04026704347995(17)261128(10)KMH21M0228	
KMH21M0198	(01)04026704347995(17)261128(10)KMH21M0198	
KMH22A0013	(01)04026704347995(17)261228(10)KMH22A0013	
KMH22A0043	(01)04026704347995(17)261228(10)KMH22A0043	
KMH22A0095	(01)04026704347995(17)261228(10)KMH22A0095	
KMH22A0326	(01)04026704347995(17)261228(10)KMH22A0326	
KMH22A0270	(01)04026704347995(17)261228(10)KMH22A0270	
KMH22B0092	(01)04026704347995(17)270128(10)KMH22B0092	

**Appendix 2 - Product Code / Lot Number / Unique Device Identifier (UDI)**

Product Code	Batch	UDI Number
19211	KMH22B0206	(01)04026704347995(17)270128(10)KMH22B0206
	KMH21L0021	(01)04026704347995(17)261028(10)KMH21L0021
	KMH21L0023	(01)04026704347995(17)261028(10)KMH21L0023
	KMH21L0024	(01)04026704347995(17)261028(10)KMH21L0024
	KMH21L0061	(01)04026704347995(17)261028(10)KMH21L0061
	KMH21L0169	(01)04026704347995(17)261028(10)KMH21L0169
	KMH21L0178	(01)04026704347995(17)261028(10)KMH21L0178
	KMH21L0182	(01)04026704347995(17)261028(10)KMH21L0182
	KMH21L0190	(01)04026704347995(17)261028(10)KMH21L0190
	KMH21L0211	(01)04026704347995(17)261028(10)KMH21L0211
	KMH21L0213	(01)04026704347995(17)261028(10)KMH21L0213
	KMH21M0002	(01)04026704347995(17)261128(10)KMH21M0002
	KMH21M0008	(01)04026704347995(17)261128(10)KMH21M0008
	KMH21M0016	(01)04026704347995(17)261128(10)KMH21M0016
	KMH21M0014	(01)04026704347995(17)261128(10)KMH21M0014
	KMH21M0015	(01)04026704347995(17)261128(10)KMH21M0015
	KMH21M0018	(01)04026704347995(17)261128(10)KMH21M0018
	KMH21M0019	(01)04026704347995(17)261128(10)KMH21M0019
	KMH21M0020	(01)04026704347995(17)261128(10)KMH21M0020
	KMH21M0031	(01)04026704347995(17)261128(10)KMH21M0031
	KMH21M0065	(01)04026704347995(17)261128(10)KMH21M0065
	KMH21M0055	(01)04026704347995(17)261128(10)KMH21M0055
	KMH21M0056	(01)04026704347995(17)261128(10)KMH21M0056
	KMH21M0057	(01)04026704347995(17)261128(10)KMH21M0057
	KMH21M0069	(01)04026704347995(17)261128(10)KMH21M0069
	KMH21M0071	(01)04026704347995(17)261128(10)KMH21M0071
	KMH21M0082	(01)04026704347995(17)261128(10)KMH21M0082
	KMH21M0106	(01)04026704347995(17)261128(10)KMH21M0106
	KMH21M0111	(01)04026704347995(17)261128(10)KMH21M0111
	KMH21M0145	(01)04026704347995(17)261128(10)KMH21M0145
	KMH21M0144	(01)04026704347995(17)261128(10)KMH21M0144
	KMH21M0150	(01)04026704347995(17)261128(10)KMH21M0150
	KMH21M0157	(01)04026704347995(17)261128(10)KMH21M0157
	KMH21M0162	(01)04026704347995(17)261128(10)KMH21M0162
	KMH21M0183	(01)04026704347995(17)261128(10)KMH21M0183
	KMH21M0185	(01)04026704347995(17)261128(10)KMH21M0185
	KMH21M0202	(01)04026704347995(17)261128(10)KMH21M0202
	KMH21M0204	(01)04026704347995(17)261128(10)KMH21M0204
	KMH21M0208	(01)04026704347995(17)261128(10)KMH21M0208
	KMH21M0215	(01)04026704347995(17)261128(10)KMH21M0215
	KMH21M0216	(01)04026704347995(17)261128(10)KMH21M0216
	KMH21M0226	(01)04026704347995(17)261128(10)KMH21M0226
	KMH21M0232	(01)04026704347995(17)261128(10)KMH21M0232
	KMH21M0231	(01)04026704347995(17)261128(10)KMH21M0231
	KMH21M0244	(01)04026704347995(17)261128(10)KMH21M0244
	KMH21M0242	(01)04026704347995(17)261128(10)KMH21M0242
	KMH21M0243	(01)04026704347995(17)261128(10)KMH21M0243
	KMH21M0249	(01)04026704347995(17)261128(10)KMH21M0249
	KMH21M0250	(01)04026704347995(17)261128(10)KMH21M0250
	KMH22A0005	(01)04026704347995(17)261228(10)KMH22A0005
KMH22A0046	(01)04026704347995(17)261228(10)KMH22A0046	
KMH22A0047	(01)04026704347995(17)261228(10)KMH22A0047	
KMH22A0051	(01)04026704347995(17)261228(10)KMH22A0051	
KMH22A0070	(01)04026704347995(17)261228(10)KMH22A0070	
KMH22A0081	(01)04026704347995(17)261228(10)KMH22A0081	
KMH22A0078	(01)04026704347995(17)261228(10)KMH22A0078	
KMH22A0082	(01)04026704347995(17)261228(10)KMH22A0082	
KMH22A0096	(01)04026704347995(17)261228(10)KMH22A0096	

**Appendix 2 - Product Code / Lot Number / Unique Device Identifier (UDI)**

Product Code	Batch	UDI Number
19211	KMH22A0097	(01)04026704347995(17)261228(10)KMH22A0097
	KMH22A0124	(01)04026704347995(17)261228(10)KMH22A0124
	KMH22A0113	(01)04026704347995(17)261228(10)KMH22A0113
	KMH22A0126	(01)04026704347995(17)261228(10)KMH22A0126
	KMH22A0129	(01)04026704347995(17)261228(10)KMH22A0129
	KMH22A0144	(01)04026704347995(17)261228(10)KMH22A0144
	KMH22A0142	(01)04026704347995(17)261228(10)KMH22A0142
	KMH22A0146	(01)04026704347995(17)261228(10)KMH22A0146
	KMH22A0151	(01)04026704347995(17)261228(10)KMH22A0151
	KMH22A0155	(01)04026704347995(17)261228(10)KMH22A0155
	KMH22A0281	(01)04026704347995(17)261228(10)KMH22A0281
	KMH22A0283	(01)04026704347995(17)261228(10)KMH22A0283
	KMH22A0286	(01)04026704347995(17)261228(10)KMH22A0286
	KMH22A0297	(01)04026704347995(17)261228(10)KMH22A0297
	KMH22A0300	(01)04026704347995(17)261228(10)KMH22A0300
	KMH22A0319	(01)04026704347995(17)261228(10)KMH22A0319
	KMH22A0327	(01)04026704347995(17)261228(10)KMH22A0327
	KMH22A0328	(01)04026704347995(17)261228(10)KMH22A0328
	KMH22A0332	(01)04026704347995(17)261228(10)KMH22A0332
	KMH22B0006	(01)04026704347995(17)270128(10)KMH22B0006
	KMH22B0001	(01)04026704347995(17)270128(10)KMH22B0001
	KMH22B0003	(01)04026704347995(17)270128(10)KMH22B0003
	KMH22B0094	(01)04026704347995(17)270128(10)KMH22B0094
	KMH22B0095	(01)04026704347995(17)270128(10)KMH22B0095
	KMH22B0096	(01)04026704347995(17)270128(10)KMH22B0096
	KMH22B0113	(01)04026704347995(17)270128(10)KMH22B0113
	KMH22B0112	(01)04026704347995(17)270128(10)KMH22B0112
	KMH22B0123	(01)04026704347995(17)270128(10)KMH22B0123
	KMH22B0127	(01)04026704347995(17)270128(10)KMH22B0127
	KMH22B0128	(01)04026704347995(17)270128(10)KMH22B0128
	KMH22B0129	(01)04026704347995(17)270128(10)KMH22B0129
	KMH22B0186	(01)04026704347995(17)270128(10)KMH22B0186
	KMH22B0194	(01)04026704347995(17)270128(10)KMH22B0194
	KMH22B0199	(01)04026704347995(17)270128(10)KMH22B0199
	KMH22B0210	(01)04026704347995(17)270128(10)KMH22B0210
	KMH22B0214	(01)04026704347995(17)270128(10)KMH22B0214
	KMH22C0058	(01)04026704347995(17)270228(10)KMH22C0058
	KMH22C0079	(01)04026704347995(17)270228(10)KMH22C0079
	KMH22C0080	(01)04026704347995(17)270228(10)KMH22C0080
	KMH21L0225	(01)04026704347995(17)261028(10)KMH21L0225
19212	202024	(01)04026704348008(17)250528(10)202024
	KMZ20K0591	(01)04026704348008(17)250928(10)KMZ20K0591
	KMZ20K1179	(01)04026704348008(17)250928(10)KMZ20K1179
	201951	(01)04026704348008(17)241128(10)201951
	KMZ20L1178	(01)04026704348008(17)251028(10)KMZ20L1178
	KMZ21A0489	(01)04026704348008(17)251228(10)KMZ21A0489
	KMZ21A1126	(01)04026704348008(17)251228(10)KMZ21A1126
	KMZ21B1183	(01)04026704348008(17)260128(10)KMZ21B1183
	KMZ21C0870	(01)04026704348008(17)260228(10)KMZ21C0870
	KMZ21D0970	(01)04026704348008(17)260328(10)KMZ21D0970
	KMZ21D1114	(01)04026704348008(17)260328(10)KMZ21D1114
	KMZ21E0016	(01)04026704348008(17)260428(10)KMZ21E0016
	KMZ21E0629	(01)04026704348008(17)260428(10)KMZ21E0629
	KMZ21E0635	(01)04026704348008(17)260428(10)KMZ21E0635
	KMZ21G0003	(01)04026704348008(17)260628(10)KMZ21G0003
	KMZ21H0247	(01)04026704348008(17)260728(10)KMZ21H0247
	202027	(01)04026704348008(17)250628(10)202027
	KMZ21J0586	(01)04026704348008(17)260828(10)KMZ21J0586

**Appendix 2 - Product Code / Lot Number / Unique Device Identifier (UDI)**

Product Code	Batch	UDI Number
19212	KMZ20H0517	(01)04026704348008(17)250728(10)KMZ20H0517
	KMZ21J0662	(01)04026704348008(17)260828(10)KMZ21J0662
	KMZ20H0852	(01)04026704348008(17)250728(10)KMZ20H0852
	KMZ21K0046	(01)04026704348008(17)260928(10)KMZ21K0046
	KMZ21K0106	(01)04026704348008(17)260928(10)KMZ21K0106
	KMZ21K0183	(01)04026704348008(17)260928(10)KMZ21K0183
	KMZ21L0888	(01)04026704348008(17)261028(10)KMZ21L0888
	KMZ20K0594	(01)04026704348008(17)250928(10)KMZ20K0594
	KMZ20K0725	(01)04026704348008(17)250928(10)KMZ20K0725
	KMZ20K0805	(01)04026704348008(17)250928(10)KMZ20K0805
	KMZ20K1320	(01)04026704348008(17)250928(10)KMZ20K1320
	KMZ20L0517	(01)04026704348008(17)251028(10)KMZ20L0517
	KMZ20L0522	(01)04026704348008(17)251028(10)KMZ20L0522
	KMZ20L0550	(01)04026704348008(17)251028(10)KMZ20L0550
	KMZ20M0251	(01)04026704348008(17)251128(10)KMZ20M0251
	KMZ20M0849	(01)04026704348008(17)251128(10)KMZ20M0849
	KMZ21A0538	(01)04026704348008(17)251228(10)KMZ21A0538
	KMZ21A0587	(01)04026704348008(17)251228(10)KMZ21A0587
	KMZ21A0622	(01)04026704348008(17)251228(10)KMZ21A0622
	KMZ21A0698	(01)04026704348008(17)251228(10)KMZ21A0698
	KMZ21A0633	(01)04026704348008(17)251228(10)KMZ21A0633
	KMZ21A0644	(01)04026704348008(17)251228(10)KMZ21A0644
	KMZ21A0652	(01)04026704348008(17)251228(10)KMZ21A0652
	KMZ21A0707	(01)04026704348008(17)251228(10)KMZ21A0707
	KMZ21A0709	(01)04026704348008(17)251228(10)KMZ21A0709
	KMZ21A0726	(01)04026704348008(17)251228(10)KMZ21A0726
	KMZ21A1007	(01)04026704348008(17)251228(10)KMZ21A1007
	KMZ21A1086	(01)04026704348008(17)251228(10)KMZ21A1086
	KMZ21B0565	(01)04026704348008(17)260128(10)KMZ21B0565
	KMZ21B0612	(01)04026704348008(17)260128(10)KMZ21B0612
	KMZ21B0780	(01)04026704348008(17)260128(10)KMZ21B0780
	KMZ21B1131	(01)04026704348008(17)260128(10)KMZ21B1131
	KMZ21B1196	(01)04026704348008(17)260128(10)KMZ21B1196
	KMZ21C0027	(01)04026704348008(17)260228(10)KMZ21C0027
	KMZ21C0059	(01)04026704348008(17)260228(10)KMZ21C0059
	KMZ21C0564	(01)04026704348008(17)260228(10)KMZ21C0564
	KMZ21C0908	(01)04026704348008(17)260228(10)KMZ21C0908
	KMZ21D0937	(01)04026704348008(17)260328(10)KMZ21D0937
	KMZ21D0942	(01)04026704348008(17)260328(10)KMZ21D0942
	KMZ21D0994	(01)04026704348008(17)260328(10)KMZ21D0994
	KMZ21D0980	(01)04026704348008(17)260328(10)KMZ21D0980
	KMZ21D1059	(01)04026704348008(17)260328(10)KMZ21D1059
	KMZ21D0995	(01)04026704348008(17)260328(10)KMZ21D0995
	KMZ21D1027	(01)04026704348008(17)260328(10)KMZ21D1027
	KMZ21D1066	(01)04026704348008(17)260328(10)KMZ21D1066
	KMZ21D1078	(01)04026704348008(17)260328(10)KMZ21D1078
	KMZ21D1092	(01)04026704348008(17)260328(10)KMZ21D1092
	KMZ21E0021	(01)04026704348008(17)260428(10)KMZ21E0021
	KMZ21E0047	(01)04026704348008(17)260428(10)KMZ21E0047
	KMZ21E0054	(01)04026704348008(17)260428(10)KMZ21E0054
KMZ21E0075	(01)04026704348008(17)260428(10)KMZ21E0075	
KMZ21E0083	(01)04026704348008(17)260428(10)KMZ21E0083	
KMZ21E0099	(01)04026704348008(17)260428(10)KMZ21E0099	
KMZ21E0118	(01)04026704348008(17)260428(10)KMZ21E0118	
KMZ21E0120	(01)04026704348008(17)260428(10)KMZ21E0120	
KMZ21E0617	(01)04026704348008(17)260428(10)KMZ21E0617	
KMZ21E0647	(01)04026704348008(17)260428(10)KMZ21E0647	
KMZ21E0662	(01)04026704348008(17)260428(10)KMZ21E0662	

**Appendix 2 - Product Code / Lot Number / Unique Device Identifier (UDI)**

Product Code	Batch	UDI Number
19212	KMZ21E0656	(01)04026704348008(17)260428(10)KMZ21E0656
	KMZ21E0666	(01)04026704348008(17)260428(10)KMZ21E0666
	KMZ21F0965	(01)04026704348008(17)260528(10)KMZ21F0965
	KMZ21H0304	(01)04026704348008(17)260728(10)KMZ21H0304
	KMZ21H0291	(01)04026704348008(17)260728(10)KMZ21H0291
	KMZ21H0297	(01)04026704348008(17)260728(10)KMZ21H0297
	KMZ21H0308	(01)04026704348008(17)260728(10)KMZ21H0308
	KMZ21H0312	(01)04026704348008(17)260728(10)KMZ21H0312
	KMZ21H0333	(01)04026704348008(17)260728(10)KMZ21H0333
	KMZ21H0338	(01)04026704348008(17)260728(10)KMZ21H0338
	KMZ21H0378	(01)04026704348008(17)260728(10)KMZ21H0378
	KMZ21H0384	(01)04026704348008(17)260728(10)KMZ21H0384
	KMZ21J0562	(01)04026704348008(17)260828(10)KMZ21J0562
	KMZ21J0563	(01)04026704348008(17)260828(10)KMZ21J0563
	KMZ21J0587	(01)04026704348008(17)260828(10)KMZ21J0587
	KMZ21J0701	(01)04026704348008(17)260828(10)KMZ21J0701
	KMZ21J0703	(01)04026704348008(17)260828(10)KMZ21J0703
	KMZ21J0713	(01)04026704348008(17)260828(10)KMZ21J0713
	KMZ21J0708	(01)04026704348008(17)260828(10)KMZ21J0708
	KMZ21J0725	(01)04026704348008(17)260828(10)KMZ21J0725
	KMZ21J0728	(01)04026704348008(17)260828(10)KMZ21J0728
	KMZ21K0075	(01)04026704348008(17)260928(10)KMZ21K0075
	KMZ21K0064	(01)04026704348008(17)260928(10)KMZ21K0064
	KMZ21K0076	(01)04026704348008(17)260928(10)KMZ21K0076
	KMZ21K0084	(01)04026704348008(17)260928(10)KMZ21K0084
	KMZ21K0098	(01)04026704348008(17)260928(10)KMZ21K0098
	KMZ21K0161	(01)04026704348008(17)260928(10)KMZ21K0161
	KMZ21K0182	(01)04026704348008(17)260928(10)KMZ21K0182
	KMZ21K0530	(01)04026704348008(17)260928(10)KMZ21K0530
	KMZ21K0597	(01)04026704348008(17)260928(10)KMZ21K0597
	KMZ21K0598	(01)04026704348008(17)260928(10)KMZ21K0598
	KMZ21K0611	(01)04026704348008(17)260928(10)KMZ21K0611
	KMZ21K0613	(01)04026704348008(17)260928(10)KMZ21K0613
	KMZ21K0945	(01)04026704348008(17)260928(10)KMZ21K0945
	KMZ21M0255	(01)04026704348008(17)261128(10)KMZ21M0255
	KMZ21M0355	(01)04026704348008(17)261128(10)KMZ21M0355
	KMZ21M0512	(01)04026704348008(17)261128(10)KMZ21M0512
	KMZ22A0124	(01)04026704348008(17)261228(10)KMZ22A0124
	KMZ21M0888	(01)04026704348008(17)261128(10)KMZ21M0888
	KMZ21M0989	(01)04026704348008(17)261128(10)KMZ21M0989
	KMZ22A0013	(01)04026704348008(17)261228(10)KMZ22A0013
	KMZ22A0264	(01)04026704348008(17)261228(10)KMZ22A0264
	KMZ22A0904	(01)04026704348008(17)261228(10)KMZ22A0904
	KMZ22B0263	(01)04026704348008(17)270128(10)KMZ22B0263
	KMZ22B0001	(01)04026704348008(17)270128(10)KMZ22B0001
	KMZ22B0596	(01)04026704348008(17)270128(10)KMZ22B0596
	KMZ22B0743	(01)04026704348008(17)270128(10)KMZ22B0743
	KMZ21L0892	(01)04026704348008(17)261028(10)KMZ21L0892
	KMZ21L0900	(01)04026704348008(17)261028(10)KMZ21L0900
	KMZ21M0261	(01)04026704348008(17)261128(10)KMZ21M0261
KMZ21M0311	(01)04026704348008(17)261128(10)KMZ21M0311	
KMZ21M0318	(01)04026704348008(17)261128(10)KMZ21M0318	
KMZ21M0332	(01)04026704348008(17)261128(10)KMZ21M0332	
KMZ21M0333	(01)04026704348008(17)261128(10)KMZ21M0333	
KMZ21M0338	(01)04026704348008(17)261128(10)KMZ21M0338	
KMZ21M0346	(01)04026704348008(17)261128(10)KMZ21M0346	
KMZ21M0363	(01)04026704348008(17)261128(10)KMZ21M0363	
KMZ21M0380	(01)04026704348008(17)261128(10)KMZ21M0380	

**Appendix 2 - Product Code / Lot Number / Unique Device Identifier (UDI)**

Product Code	Batch	UDI Number
19212	KMZ21M0470	(01)04026704348008(17)261128(10)KMZ21M0470
	KMZ21M0530	(01)04026704348008(17)261128(10)KMZ21M0530
	KMZ21M0639	(01)04026704348008(17)261128(10)KMZ21M0639
	KMZ21M0591	(01)04026704348008(17)261128(10)KMZ21M0591
	KMZ21M0707	(01)04026704348008(17)261128(10)KMZ21M0707
	KMZ21M0724	(01)04026704348008(17)261128(10)KMZ21M0724
	KMZ21M0896	(01)04026704348008(17)261128(10)KMZ21M0896
	KMZ21M0906	(01)04026704348008(17)261128(10)KMZ21M0906
	KMZ21M0914	(01)04026704348008(17)261128(10)KMZ21M0914
	KMZ21M0943	(01)04026704348008(17)261128(10)KMZ21M0943
	KMZ21M0954	(01)04026704348008(17)261128(10)KMZ21M0954
	KMZ21M1006	(01)04026704348008(17)261128(10)KMZ21M1006
	KMZ21M1053	(01)04026704348008(17)261128(10)KMZ21M1053
	KMZ21M1037	(01)04026704348008(17)261128(10)KMZ21M1037
	KMZ22A0016	(01)04026704348008(17)261228(10)KMZ22A0016
	KMZ22A0021	(01)04026704348008(17)261228(10)KMZ22A0021
	KMZ22A0025	(01)04026704348008(17)261228(10)KMZ22A0025
	KMZ22A0068	(01)04026704348008(17)261228(10)KMZ22A0068
	KMZ22A0087	(01)04026704348008(17)261228(10)KMZ22A0087
	KMZ22A0132	(01)04026704348008(17)261228(10)KMZ22A0132
	KMZ22A0280	(01)04026704348008(17)261228(10)KMZ22A0280
	KMZ22A0439	(01)04026704348008(17)261228(10)KMZ22A0439
	KMZ22A0661	(01)04026704348008(17)261228(10)KMZ22A0661
	KMZ22A0713	(01)04026704348008(17)261228(10)KMZ22A0713
	KMZ22A0730	(01)04026704348008(17)261228(10)KMZ22A0730
	KMZ22A0744	(01)04026704348008(17)261228(10)KMZ22A0744
	KMZ22A0755	(01)04026704348008(17)261228(10)KMZ22A0755
	KMZ22A0767	(01)04026704348008(17)261228(10)KMZ22A0767
	KMZ22A0777	(01)04026704348008(17)261228(10)KMZ22A0777
	KMZ22A1002	(01)04026704348008(17)261228(10)KMZ22A1002
	KMZ22A1017	(01)04026704348008(17)261228(10)KMZ22A1017
	KMZ22B0016	(01)04026704348008(17)270128(10)KMZ22B0016
	KMZ22B0064	(01)04026704348008(17)270128(10)KMZ22B0064
	KMZ22B0022	(01)04026704348008(17)270128(10)KMZ22B0022
	KMZ22B0071	(01)04026704348008(17)270128(10)KMZ22B0071
	KMZ22B0080	(01)04026704348008(17)270128(10)KMZ22B0080
	KMZ22B0090	(01)04026704348008(17)270128(10)KMZ22B0090
	KMZ22B0097	(01)04026704348008(17)270128(10)KMZ22B0097
	KMZ22B0106	(01)04026704348008(17)270128(10)KMZ22B0106
	KMZ22B0631	(01)04026704348008(17)270128(10)KMZ22B0631
KMZ22B0636	(01)04026704348008(17)270128(10)KMZ22B0636	
KMZ22B0651	(01)04026704348008(17)270128(10)KMZ22B0651	
KMZ22C0005	(01)04026704348008(17)270228(10)KMZ22C0005	
KMZ22C0117	(01)04026704348008(17)270228(10)KMZ22C0117	
19261	KMH21B0193	(01)04026704388486(17)260128(10)KMH21B0193
	KMH20L0234	(01)04026704388486(17)251028(10)KMH20L0234
	KMH21J0246	(01)04026704388486(17)260828(10)KMH21J0246
	KMH21L0294	(01)04026704388486(17)261028(10)KMH21L0294
	KMH22A0230	(01)04026704388486(17)261228(10)KMH22A0230
19272	KMZ21E0613	(01)04026704388509(17)260428(10)KMZ21E0613
	KMZ21F0939	(01)04026704388509(17)260528(10)KMZ21F0939
19012T	KMZ21A0712	(01)04026704400423(17)251228(10)KMZ21A0712
	KMZ21A0909	(01)04026704400423(17)251228(10)KMZ21A0909
	KMZ21B0141	(01)04026704400423(17)260128(10)KMZ21B0141
	KMZ21G0174	(01)04026704400423(17)260628(10)KMZ21G0174
	KMZ21G0180	(01)04026704400423(17)260628(10)KMZ21G0180
	KMZ21G0181	(01)04026704400423(17)260628(10)KMZ21G0181
	KMZ21H0277	(01)04026704400423(17)260728(10)KMZ21H0277

**Appendix 2 - Product Code / Lot Number / Unique Device Identifier (UDI)**

Product Code	Batch	UDI Number
19012T	KMZ21K0208	(01)04026704400423(17)260928(10)KMZ21K0208
	KMZ21K0648	(01)04026704400423(17)260928(10)KMZ21K0648
	KMZ21K0707	(01)04026704400423(17)260928(10)KMZ21K0707
	KMZ21L0309	(01)04026704400423(17)261028(10)KMZ21L0309
	KMZ21K0856	(01)04026704400423(17)260928(10)KMZ21K0856
	KMZ21M0061	(01)04026704400423(17)261128(10)KMZ21M0061
	KMZ21A0427	(01)04026704400423(17)251228(10)KMZ21A0427
	KMZ21A0522	(01)04026704400423(17)251228(10)KMZ21A0522
	KMZ21A0535	(01)04026704400423(17)251228(10)KMZ21A0535
	KMZ21A0617	(01)04026704400423(17)251228(10)KMZ21A0617
	KMZ21A0926	(01)04026704400423(17)251228(10)KMZ21A0926
	KMZ21J0386	(01)04026704400423(17)260828(10)KMZ21J0386
	KMZ21M0813	(01)04026704400423(17)261128(10)KMZ21M0813
	KMZ22B0105	(01)04026704400423(17)270128(10)KMZ22B0105
	KMZ22B0118	(01)04026704400423(17)270128(10)KMZ22B0118
KMZ22B0117	(01)04026704400423(17)270128(10)KMZ22B0117	
191667-000100	KMH21B0197	(01)04026704561544(17)260128(10)KMH21B0197
	KMH22A0181	(01)04026704561544(17)261228(10)KMH22A0181
	KMH22A0185	(01)04026704561544(17)261228(10)KMH22A0185
19211T	KMH20L0362	(01)04026704400478(17)251028(10)KMH20L0362
	KMH20M0276	(01)04026704400478(17)251128(10)KMH20M0276
	KMH21A0397	(01)04026704400478(17)251228(10)KMH21A0397
	KMH21C0228	(01)04026704400478(17)260228(10)KMH21C0228
	KMH21K0122	(01)04026704400478(17)260928(10)KMH21K0122
	KMH21K0130	(01)04026704400478(17)260928(10)KMH21K0130
19212T	KMZ20K0917	(01)04026704400485(17)250928(10)KMZ20K0917
	KMZ20L0912	(01)04026704400485(17)251028(10)KMZ20L0912
	KMZ21A1017	(01)04026704400485(17)251228(10)KMZ21A1017
	KMZ21B0688	(01)04026704400485(17)260128(10)KMZ21B0688
	KMZ21B1156	(01)04026704400485(17)260128(10)KMZ21B1156
	KMZ21C0760	(01)04026704400485(17)260228(10)KMZ21C0760
	KMZ21C1163	(01)04026704400485(17)260228(10)KMZ21C1163
	KMZ21D1051	(01)04026704400485(17)260328(10)KMZ21D1051
	KMZ21E0245	(01)04026704400485(17)260428(10)KMZ21E0245
	KMZ21E0420	(01)04026704400485(17)260428(10)KMZ21E0420
	KMZ21F0202	(01)04026704400485(17)260528(10)KMZ21F0202
	KMZ21F0773	(01)04026704400485(17)260528(10)KMZ21F0773
	KMZ20K1295	(01)04026704400485(17)250928(10)KMZ20K1295
	KMZ20K1311	(01)04026704400485(17)250928(10)KMZ20K1311
	KMZ20L0944	(01)04026704400485(17)251028(10)KMZ20L0944
	KMZ21A0533	(01)04026704400485(17)251228(10)KMZ21A0533
	KMZ21A0867	(01)04026704400485(17)251228(10)KMZ21A0867
	KMZ21B1188	(01)04026704400485(17)260128(10)KMZ21B1188
	KMZ21C0712	(01)04026704400485(17)260228(10)KMZ21C0712
	KMZ21C0882	(01)04026704400485(17)260228(10)KMZ21C0882
	KMZ21D1020	(01)04026704400485(17)260328(10)KMZ21D1020
	KMZ21E0253	(01)04026704400485(17)260428(10)KMZ21E0253
	KMZ21F0786	(01)04026704400485(17)260528(10)KMZ21F0786
	KMZ21G0060	(01)04026704400485(17)260628(10)KMZ21G0060
	KMZ21G0353	(01)04026704400485(17)260628(10)KMZ21G0353
	KMZ21J0403	(01)04026704400485(17)260828(10)KMZ21J0403
	KMZ21J0407	(01)04026704400485(17)260828(10)KMZ21J0407
	KMZ21K0113	(01)04026704400485(17)260928(10)KMZ21K0113
	KMZ21K0223	(01)04026704400485(17)260928(10)KMZ21K0223
	KMZ21L0114	(01)04026704400485(17)261028(10)KMZ21L0114
	KMZ21L0321	(01)04026704400485(17)261028(10)KMZ21L0321
	KMZ21G0360	(01)04026704400485(17)260628(10)KMZ21G0360
	KMZ21G0392	(01)04026704400485(17)260628(10)KMZ21G0392

**Appendix 2 - Product Code / Lot Number / Unique Device Identifier (UDI)**

Product Code	Batch	UDI Number
19212T	KMZ21J0428	(01)04026704400485(17)260828(10)KMZ21J0428
	KMZ21J0435	(01)04026704400485(17)260828(10)KMZ21J0435
	KMZ21K0028	(01)04026704400485(17)260928(10)KMZ21K0028
	KMZ21K0191	(01)04026704400485(17)260928(10)KMZ21K0191
	KMZ21K0230	(01)04026704400485(17)260928(10)KMZ21K0230
	KMZ21K0239	(01)04026704400485(17)260928(10)KMZ21K0239
	KMZ21K0260	(01)04026704400485(17)260928(10)KMZ21K0260
	KMZ21K0329	(01)04026704400485(17)260928(10)KMZ21K0329
	KMZ21K0610	(01)04026704400485(17)260928(10)KMZ21K0610
	KMZ21K0552	(01)04026704400485(17)260928(10)KMZ21K0552
	KMZ21K0562	(01)04026704400485(17)260928(10)KMZ21K0562
	KMZ21K0579	(01)04026704400485(17)260928(10)KMZ21K0579
	KMZ21M0066	(01)04026704400485(17)261128(10)KMZ21M0066
	KMZ21M0476	(01)04026704400485(17)261128(10)KMZ21M0476
	KMZ22A0035	(01)04026704400485(17)261228(10)KMZ22A0035
	KMZ22A0059	(01)04026704400485(17)261228(10)KMZ22A0059
	KMZ22A0471	(01)04026704400485(17)261228(10)KMZ22A0471
	KMZ22A0483	(01)04026704400485(17)261228(10)KMZ22A0483
	KMZ22B0043	(01)04026704400485(17)270128(10)KMZ22B0043
	KMZ22B0439	(01)04026704400485(17)270128(10)KMZ22B0439
	KMZ22B0871	(01)04026704400485(17)270128(10)KMZ22B0871
	KMZ21L0209	(01)04026704400485(17)261028(10)KMZ21L0209
	KMZ21L0536	(01)04026704400485(17)261028(10)KMZ21L0536
	KMZ21M0572	(01)04026704400485(17)261128(10)KMZ21M0572
	KMZ21M0486	(01)04026704400485(17)261128(10)KMZ21M0486
	KMZ21M0526	(01)04026704400485(17)261128(10)KMZ21M0526
	KMZ21M0586	(01)04026704400485(17)261128(10)KMZ21M0586
	KMZ21M0601	(01)04026704400485(17)261128(10)KMZ21M0601
	KMZ21M0606	(01)04026704400485(17)261128(10)KMZ21M0606
	KMZ21M0658	(01)04026704400485(17)261128(10)KMZ21M0658
	KMZ21M0742	(01)04026704400485(17)261128(10)KMZ21M0742
	KMZ21M0775	(01)04026704400485(17)261128(10)KMZ21M0775
	KMZ22A0076	(01)04026704400485(17)261228(10)KMZ22A0076
	KMZ22A0085	(01)04026704400485(17)261228(10)KMZ22A0085
KMZ22A0482	(01)04026704400485(17)261228(10)KMZ22A0482	
KMZ22B0044	(01)04026704400485(17)270128(10)KMZ22B0044	
KMZ22B0123	(01)04026704400485(17)270128(10)KMZ22B0123	
19261T	KMH20L0409	(01)04026704400492(17)251028(10)KMH20L0409
	KMH21B0343	(01)04026704400492(17)260128(10)KMH21B0343
	KMH21C0177	(01)04026704400492(17)260228(10)KMH21C0177
	KMH21C0387	(01)04026704400492(17)260228(10)KMH21C0387
	KMH21G0426	(01)04026704400492(17)260628(10)KMH21G0426
	KMH21M0182	(01)04026704400492(17)261128(10)KMH21M0182
KMH22B0151	(01)04026704400492(17)270128(10)KMH22B0151	
19262T	KMZ20M0661	(01)04026704400508(17)251128(10)KMZ20M0661
	KMZ21B0811	(01)04026704400508(17)260128(10)KMZ21B0811
	KMZ21B0918	(01)04026704400508(17)260128(10)KMZ21B0918
	KMZ21B1172	(01)04026704400508(17)260128(10)KMZ21B1172
	KMZ21C0929	(01)04026704400508(17)260228(10)KMZ21C0929
	KMZ21G0728	(01)04026704400508(17)260628(10)KMZ21G0728
	KMZ21G0735	(01)04026704400508(17)260628(10)KMZ21G0735
	KMZ21K0510	(01)04026704400508(17)260928(10)KMZ21K0510
	KMZ21K0592	(01)04026704400508(17)260928(10)KMZ21K0592
	KMZ21K0952	(01)04026704400508(17)260928(10)KMZ21K0952
	KMZ21L0218	(01)04026704400508(17)261028(10)KMZ21L0218
	KMZ20L0968	(01)04026704400508(17)251028(10)KMZ20L0968
	KMZ20L0977	(01)04026704400508(17)251028(10)KMZ20L0977
	KMZ20L1275	(01)04026704400508(17)251028(10)KMZ20L1275



**Appendix 2 - Product Code / Lot Number / Unique Device Identifier (UDI)**

Product Code	Batch	UDI Number
19262T	KMZ21C0884	(01)04026704400508(17)260228(10)KMZ21C0884
	KMZ21C0910	(01)04026704400508(17)260228(10)KMZ21C0910
	KMZ21G0739	(01)04026704400508(17)260628(10)KMZ21G0739
	KMZ21K0539	(01)04026704400508(17)260928(10)KMZ21K0539
	KMZ21K0524	(01)04026704400508(17)260928(10)KMZ21K0524
	KMZ21K0572	(01)04026704400508(17)260928(10)KMZ21K0572
	KMZ21K0578	(01)04026704400508(17)260928(10)KMZ21K0578
	KMZ22B0125	(01)04026704400508(17)270128(10)KMZ22B0125
	KMZ22B0345	(01)04026704400508(17)270128(10)KMZ22B0345
	KMZ22B0809	(01)04026704400508(17)270128(10)KMZ22B0809
	KMZ22B0869	(01)04026704400508(17)270128(10)KMZ22B0869
	KMZ22B0828	(01)04026704400508(17)270128(10)KMZ22B0828
	KMZ22B0838	(01)04026704400508(17)270128(10)KMZ22B0838
	19272T	KMZ20K0838
KMZ20L0885		(01)04026704400515(17)251028(10)KMZ20L0885
KMZ20L1056		(01)04026704400515(17)251028(10)KMZ20L1056
KMZ20M0777		(01)04026704400515(17)251128(10)KMZ20M0777
KMZ20M0387		(01)04026704400515(17)251128(10)KMZ20M0387
KMZ21A0734		(01)04026704400515(17)251228(10)KMZ21A0734
KMZ21B0533		(01)04026704400515(17)260128(10)KMZ21B0533
KMZ21B0181		(01)04026704400515(17)260128(10)KMZ21B0181
KMZ21B1171		(01)04026704400515(17)260128(10)KMZ21B1171
KMZ21C0388		(01)04026704400515(17)260228(10)KMZ21C0388
KMZ21C0883		(01)04026704400515(17)260228(10)KMZ21C0883
KMZ21C1160		(01)04026704400515(17)260228(10)KMZ21C1160
KMZ21G1038		(01)04026704400515(17)260628(10)KMZ21G1038
KMZ21G0751		(01)04026704400515(17)260628(10)KMZ21G0751
KMZ21G0971		(01)04026704400515(17)260628(10)KMZ21G0971
KMZ21K0266		(01)04026704400515(17)260928(10)KMZ21K0266
KMZ21J0446		(01)04026704400515(17)260828(10)KMZ21J0446
KMZ21J0534		(01)04026704400515(17)260828(10)KMZ21J0534
KMZ21L0541		(01)04026704400515(17)261028(10)KMZ21L0541
KMZ21L0220		(01)04026704400515(17)261028(10)KMZ21L0220
KMZ20K0946		(01)04026704400515(17)250928(10)KMZ20K0946
KMZ20K1260		(01)04026704400515(17)250928(10)KMZ20K1260
KMZ20K1294		(01)04026704400515(17)250928(10)KMZ20K1294
KMZ20K1310		(01)04026704400515(17)250928(10)KMZ20K1310
KMZ20L0724		(01)04026704400515(17)251028(10)KMZ20L0724
KMZ20L0537		(01)04026704400515(17)251028(10)KMZ20L0537
KMZ20L0706		(01)04026704400515(17)251028(10)KMZ20L0706
KMZ20L0886		(01)04026704400515(17)251028(10)KMZ20L0886
KMZ20L0945		(01)04026704400515(17)251028(10)KMZ20L0945
KMZ20L1091		(01)04026704400515(17)251028(10)KMZ20L1091
KMZ20L1103		(01)04026704400515(17)251028(10)KMZ20L1103
KMZ20L1104		(01)04026704400515(17)251028(10)KMZ20L1104
KMZ21A0693		(01)04026704400515(17)251228(10)KMZ21A0693
KMZ21A0618		(01)04026704400515(17)251228(10)KMZ21A0618
KMZ21A0635		(01)04026704400515(17)251228(10)KMZ21A0635
KMZ21A0706		(01)04026704400515(17)251228(10)KMZ21A0706
KMZ21A0714		(01)04026704400515(17)251228(10)KMZ21A0714
KMZ21A0725		(01)04026704400515(17)251228(10)KMZ21A0725
KMZ21B0101		(01)04026704400515(17)260128(10)KMZ21B0101
KMZ21B0163		(01)04026704400515(17)260128(10)KMZ21B0163
KMZ21B0689		(01)04026704400515(17)260128(10)KMZ21B0689
KMZ21C0270		(01)04026704400515(17)260228(10)KMZ21C0270
KMZ21C0398		(01)04026704400515(17)260228(10)KMZ21C0398
KMZ21C1159		(01)04026704400515(17)260228(10)KMZ21C1159
KMZ21F0777		(01)04026704400515(17)260528(10)KMZ21F0777

**Appendix 2 - Product Code / Lot Number / Unique Device Identifier (UDI)**

Product Code	Batch	UDI Number
19272T	KMZ21G0399	(01)04026704400515(17)260628(10)KMZ21G0399
	KMZ21G0403	(01)04026704400515(17)260628(10)KMZ21G0403
	KMZ21G0404	(01)04026704400515(17)260628(10)KMZ21G0404
	KMZ21G0786	(01)04026704400515(17)260628(10)KMZ21G0786
	KMZ21G0765	(01)04026704400515(17)260628(10)KMZ21G0765
	KMZ21G0838	(01)04026704400515(17)260628(10)KMZ21G0838
	KMZ21G0806	(01)04026704400515(17)260628(10)KMZ21G0806
	KMZ21G0911	(01)04026704400515(17)260628(10)KMZ21G0911
	KMZ21G0912	(01)04026704400515(17)260628(10)KMZ21G0912
	KMZ21G0917	(01)04026704400515(17)260628(10)KMZ21G0917
	KMZ21G0954	(01)04026704400515(17)260628(10)KMZ21G0954
	KMZ21G0955	(01)04026704400515(17)260628(10)KMZ21G0955
	KMZ21G0956	(01)04026704400515(17)260628(10)KMZ21G0956
	KMZ21G0964	(01)04026704400515(17)260628(10)KMZ21G0964
	KMZ21G0992	(01)04026704400515(17)260628(10)KMZ21G0992
	KMZ21G0982	(01)04026704400515(17)260628(10)KMZ21G0982
	KMZ21G0986	(01)04026704400515(17)260628(10)KMZ21G0986
	KMZ21G0993	(01)04026704400515(17)260628(10)KMZ21G0993
	KMZ21G0995	(01)04026704400515(17)260628(10)KMZ21G0995
	KMZ21G1001	(01)04026704400515(17)260628(10)KMZ21G1001
	KMZ21G1008	(01)04026704400515(17)260628(10)KMZ21G1008
	KMZ21G1017	(01)04026704400515(17)260628(10)KMZ21G1017
	KMZ21H0257	(01)04026704400515(17)260728(10)KMZ21H0257
	KMZ21H0299	(01)04026704400515(17)260728(10)KMZ21H0299
	KMZ21H0405	(01)04026704400515(17)260728(10)KMZ21H0405
	KMZ21J0447	(01)04026704400515(17)260828(10)KMZ21J0447
	KMZ21J0546	(01)04026704400515(17)260828(10)KMZ21J0546
	KMZ21J0554	(01)04026704400515(17)260828(10)KMZ21J0554
	KMZ21J0555	(01)04026704400515(17)260828(10)KMZ21J0555
	KMZ21K0164	(01)04026704400515(17)260928(10)KMZ21K0164
	KMZ21J0579	(01)04026704400515(17)260828(10)KMZ21J0579
	KMZ21J0580	(01)04026704400515(17)260828(10)KMZ21J0580
	KMZ21K0045	(01)04026704400515(17)260928(10)KMZ21K0045
	KMZ21K0114	(01)04026704400515(17)260928(10)KMZ21K0114
	KMZ21K0163	(01)04026704400515(17)260928(10)KMZ21K0163
	KMZ21K0308	(01)04026704400515(17)260928(10)KMZ21K0308
	KMZ21K0274	(01)04026704400515(17)260928(10)KMZ21K0274
	KMZ21K0295	(01)04026704400515(17)260928(10)KMZ21K0295
	KMZ21K0299	(01)04026704400515(17)260928(10)KMZ21K0299
	KMZ21K0314	(01)04026704400515(17)260928(10)KMZ21K0314
	KMZ21K0323	(01)04026704400515(17)260928(10)KMZ21K0323
	KMZ21K0318	(01)04026704400515(17)260928(10)KMZ21K0318
	KMZ21L0047	(01)04026704400515(17)261028(10)KMZ21L0047
	KMZ21M0781	(01)04026704400515(17)261128(10)KMZ21M0781
	KMZ22A0635	(01)04026704400515(17)261228(10)KMZ22A0635
	KMZ22B0130	(01)04026704400515(17)270128(10)KMZ22B0130
	KMZ22B0509	(01)04026704400515(17)270128(10)KMZ22B0509
	KMZ21L0074	(01)04026704400515(17)261028(10)KMZ21L0074
	KMZ21L0077	(01)04026704400515(17)261028(10)KMZ21L0077
	KMZ21L0079	(01)04026704400515(17)261028(10)KMZ21L0079
KMZ21L0219	(01)04026704400515(17)261028(10)KMZ21L0219	
KMZ21L0238	(01)04026704400515(17)261028(10)KMZ21L0238	
KMZ21L0239	(01)04026704400515(17)261028(10)KMZ21L0239	
KMZ21L0244	(01)04026704400515(17)261028(10)KMZ21L0244	
KMZ21L0295	(01)04026704400515(17)261028(10)KMZ21L0295	
KMZ21L0296	(01)04026704400515(17)261028(10)KMZ21L0296	
KMZ21L0307	(01)04026704400515(17)261028(10)KMZ21L0307	
KMZ21L0308	(01)04026704400515(17)261028(10)KMZ21L0308	

**Appendix 2 - Product Code / Lot Number / Unique Device Identifier (UDI)**

Product Code	Batch	UDI Number
19272T	KMZ21L0316	(01)04026704400515(17)261028(10)KMZ21L0316
	KMZ21L0333	(01)04026704400515(17)261028(10)KMZ21L0333
	KMZ21L0334	(01)04026704400515(17)261028(10)KMZ21L0334
	KMZ21L0337	(01)04026704400515(17)261028(10)KMZ21L0337
	KMZ21L0344	(01)04026704400515(17)261028(10)KMZ21L0344
	KMZ21L0367	(01)04026704400515(17)261028(10)KMZ21L0367
	KMZ21L0368	(01)04026704400515(17)261028(10)KMZ21L0368
	KMZ21L0382	(01)04026704400515(17)261028(10)KMZ21L0382
	KMZ21L0389	(01)04026704400515(17)261028(10)KMZ21L0389
	KMZ21L0878	(01)04026704400515(17)261028(10)KMZ21L0878
	KMZ21L0880	(01)04026704400515(17)261028(10)KMZ21L0880
	KMZ21L0901	(01)04026704400515(17)261028(10)KMZ21L0901
	KMZ21L0924	(01)04026704400515(17)261028(10)KMZ21L0924
	KMZ21M0050	(01)04026704400515(17)261128(10)KMZ21M0050
	KMZ21M0756	(01)04026704400515(17)261128(10)KMZ21M0756
	KMZ21M0786	(01)04026704400515(17)261128(10)KMZ21M0786
	KMZ21M0790	(01)04026704400515(17)261128(10)KMZ21M0790
	KMZ21M0805	(01)04026704400515(17)261128(10)KMZ21M0805
	KMZ22B0131	(01)04026704400515(17)270128(10)KMZ22B0131
	KMZ22B0192	(01)04026704400515(17)270128(10)KMZ22B0192
KMZ22B0524	(01)04026704400515(17)270128(10)KMZ22B0524	
KMZ21G0770	(01)04026704400515(17)260628(10)KMZ21G0770	

- end of document -