

Drägerwerk AG & Co. KGaA, 23542 Lübeck

To our customers of the
Dräger Anesthesia Workstations Atlan

July 2026

Second Addendum Important Safety Notice!

Possible failure of automatic ventilation
due to a fault in the motor electronics

**Incorrect manufacturing date string within UDIs provided with our
Atlan - Addendum Important Safety Notice!, dated April 2026**

Dear Sir or Madam

Due to a technical issue, the UDIs we provided with our last Field Safety Notice contain incorrect string section (11) for the manufacturing date.

UDI
(01)04048675556343(11)240109(21)ASSN-0311(93)8621300-15

Pic. 1 – incorrect string section (11)

To identify the affected devices properly, you may simply compare the serial number string section (21) with your serial numbers on the labelling of your devices (refer to Pic. 3).

UDI
(01)04048675556343(11)240109(21)ASSN-0311(93)8621300-15


Pic. 2 – serial number string section (21)

Drägerwerk AG & Co. KGaA
Moislinger Allee 53-55
23558 Lübeck, Germany
Postal address:
23542 Lübeck, Germany
Tel +49 451 882-0
Fax +49 451 882-2080
info@draeger.com
www.draeger.com
VAT no. DE135082211

Bank details:
Commerzbank AG, Lübeck
IBAN: DE95 2304 0022 0014 6795 00
Swift-Code: COBA DE FF 230
Sparkasse zu Lübeck
IBAN: DE15 2305 0101 0001 0711 17
Swift-Code: NOLADE21SPL

Registered office: Lübeck
Commercial register:
Local court Lübeck HRB 7903 HL
General partner: Drägerwerk Verwaltungs AG
Registered office: Lübeck
Commercial register:
Local court Lübeck HRB 7395 HL

Chairman of the Supervisory Board
for Drägerwerk AG & Co. KGaA
and Drägerwerk Verwaltungs AG:
Stefan Lauer
Executive Board:
Stefan Dräger (chairman)
Rainer Klug
Gert-Hartwig Lescow
Dr. Reiner Piske
Anton Schrofner



Pic. 3 – Device label with serial number [SN]

We apologize for the inconvenience resulting from this error.

Makoto Onda

Makoto Onda 2026-07-14

Global Product Manager, Atlan family

Business Unit Therapy, Anesthesiology