

RA2026-4417798

URGENT FIELD SAFETY NOTICE

MEDICAL DEVICE RECALL

August 2026

Product Field Action #: RA2026-4417798

Product Family Names: Orthopaedic implant components manufactured from UHMWPE (polyethylene), including legacy tibial inserts, patellas, acetabular liners/cups, shoulder and extremity implant components

Identification of the Affected Products: See **Appendix A: Part/Lot Number Attachment: RA2026-4417798**

Dear Customer,

Stryker has initiated a voluntary, lot specific recall for the devices listed in the **Appendix A: Part/Lot Number Attachment: RA2026-4417798**. The intent of this letter is to list all known hazards and harms potentially associated with the below noted issue and list the risk mitigation factors associated with the use of the product.

Issue:

Stryker has become aware that the devices listed in the **Appendix A: Part/Lot Number Attachment: RA2026-4417798** were manufactured with Ultra High Molecular Weight Polyethylene (UHMWPE) raw material aged in excess of five (5) years.

Industry guidance for storage duration of UHMWPE raw material used in the manufacture of implant bearing/articulating material of prosthetic joints is less than five (5) years. Product manufactured using UHMWPE raw material beyond the recommended storage duration has the potential for elevated levels of oxidation. Oxidation of UHMWPE material can impact material properties, potentially resulting in increased wear.

Potential Hazards:

- Excessive Wear Debris
- Material Fragments
- Fractured Device
- Disassociated Device

Potential Harms:

- Revision Surgery
- Pain
- Inflammation

Risk Mitigations:

Product with a high oxidation index can become discolored. However, not all product with a high oxidation index becomes discolored. In instances where the product is discolored, the issue may be recognized by the user. Where present, this visual observation may assist users in identifying potentially impacted devices prior to use. Based on medical assessment, the occurrence of the potential harms described above is considered low and unlikely. Post-market surveillance and complaint data have not identified adverse trends associated with the devices in the scope of this issue. Review of complaint and

vigilance data has not identified evidence suggesting elevated oxidation levels are associated with an increased clinical risk in the affected device population.

Recommendations for patients already implanted with an impacted device:

Patients treated with an impacted product should continue to be followed per the normal protocol established by their surgeon(s). There are no recommended changes to the frequency of the standard follow-up care protocol.

Actions Needed:

Our records indicate that you may have received the affected product(s). It is Stryker's responsibility as the manufacturer to ensure that customers who may have received these affected products also receive this important communication. We therefore request that you read this notice carefully and complete the following actions.

1. Please inform healthcare practitioners who have used devices in scope of this Urgent Medical Device Recall listed in the **Appendix A: Part/Lot Number Attachment: RA2026-4417798** and forward this notice to all individuals who need to be made aware.
2. Sign and return the enclosed Business Reply Form by email to **<XXXX@stryker.com>** to confirm receipt of this notification/documenting product segregation.
 - a. Response is required, even if you may not have any physical inventory on site anymore. It may be that you no longer have any physical inventory on site. Completing this form will allow us to update our records and will also eliminate the need for us to send further unnecessary communication on this matter. Therefore, please complete the form even if you no longer have any of the subject devices in your physical inventory.
3. Upon receipt of the complete Business Reply Form, Stryker will contact you to arrange for the return of your affected product(s). A replacement will be provided upon receipt of the recalled product.
4. Maintain awareness of this communication internally until all required actions have been completed within your facility.
5. If you have further distributed the affected product, please notify the applicable parties at once about this recall. You may copy and distribute this notification letter.
 - a) If possible, inform us if any of the subject devices have been distributed to other organizations, including contact details.
 - b) If you are a distributor, note that you are responsible for notifying your affected customers

Please inform us of any adverse events and/or report them to the Health/Competent Authorities in accordance with current regulations.

Your designated contact person for this action is given below. Should you have any queries concerning this matter please do not hesitate to contact them directly.

Name:	Position:	Email:
-------	-----------	--------

In line with the recommendations of the Medical Device Coordination Group Guidance document Ref MDCG 2023-3 and EU MDR 2017/745, we can confirm that this FSCA has been notified appropriately to the National Competent Authority for your country.

On behalf of Stryker, we thank you sincerely for your help and support in completing this action and regret any inconvenience caused. We would like to reassure you that Stryker is committed to ensuring that only conforming devices, meeting our high internal quality standards and your expectations, remain on the market.

Sincerely,

RA2026-4417798**Business Reply Form- response required****MEDICAL DEVICE RECALL****August 2026****Product Field Action #:** RA2026-4417798

Product Names: Orthopaedic implant components manufactured from UHMWPE barstock, including knee system components (e.g., tibial inserts), acetabular liners/cups, shoulder and extremity implant components

I have received the **Urgent Medical Device Recall** letter from Stryker dated **August XX, 2026**, stating that the company has initiated a voluntary recall of the above-referenced affected products.

We have not located any of these devices in our inventory <i>(please add check mark to box):</i>			
We have located and returned the following devices (Document total # returned):			
Product	Catalog Number	Lot Number	Total # Returned:

Hospital/Surgeon/Branch Name_____
Date_____
Hospital/Surgeon/Branch Address_____
Hospital/Surgeon/Branch Rep
(Signature)

**PLEASE COMPLETE THIS FORM WITHIN XX BUSINESS DAYS AND RETURN IT
BY USING THE EMAIL LISTED BELOW: xxx@stryker.com**

Appendix A

Part/Lot Number Attachment: RA2026-4417798			
Part Number	Lot Number	Description	UDI DI
2001-2246	13MFT	OMNIFIT SERIES I ACET. INSERT 10° 22mm ID	07613327037883
2099-2250	1K5VWW	CONSTRAINED ACETABULAR INSERTS 22MM 50MM	07613327031430
2099-2250	2367AK	CONSTRAINED ACETABULAR INSERTS 22MM 50MM	07613327031430
2099-2250	3X0NJD	CONSTRAINED ACETABULAR INSERTS 22MM 50MM	07613327031430
2099-2250	4W4MPL	CONSTRAINED ACETABULAR INSERTS 22MM 50MM	07613327031430
2099-2250	LV80NE	CONSTRAINED ACETABULAR INSERTS 22MM 50MM	07613327031430
2099-2252	NV561W	CONSTRAINED ACETABULAR INSERTS 22MM 52MM	07613327031447
2099-2252	PM326V	CONSTRAINED ACETABULAR INSERTS 22MM 52MM	07613327031447
2099-2256	3H5RRK	CONSTRAINED ACETABULAR INSERTS 22MM 56MM	07613327031454
2099-2256	AA7JK0	CONSTRAINED ACETABULAR INSERTS 22MM 56MM	07613327031454
2099-2861	MLLV6	CONSTRAINED ACETABULAR INSERTS 28MM 61MM	07613327031591
3044-0026	14SCVA	SCORPIO RECESSED PATELLA 26MM	07613327033779
3044-0026	14SCVD	SCORPIO RECESSED PATELLA 26MM	07613327033779
3044-0026	14SCVF	SCORPIO RECESSED PATELLA 26MM	07613327033779
3045-0026	12WAHA	SCORPIO RECESSED METAL BACKED PATELLA	N/A
3045-0026	12WAHB	SCORPIO RECESSED METAL BACKED PATELLA	N/A
3045-0026	HOHE	SCORPIO RECESSED METAL BACKED PATELLA	N/A
3045-0026	HOVX	SCORPIO RECESSED METAL BACKED PATELLA	N/A
3045-0028	13VAW	SCORPIO RECESSED METAL BACKED PATELLA	N/A
3045-0028	13VAWA	SCORPIO RECESSED METAL BACKED PATELLA	N/A
3045-0030	BIYIE	SCORPIO RECESSED METAL BACKED PATELLA	N/A
3045-0030	HNZY	SCORPIO RECESSED METAL BACKED PATELLA	N/A
3045-0032	BIYIH	SCORPIO RECESSED METAL BACKED PATELLA	N/A
3045-0032	HOYK	SCORPIO RECESSED METAL BACKED PATELLA	N/A
3052-1115	MJ7YY6	SERIES II TIBIAL BEAR INSERT NO 11 15MM	07613327034240
5361-6105	AF2S09U	SHOULDER - GLENOID #5	07613327138214
5361-6105	AF2S09Z	SHOULDER - GLENOID #5	07613327138214
5361-6105	P1H41	SHOULDER - GLENOID #5	07613327138214
5361-6105	P4T02	SHOULDER - GLENOID #5	07613327138214
5361-6105	P4T02X1	SHOULDER - GLENOID #5	07613327138214
5361-6105	P4T02X2	SHOULDER - GLENOID #5	07613327138214
5361-6105	P4T02X3	SHOULDER - GLENOID #5	07613327138214
5361-6105	P4T02X4	SHOULDER - GLENOID #5	07613327138214
5361-6105	P4T02X5	SHOULDER - GLENOID #5	07613327138214
5361-6105	P4T02X6	SHOULDER - GLENOID #5	07613327138214
5361-6107	AF3E00	SHOULDER - GLENOID #7	07613327138207
5361-6107	AF3E00B	SHOULDER - GLENOID #7	07613327138207
5361-6107	AF3E00F	SHOULDER - GLENOID #7	07613327138207

5361-6107	AF3E00G	SHOULDER - GLENOID #7	07613327138207
5361-6107	AF3E00Q	SHOULDER - GLENOID #7	07613327138207
5361-6107	AF3E00S	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4A238	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4A238X1	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4A238X2	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4A238X3	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4A238X4	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4A238X5	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4A238X6	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4A238X7	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4A238X8	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4L09	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4L09X1	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4L09X10	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4L09X2	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4L09X3	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4L09X4	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4L09X5	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4L09X6	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4L09X7	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4L09X8	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4L09X9	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4N06	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4N06X1	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4N06X2	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4N06X3	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4N06X4	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4N06X5	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4N06X6	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4N06X7	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4T03	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4T03X1	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4T03X2	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4T03X3	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4T03X4	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4T03X5	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4T03X6	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4T03X7	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4T03X8	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4V01	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4V01X1	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4V01X2	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4V01X3	SHOULDER - GLENOID #7	07613327138207

5361-6107	P4V01X4	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4V01X5	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4V01X6	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4V01X7	SHOULDER - GLENOID #7	07613327138207
5361-6107	P4V01X8	SHOULDER - GLENOID #7	07613327138207
5361-6109	AF3E02I	SHOULDER - GLENOID #9	07613327138153
5361-6109	AF3E02S	SHOULDER - GLENOID #9	07613327138153
5361-6109	AF3E02Z	SHOULDER - GLENOID #9	07613327138153
5361-6109	P1H39Q	SHOULDER - GLENOID #9	07613327138153
5361-6111	AFZT00U	SHOULDER - GLENOID #11	07613327138191
5361-6111	AFZT00Z	SHOULDER - GLENOID #11	07613327138191
5361-6111	P2N04B	SHOULDER - GLENOID #11	07613327138191
5361-6111	P3K57	SHOULDER - GLENOID #11	07613327138191
5361-6111	P3K57B	SHOULDER - GLENOID #11	07613327138191
5361-6111	P3K57F	SHOULDER - GLENOID #11	07613327138191
5361-6111	PYE43B	SHOULDER - GLENOID #11	07613327138191
5361-6111	PZA70F	SHOULDER - GLENOID #11	07613327138191
5361-6111	PZA70U	SHOULDER - GLENOID #11	07613327138191
5361-6205	1NMMPE	SHOULDER - KEELED GLENOID #5	07613327138177
5361-6205	1NNMPE	SHOULDER - KEELED GLENOID #5	07613327138177
5361-6205	P3K58	SHOULDER - KEELED GLENOID #5	07613327138177
5361-6205	P3K58B	SHOULDER - KEELED GLENOID #5	07613327138177
5361-6205	P3K58F	SHOULDER - KEELED GLENOID #5	07613327138177
5361-6207	AF2S10B	SHOULDER - KEELED GLENOID #7	07613327138269
5361-6207	AF2S10F	SHOULDER - KEELED GLENOID #7	07613327138269
5361-6207	AF2S10G	SHOULDER - KEELED GLENOID #7	07613327138269
5361-6207	P4A241	SHOULDER - KEELED GLENOID #7	07613327138269
5361-6207	P4A241X1	SHOULDER - KEELED GLENOID #7	07613327138269
5361-6207	P4A241X2	SHOULDER - KEELED GLENOID #7	07613327138269
5361-6207	P4A241X3	SHOULDER - KEELED GLENOID #7	07613327138269
5361-6207	P4V03	SHOULDER - KEELED GLENOID #7	07613327138269
5361-6207	P4V03X1	SHOULDER - KEELED GLENOID #7	07613327138269
5361-6207	P4V03X2	SHOULDER - KEELED GLENOID #7	07613327138269
5361-6207	P4V03X3	SHOULDER - KEELED GLENOID #7	07613327138269
5361-6207	P4V03X4	SHOULDER - KEELED GLENOID #7	07613327138269
5361-6207	P4V03X5	SHOULDER - KEELED GLENOID #7	07613327138269
5361-6207	PZKDFX1	SHOULDER - KEELED GLENOID #7	07613327138269
5361-6209	AF2K03Z	SHOULDER - KEELED GLENOID #9	07613327138238
5361-6209	AFZS10S	SHOULDER - KEELED GLENOID #9	07613327138238
5361-6209	AFZS10U	SHOULDER - KEELED GLENOID #9	07613327138238
5361-6211	AFZS09F	SHOULDER - KEELED GLENOID #11	07613327138283
5361-6211	AFZS09G	SHOULDER - KEELED GLENOID #11	07613327138283
5361-6211	AFZS09S	SHOULDER - KEELED GLENOID #11	07613327138283
5361-6211	P2N08	SHOULDER - KEELED GLENOID #11	07613327138283

5361-6211	PYC30X1	SHOULDER - KEELED GLENOID #11	07613327138283
5361-6211	PYW24X1	SHOULDER - KEELED GLENOID #11	07613327138283
5361-6211	PZKDX1	SHOULDER - KEELED GLENOID #11	07613327138283
5612-P-411	360958	TRIATHLON HINGE INSERT SIZE 4 11MM	07613327601893
5612-P-413	354573	TRIATHLON HINGE INSERT SIZE 4 13MM	07613327602012
6074-4022	2M5NTY	SYSTEM ONE BIPOLAR 40MM	07613327031775
6074-4022	LH25R5	SYSTEM ONE BIPOLAR 40MM	07613327031775
6074-4022	M36N83	SYSTEM ONE BIPOLAR 40MM	07613327031775
6074-4322	W40E3T	SYSTEM ONE BIPOLAR 43MM	07613327031805
6074-4422	R84K20	SYSTEM ONE BIPOLAR 44MM	07613327031812
61-3258	191723	ACETABULAR CUP UHMWPE	07613327038538
61-3258	191867	ACETABULAR CUP UHMWPE	07613327038538
620-00-28D	MET76R	TRIDENT 0° POLY INSERT 28mm ID	07613327313581
620-00-28D	MHANP2	TRIDENT 0° POLY INSERT 28mm ID	07613327313581
620-00-28D	MHD1N0	TRIDENT 0° POLY INSERT 28mm ID	07613327313581
6299-N-332	19ECK	PCA CONST LINER 32MM NEUTRAL C - SHELL	07613327138870
6302-H-632	15NKT	SYS 12 CONST LINER 32ID 10 DEG P6 +4MM	07613327138979
6636-0-419	13SFLA	EIUS UNI TIB SMALL RM/LL 9MM	N/A
690-00-22D	123J5R	TRIDENT 0 DEG CONSTRAINED INSERT 22D	07613327025538
690-00-22D	194RN0	TRIDENT 0 DEG CONSTRAINED INSERT 22D	07613327025538
690-00-22D	274HMK	TRIDENT 0 DEG CONSTRAINED INSERT 22D	07613327025538
690-00-22D	3M77NN	TRIDENT 0 DEG CONSTRAINED INSERT 22D	07613327025538
690-00-22D	5J8P89	TRIDENT 0 DEG CONSTRAINED INSERT 22D	07613327025538
690-00-22D	5L82T0	TRIDENT 0 DEG CONSTRAINED INSERT 22D	07613327025538
690-00-22D	986HVM	TRIDENT 0 DEG CONSTRAINED INSERT 22D	07613327025538
690-00-22D	A23P1X	TRIDENT 0 DEG CONSTRAINED INSERT 22D	07613327025538
690-00-22D	E923WW	TRIDENT 0 DEG CONSTRAINED INSERT 22D	07613327025538
690-00-22D	MLMEW8	TRIDENT 0 DEG CONSTRAINED INSERT 22D	07613327025538
690-00-22D	MLMTVP	TRIDENT 0 DEG CONSTRAINED INSERT 22D	07613327025538
690-00-22D	MLNPX0	TRIDENT 0 DEG CONSTRAINED INSERT 22D	07613327025538
690-00-22D	VX6813	TRIDENT 0 DEG CONSTRAINED INSERT 22D	07613327025538
690-00-22E	4100R9	TRIDENT 0 DEG CONSTRAINED INSERT 22E	07613327025545
690-00-22E	5V78KA	TRIDENT 0 DEG CONSTRAINED INSERT 22E	07613327025545
690-00-22E	711K4K	TRIDENT 0 DEG CONSTRAINED INSERT 22E	07613327025545
690-00-22E	EY7HL9	TRIDENT 0 DEG CONSTRAINED INSERT 22E	07613327025545
690-00-22E	KP576W	TRIDENT 0 DEG CONSTRAINED INSERT 22E	07613327025545
690-00-22E	L620XW	TRIDENT 0 DEG CONSTRAINED INSERT 22E	07613327025545
690-00-22E	MD8JH9	TRIDENT 0 DEG CONSTRAINED INSERT 22E	07613327025545
690-00-22E	T61W43	TRIDENT 0 DEG CONSTRAINED INSERT 22E	07613327025545
690-00-22E	XA0L27	TRIDENT 0 DEG CONSTRAINED INSERT 22E	07613327025545
690-00-28F	284X4X	TRIDENT 0 DEG CONSTRAINED INSERT 28F	07613327025552
690-00-28F	2K2DW1	TRIDENT 0 DEG CONSTRAINED INSERT 28F	07613327025552
690-00-28F	E47WVY	TRIDENT 0 DEG CONSTRAINED INSERT 28F	07613327025552
690-00-28F	MLL8JJ	TRIDENT 0 DEG CONSTRAINED INSERT 28F	07613327025552

690-00-28F	MLLKAY	TRIDENT 0 DEG CONSTRAINED INSERT 28F	07613327025552
690-00-28F	MLLP6K	TRIDENT 0 DEG CONSTRAINED INSERT 28F	07613327025552
690-00-28F	MLLYT7	TRIDENT 0 DEG CONSTRAINED INSERT 28F	07613327025552
690-00-28F	YE6959	TRIDENT 0 DEG CONSTRAINED INSERT 28F	07613327025552
690-00-28G	6D8274	TRIDENT 0 DEG CONSTRAINED INSERT 28G	07613327025583
690-00-28G	9N8EWW	TRIDENT 0 DEG CONSTRAINED INSERT 28G	07613327025583
690-00-32I	3A2P68	TRIDENT 0 DEG CONSTRAINED INSERT 32I	07613327025590
690-00-32I	4A40TV	TRIDENT 0 DEG CONSTRAINED INSERT 32I	07613327025590
690-10-22E	29385D	TRIDENT CONSTRAINED INSERT	07613327024289
690-10-22E	906WEN	TRIDENT CONSTRAINED INSERT	07613327024289
690-10-22E	DJ7TP1	TRIDENT CONSTRAINED INSERT	07613327024289
690-10-22E	EK4TK5	TRIDENT CONSTRAINED INSERT	07613327024289
690-10-22E	LA22KA	TRIDENT CONSTRAINED INSERT	07613327024289
690-10-22E	WL248Y	TRIDENT CONSTRAINED INSERT	07613327024289
690-10-22E	XV01J8	TRIDENT CONSTRAINED INSERT	07613327024289
690-10-22F	A6592H	TRIDENT CONSTRAINED INSERT	07613327024357
690-10-22F	L7140Y	TRIDENT CONSTRAINED INSERT	07613327024357
690-10-22F	P103P0	TRIDENT CONSTRAINED INSERT	07613327024357
690-10-22F	W745P8	TRIDENT CONSTRAINED INSERT	07613327024357
690-10-28G	MLL4Y7	TRIDENT CONSTRAINED INSERT	07613327024326
690-10-28H	MLKNT7	TRIDENT CONSTRAINED INSERT	07613327024319
690-10-28H	MLL1YN	TRIDENT CONSTRAINED INSERT	07613327024319
690-10-28H	MLLYXD	TRIDENT CONSTRAINED INSERT	07613327024319
690-10-28I	MLL8L9	TRIDENT CONSTRAINED INSERT	07613327024364
690-10-28I	MLLXWY	TRIDENT CONSTRAINED INSERT	07613327024364
690-10-28K	MHNRM4	TRIDENT CONSTR ACET INSERT	07613327024449
690-10-28K	MKRH1M	TRIDENT CONSTR ACET INSERT	07613327024449
690-10-28K	MLAAV7	TRIDENT CONSTR ACET INSERT	07613327024449
69-2244	990V1W	ALL POLY CONSTRAINED INSERT 44MM	07613327025446
69-2244	D56JV4	ALL POLY CONSTRAINED INSERT 44MM	07613327025446
69-2244	D74X5E	ALL POLY CONSTRAINED INSERT 44MM	07613327025446
69-2244	DX56K2	ALL POLY CONSTRAINED INSERT 44MM	07613327025446
69-2244	E67D6N	ALL POLY CONSTRAINED INSERT 44MM	07613327025446
69-2244	EP10WL	ALL POLY CONSTRAINED INSERT 44MM	07613327025446
69-2244	KP14HN	ALL POLY CONSTRAINED INSERT 44MM	07613327025446
69-2244	NN8NM2	ALL POLY CONSTRAINED INSERT 44MM	07613327025446
69-2244	P36TXM	ALL POLY CONSTRAINED INSERT 44MM	07613327025446
69-2244	P95L8X	ALL POLY CONSTRAINED INSERT 44MM	07613327025446
69-2244	RE4EK0	ALL POLY CONSTRAINED INSERT 44MM	07613327025446
69-2244	X82K0D	ALL POLY CONSTRAINED INSERT 44MM	07613327025446
69-2246	8X6EL1	ALL POLY CONSTRAINED INSERT 46MM	07613327025453
69-2248	9J3DPH	ALL POLY CONSTRAINED INSERT 48MM	07613327025460
69-2248	A41E7D	ALL POLY CONSTRAINED INSERT 48MM	07613327025460
69-2248	AD39A3	ALL POLY CONSTRAINED INSERT 48MM	07613327025460

69-2248	RP673W	ALL POLY CONSTRAINED INSERT 48MM	07613327025460
69-2850	MLL1M3	ALL POLY CONSTRAINED INSERT 50MM	07613327025477
69-2850	MLLXT1	ALL POLY CONSTRAINED INSERT 50MM	07613327025477
69-2850	MLLY6H	ALL POLY CONSTRAINED INSERT 50MM	07613327025477
69-2852	MLL0NM	ALL POLY CONSTRAINED INSERT 52MM	07613327025484
69-2852	MLLAHY	ALL POLY CONSTRAINED INSERT 52MM	07613327025484
69-2854	MLKY8A	ALL POLY CONSTRAINED INSERT 54MM	07613327025491
69-2854	MLL8E1	ALL POLY CONSTRAINED INSERT 54MM	07613327025491
72-4-9114	8H39KP	SCORPIO TS TIB INS	07613327034875
73-0110	HHPW	SCORPIO M-DOME PATELLA	07613327050288
73-0710	8P83	SCORPIO M-DOME PATELLA	07613327050295
73-0710	DL4V	SCORPIO M-DOME PATELLA	07613327050295
73-0910	14NACD	SCORPIO M-DOME PATELLA	07613327050332
73-2610	193307A	SCORPIO INSET PATELLA	07613327014082
73-2610	193307E	SCORPIO INSET PATELLA	07613327014082
73-2710	14HDSA	SCORPIO C-DOME PATELLA	07613327033847
73-2710	NPPV	SCORPIO C-DOME PATELLA	07613327033847
73-2910	14AGFD	SCORPIO C-DOME PATELLA	07613327033854
73-3110	1L45	SCORPIO U-DOME PATELLA	07613327033861
73-3110	4JM5	SCORPIO U-DOME PATELLA	07613327033861
73-3110	7EL5	SCORPIO U-DOME PATELLA	07613327033861
73-3110	AME3	SCORPIO U-DOME PATELLA	07613327033861
73-3110	MK74	SCORPIO U-DOME PATELLA	07613327033861
73-3110	NTV7	SCORPIO U-DOME PATELLA	07613327033861
73-3110	PYHM	SCORPIO U-DOME PATELLA	07613327033861
73-3110	RE7J	SCORPIO U-DOME PATELLA	07613327033861
73-3110	VEWJ	SCORPIO U-DOME PATELLA	07613327033861
73-3508	HMM5	SCORPIO U-DOME PATELLA	07613327033885
73-3508	NLRE	SCORPIO U-DOME PATELLA	07613327033885
73-3708	3467	SCORPIO U-DOME PATELLA	07613327033892
73-3708	3D8D	SCORPIO U-DOME PATELLA	07613327033892
73-3708	PMRHX1	SCORPIO U-DOME PATELLA	07613327033892
73-3708	T6YJ	SCORPIO U-DOME PATELLA	07613327033892
73-3708	T9ME	SCORPIO U-DOME PATELLA	07613327033892
73-3708	VWDA	SCORPIO U-DOME PATELLA	07613327033892
73-3708	YAM7	SCORPIO U-DOME PATELLA	07613327033892
73-3710	3W25	SCORPIO U-DOME PATELLA	07613327033908
73-3710	DLYP	SCORPIO U-DOME PATELLA	07613327033908
73-3910	14NEYA	SCORPIO U-DOME PATELLA	07613327033915
73-3910	1XED	SCORPIO U-DOME PATELLA	07613327033915
73-3910	9Y3T	SCORPIO U-DOME PATELLA	07613327033915
73-3910	D1M9	SCORPIO U-DOME PATELLA	07613327033915
73-3910	D2YP	SCORPIO U-DOME PATELLA	07613327033915
73-3910	EX8T	SCORPIO U-DOME PATELLA	07613327033915

73-3910	KRT3	SCORPIO U-DOME PATELLA	07613327033915
73-3910	KRT3X1	SCORPIO U-DOME PATELLA	07613327033915
73-3910	T325	SCORPIO U-DOME PATELLA	07613327033915
73-3910	V8X8	SCORPIO U-DOME PATELLA	07613327033915
73-3910	X6YE	SCORPIO U-DOME PATELLA	07613327033915
C-M062-6-800	RX06268A	SELF-PRESSURIZING GLENOID SIZE 8	N/A
C-M062-7-500	RX06275A	SELF-PRESSURIZING GLENOID SIZE 8	N/A
C-M063-2-600	RX06326A	SELF-PRESSURIZING GLENOID SIZE 8	N/A
C-M063-5-000	RX06350A	SELF-PRESSURIZING GLENOID SIZE 9	N/A
C-M063-5-000	RX06350AA	SELF-PRESSURIZING GLENOID SIZE 9	N/A
C-M063-5-800	RX06358A	SELF-PRESSURIZING GLENOID SIZE 9	N/A
C-M064-1-300	RX06413A	DURACON XL ANGLED TIB INS 11/14MM	N/A
C-M064-4-300	RX06443A	SELF-PRESSURIZING GLENOID SIZE 9	N/A
C-M066-0-901	RX06609B	TRIDENT ALL POLY CUP 44 X 66MM	N/A
C-M066-3-000	RX06630A	DURACON TIB INSERT LG 12MM MED 14MM LAT	N/A
C-M067-3-103	MHP0AX	MRH KNEE TIBIAL INSERT FOR THIN POST	N/A
C-M067-5-002	RX06750C	TRIATHLON NO. 9 TIBIAL INSERT 9MM	N/A
C-M067-5-003	RX06750D	TRIATHLON NO. 9 TIBIAL INSERT 11MM	N/A
C-M071-6-601	G3055684	KALI TIBIAL INSERT LARGE RIGHT 8mm	N/A
C-M071-6-700	G3043883	ARC2F ACETABULAR INSERT 48/50/52mm OD x 28mm ID	N/A
C-M071-8-400	G3043878	ARC2F ACETABULAR INSERT 50mm OD x 28mm ID 10DEGREE HOOD	N/A
C-M079-2-801	G3264034	ABG I ACETABULAR INSERT HOODED 28mm ID X 48mm	N/A
C-M082-7-303	L6237	EXETER WINGLESS CENTRALISER FOR C-M082-7-302	N/A
C-M104-2-002	192371	TRIDENT CONSTRAINED ACETABULAR INSERT SIZE A X 22mm ID UH1	N/A
C-M104-2-003	192370	TRIDENT CONSTRAINED ACETABULAR INSERT SIZE B X 22mm ID UH1	N/A
C-M107-5-301	193607	EXETER WINGLESS CENTRALIZER XSM	N/A
C-M107-7-801	193479	MRS RADIAL BUSHING FORITH C-M107-7-800	N/A
C-M110-1-201	573P37	BEARING AND AXLE BUSHING ASSY STANDARD ELBOW	07613327504811
C-M112-2-500	P22E159	AXLE PIN ASSY LARGE ELBOW COBALT CHROME	N/A
C-M113-1-301	MH621A	BEARING AND AXLE BUSHING ASSY STD ELBOW	07613327590074
C-M113-5-600	194593	KRH REV TIB INS 16mm THK .640 DIA STEM	N/A
C-M113-5-603	194594	KRH REV TIB INS 16mm THK .655 DIA STEM	N/A
C-MM100253-02	N8873H	BUSHING AND BEARING USE C-MM100253-00	N/A
C-MM100411-02	194961	GMRS SML TIBIA 8MM WITH CASE 100411	N/A
C-MM100411-03	194962	GMRS SML TIBIA 11MM WITH CASE 100411	N/A
C-MM100411-04	194963	GMRS SML TIBIA 14MM WITH CASE 100411	N/A
C-MM100507-04	LNH795	GMRS TIB INSRT SML 10MM WITH CASE 100507	N/A
C-MM100507-06	LNH797	GMRS TIB INSRT SML 16MM WITH CASE 100507	N/A
C-MM100520-02	195080	GMRS SML POLY TIBIA 8MM CASE 100520	N/A
C-MM100520-03	195081	GMRS SML POLY TIBIA 11MM CASE 100520	N/A
C-MM100520-04	195082	GMRS SML POLY TIBIA 16MM CASE 100520	N/A
C-MM100616-01	195260	GMRS EXTRA SM ALL-POLY TIBIA SIZE 8mm	N/A
C-MM100616-02	195261	GMRS EXTRA SM ALL-POLY TIBIA SIZE 11mm	N/A

C-MM100616-03	195262	GMRS EXTRA SM ALL-POLY TIBIA SIZE 16mm	N/A
UH1-38-22	2A01TA	UHR BIPOLAR 22X38MM	07613327016789
UH1-38-22	430715	UHR BIPOLAR 22X38MM	07613327016789
UH1-38-22	49510K	UHR BIPOLAR 22X38MM	07613327016789
UH1-38-22	7E5X23	UHR BIPOLAR 22X38MM	07613327016789
UH1-38-22	852J98	UHR BIPOLAR 22X38MM	07613327016789
UH1-38-22	D23664	UHR BIPOLAR 22X38MM	07613327016789
UH1-38-22	NE659N	UHR BIPOLAR 22X38MM	07613327016789
UH1-40-22	128LXW	UHR BIPOLAR 22X40MM	07613327016796
UH1-40-22	16182H	UHR BIPOLAR 22X40MM	07613327016796
UH1-40-22	1J5955	UHR BIPOLAR 22X40MM	07613327016796
UH1-40-22	1T8203	UHR BIPOLAR 22X40MM	07613327016796
UH1-40-22	2137N4	UHR BIPOLAR 22X40MM	07613327016796
UH1-40-22	2541YM	UHR BIPOLAR 22X40MM	07613327016796
UH1-40-22	2Y52Y2	UHR BIPOLAR 22X40MM	07613327016796
UH1-40-22	4E19VD	UHR BIPOLAR 22X40MM	07613327016796
UH1-40-22	4K7XYL	UHR BIPOLAR 22X40MM	07613327016796
UH1-40-22	577M4L	UHR BIPOLAR 22X40MM	07613327016796
UH1-40-22	5D3AVR	UHR BIPOLAR 22X40MM	07613327016796
UH1-40-22	746AV4	UHR BIPOLAR 22X40MM	07613327016796
UH1-40-22	7V8J0V	UHR BIPOLAR 22X40MM	07613327016796
UH1-40-22	828TAR	UHR BIPOLAR 22X40MM	07613327016796
UH1-40-22	875M72	UHR BIPOLAR 22X40MM	07613327016796
UH1-40-22	8V63LX	UHR BIPOLAR 22X40MM	07613327016796
UH1-40-22	9A3W5N	UHR BIPOLAR 22X40MM	07613327016796
UH1-40-22	A91P07	UHR BIPOLAR 22X40MM	07613327016796
UH1-40-22	DL4D38	UHR BIPOLAR 22X40MM	07613327016796
UH1-40-22	LW709R	UHR BIPOLAR 22X40MM	07613327016796
UH1-40-22	RA6A1H	UHR BIPOLAR 22X40MM	07613327016796
UH1-40-22	VT4K59	UHR BIPOLAR 22X40MM	07613327016796
UH1-40-22	X95RTV	UHR BIPOLAR 22X40MM	07613327016796
UH1-41-22	343653	UHR BIPOLAR 22X41MM	07613327016819
UH1-41-22	4P1TR7	UHR BIPOLAR 22X41MM	07613327016819
UH1-41-22	521MLK	UHR BIPOLAR 22X41MM	07613327016819
UH1-41-22	7Y8HEA	UHR BIPOLAR 22X41MM	07613327016819
UH1-41-22	975A6Y	UHR BIPOLAR 22X41MM	07613327016819
UH1-41-22	AR0KX5	UHR BIPOLAR 22X41MM	07613327016819
UH1-41-22	E788E6	UHR BIPOLAR 22X41MM	07613327016819
UH1-41-22	KL1YTE	UHR BIPOLAR 22X41MM	07613327016819
UH1-41-22	KN0086	UHR BIPOLAR 22X41MM	07613327016819
UH1-41-22	P043ME	UHR BIPOLAR 22X41MM	07613327016819
UH1-41-22	PL4L5M	UHR BIPOLAR 22X41MM	07613327016819
UH1-41-26	157V1T	UHR BIPOLAR 26X41mm	07613327016802
UH1-41-26	1D33YX	UHR BIPOLAR 26X41mm	07613327016802

UH1-41-26	1K1E6M	UHR BIPOLAR 26X41mm	07613327016802
UH1-41-26	1R5LWT	UHR BIPOLAR 26X41mm	07613327016802
UH1-41-26	218RED	UHR BIPOLAR 26X41mm	07613327016802
UH1-41-26	254MT3	UHR BIPOLAR 26X41mm	07613327016802
UH1-41-26	287D4A	UHR BIPOLAR 26X41mm	07613327016802
UH1-41-26	2K7HXW	UHR BIPOLAR 26X41mm	07613327016802
UH1-41-26	2X5MTA	UHR BIPOLAR 26X41mm	07613327016802
UH1-41-26	3J1A39	UHR BIPOLAR 26X41mm	07613327016802
UH1-41-26	3K3P0T	UHR BIPOLAR 26X41mm	07613327016802
UH1-41-26	473K4V	UHR BIPOLAR 26X41mm	07613327016802
UH1-41-26	5A1T5N	UHR BIPOLAR 26X41mm	07613327016802
UH1-41-26	5D6KJ5	UHR BIPOLAR 26X41mm	07613327016802
UH1-41-26	5V7H82	UHR BIPOLAR 26X41mm	07613327016802
UH1-41-26	6J8XKT	UHR BIPOLAR 26X41mm	07613327016802
UH1-41-26	6W1MKA	UHR BIPOLAR 26X41mm	07613327016802
UH1-41-26	8R29TH	UHR BIPOLAR 26X41mm	07613327016802
UH1-41-26	9T7N8Y	UHR BIPOLAR 26X41mm	07613327016802
UH1-41-26	AL5K8M	UHR BIPOLAR 26X41mm	07613327016802
UH1-41-26	AW02V8	UHR BIPOLAR 26X41mm	07613327016802
UH1-41-26	DT8JEW	UHR BIPOLAR 26X41mm	07613327016802
UH1-41-26	H04PHL	UHR BIPOLAR 26X41mm	07613327016802
UH1-41-26	JR6R9L	UHR BIPOLAR 26X41mm	07613327016802
UH1-41-26	LN70PY	UHR BIPOLAR 26X41mm	07613327016802
UH1-41-26	MA6M6E	UHR BIPOLAR 26X41mm	07613327016802
UH1-41-26	PE7T6P	UHR BIPOLAR 26X41mm	07613327016802
UH1-41-26	PY8PLH	UHR BIPOLAR 26X41mm	07613327016802
UH1-41-26	RY49EW	UHR BIPOLAR 26X41mm	07613327016802
UH1-41-26	T1497X	UHR BIPOLAR 26X41mm	07613327016802
UH1-41-26	V70PNP	UHR BIPOLAR 26X41mm	07613327016802
UH1-41-26	VJ6WXM	UHR BIPOLAR 26X41mm	07613327016802
UH1-41-26	X82M1D	UHR BIPOLAR 26X41mm	07613327016802
UH1-42-22	1R2VD1	UHR BIPOLAR 22X42MM	07613327016833
UH1-42-22	2E84HD	UHR BIPOLAR 22X42MM	07613327016833
UH1-42-22	2X333W	UHR BIPOLAR 22X42MM	07613327016833
UH1-42-22	4L5YMY	UHR BIPOLAR 22X42MM	07613327016833
UH1-42-22	6J2W2P	UHR BIPOLAR 22X42MM	07613327016833
UH1-42-22	7238R2	UHR BIPOLAR 22X42MM	07613327016833
UH1-42-22	7883W1	UHR BIPOLAR 22X42MM	07613327016833
UH1-42-22	7H2L8M	UHR BIPOLAR 22X42MM	07613327016833
UH1-42-22	8K315D	UHR BIPOLAR 22X42MM	07613327016833
UH1-42-22	AK2P91	UHR BIPOLAR 22X42MM	07613327016833
UH1-42-22	ET3HWR	UHR BIPOLAR 22X42MM	07613327016833
UH1-42-22	L16J85	UHR BIPOLAR 22X42MM	07613327016833
UH1-42-22	L93PY2	UHR BIPOLAR 22X42MM	07613327016833

UH1-42-22	ME1485	UHR BIPOLAR 22X42MM	07613327016833
UH1-42-22	MY1N20	UHR BIPOLAR 22X42MM	07613327016833
UH1-42-22	P987D0	UHR BIPOLAR 22X42MM	07613327016833
UH1-42-22	R92PDW	UHR BIPOLAR 22X42MM	07613327016833
UH1-42-22	TM0KR0	UHR BIPOLAR 22X42MM	07613327016833
UH1-42-22	X884ET	UHR BIPOLAR 22X42MM	07613327016833
UH1-42-22	XD3961	UHR BIPOLAR 22X42MM	07613327016833
UH1-42-22	YK39TT	UHR BIPOLAR 22X42MM	07613327016833
UH1-42-26	186JE6	UHR BIPOLAR 26X42mm	07613327016826
UH1-42-26	5Y1RL7	UHR BIPOLAR 26X42mm	07613327016826
UH1-42-26	618AMX	UHR BIPOLAR 26X42mm	07613327016826
UH1-42-26	7J0XET	UHR BIPOLAR 26X42mm	07613327016826
UH1-42-26	MLKT5A	UHR BIPOLAR 26X42mm	07613327016826
UH1-42-26	MLL022	UHR BIPOLAR 26X42mm	07613327016826
UH1-42-26	MLL05K	UHR BIPOLAR 26X42mm	07613327016826
UH1-42-26	R50NW2	UHR BIPOLAR 26X42mm	07613327016826
UH1-42-28	1R5MDV	UHR BIPOLAR 28X42mm	07613327016857
UH1-42-28	275NEA	UHR BIPOLAR 28X42mm	07613327016857
UH1-42-28	2V6KMY	UHR BIPOLAR 28X42mm	07613327016857
UH1-42-28	4215LV	UHR BIPOLAR 28X42mm	07613327016857
UH1-42-28	MLKEHT	UHR BIPOLAR 28X42mm	07613327016857
UH1-42-28	MLLJ5X	UHR BIPOLAR 28X42mm	07613327016857
UH1-42-28	MLLJK9	UHR BIPOLAR 28X42mm	07613327016857
UH1-42-28	MLLV75	UHR BIPOLAR 28X42mm	07613327016857
UH1-43-22	11138V	UHR BIPOLAR 22X43MM	07613327016840
UH1-43-22	23005R	UHR BIPOLAR 22X43MM	07613327016840
UH1-43-22	481KR9	UHR BIPOLAR 22X43MM	07613327016840
UH1-43-22	516JLV	UHR BIPOLAR 22X43MM	07613327016840
UH1-43-22	761T5M	UHR BIPOLAR 22X43MM	07613327016840
UH1-43-22	9W3NWW	UHR BIPOLAR 22X43MM	07613327016840
UH1-43-22	DA0W51	UHR BIPOLAR 22X43MM	07613327016840
UH1-43-22	EK6069	UHR BIPOLAR 22X43MM	07613327016840
UH1-43-22	KP0AAD	UHR BIPOLAR 22X43MM	07613327016840
UH1-43-22	MD1Y5V	UHR BIPOLAR 22X43MM	07613327016840
UH1-43-22	NT68KL	UHR BIPOLAR 22X43MM	07613327016840
UH1-43-22	RT8N8H	UHR BIPOLAR 22X43MM	07613327016840
UH1-43-22	RY828P	UHR BIPOLAR 22X43MM	07613327016840
UH1-43-26	192J0R	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	1D45L5	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	1H5V53	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	2P2A99	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	2V6D0R	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	3248V2	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	334851	UHR BIPOLAR 26X43mm	07613327016864

UH1-43-26	386P0V	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	3N33VT	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	402RP3	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	4644KD	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	4J5Y6H	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	4L3JPP	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	5H1A45	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	5Y2422	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	626J60	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	6E33NR	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	6J30X2	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	6M2KK8	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	6V092N	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	761EJP	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	7820L1	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	791NPH	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	7D0H4K	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	7K28Y8	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	7K63EE	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	7N7WX1	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	817AJJ	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	8J6TAL	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	8L6R73	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	8T3XDK	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	983JDR	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	998H5N	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	9E631T	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	A93KWD	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	AY40L3	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	DJ105H	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	E521L6	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	EX5DA9	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	HW7XYK	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	J94M7R	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	JT565R	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	JX0L8R	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	K939KP	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	KA30HV	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	KE5770	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	KP1Y0E	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	KV3L3H	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	MD76WR	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	MJ1DY9	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	MK2MN3	UHR BIPOLAR 26X43mm	07613327016864

UH1-43-26	MLL54K	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	MLL65N	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	MLLV64	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	MR6D42	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	MW2J8W	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	N314DK	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	NP48X0	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	P147V0	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	P502R8	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	P92DV1	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	RL1TL1	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	RN8NM1	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	T66JK4	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	TA20N9	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	TT3EJR	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	TV197Y	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	TY7YV5	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	V68NXH	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	V91EK1	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	VW74NJ	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	VX5T02	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	VY0525	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	WA0VRE	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	WT7P1A	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	X48P1Y	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	Y63LYJ	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	Y648WL	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	YE3W2X	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	YJ6PHM	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	YL3HVV	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	YL6L7K	UHR BIPOLAR 26X43mm	07613327016864
UH1-43-26	YX8NTJ	UHR BIPOLAR 26X43mm	07613327016864
UH1-44-26	336652	UHR BIPOLAR 26X44mm	07613327016888
UH1-44-26	581XTY	UHR BIPOLAR 26X44mm	07613327016888
UH1-44-26	5H1YRH	UHR BIPOLAR 26X44mm	07613327016888
UH1-44-26	640Y4R	UHR BIPOLAR 26X44mm	07613327016888
UH1-44-26	7Y7H14	UHR BIPOLAR 26X44mm	07613327016888
UH1-44-26	9J7T4L	UHR BIPOLAR 26X44mm	07613327016888
UH1-44-26	A881HX	UHR BIPOLAR 26X44mm	07613327016888
UH1-44-26	JN7TXM	UHR BIPOLAR 26X44mm	07613327016888
UH1-44-26	KT4P91	UHR BIPOLAR 26X44mm	07613327016888
UH1-44-26	LJ44JX	UHR BIPOLAR 26X44mm	07613327016888
UH1-44-26	M52APN	UHR BIPOLAR 26X44mm	07613327016888
UH1-44-26	MLL2PT	UHR BIPOLAR 26X44mm	07613327016888

UH1-44-26	MLL5X1	UHR BIPOLAR 26X44mm	07613327016888
UH1-44-26	MLLD4R	UHR BIPOLAR 26X44mm	07613327016888
UH1-44-26	MLLPK1	UHR BIPOLAR 26X44mm	07613327016888
UH1-44-26	MLLV2Y	UHR BIPOLAR 26X44mm	07613327016888
UH1-44-26	NR773K	UHR BIPOLAR 26X44mm	07613327016888
UH1-44-26	NV0NVX	UHR BIPOLAR 26X44mm	07613327016888
UH1-44-26	PD7P36	UHR BIPOLAR 26X44mm	07613327016888
UH1-44-26	R74D1D	UHR BIPOLAR 26X44mm	07613327016888
UH1-44-26	T31P16	UHR BIPOLAR 26X44mm	07613327016888
UH1-44-26	VK7JL5	UHR BIPOLAR 26X44mm	07613327016888
UH1-44-26	W91ENE	UHR BIPOLAR 26X44mm	07613327016888
UH1-44-26	XD45XY	UHR BIPOLAR 26X44mm	07613327016888
UH1-44-26	XN8EWL	UHR BIPOLAR 26X44mm	07613327016888
UH1-44-26	Y385N4	UHR BIPOLAR 26X44mm	07613327016888
UH1-44-26	YN1AJN	UHR BIPOLAR 26X44mm	07613327016888
UH1-44-28	106X12	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	121A72	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	123MX4	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	138KNV	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	1401YE	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	1429X2	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	154N34	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	1740KE	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	1830AR	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	186WEL	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	1D8LYA	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	1E4HK2	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	1E6W32	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	1E7RRN	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	1H1P69	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	1K10H5	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	1K58A6	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	1K63LK	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	1M1W89	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	1M6JWR	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	1N6VPP	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	1R8333	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	1V4LXR	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	1X8R7X	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	1Y2V83	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	1Y7DPY	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	201MY5	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	2116XX	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	2286V1	UHR BIPOLAR 28X44mm	07613327016895

UH1-44-28	25090N	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	251WMR	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	2550NN	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	261Y8Y	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	287N2R	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	2D4N58	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	2X5T7Y	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	335351	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	335EXD	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	391JXR	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	3K6YR2	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	444636	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	4E45R8	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	4K15YY	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	4L17DL	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	4L29A4	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	4M7J5A	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	4N2091	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	4P3L83	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	4T85VA	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	4W227H	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	518WTR	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	521K2M	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	5246AR	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	5325D6	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	534JYN	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	555249	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	558TVA	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	581YM9	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	5D1Y5V	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	5H0EXH	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	5K0AMY	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	5P0H6D	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	5R3H24	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	5V2RJ0	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	600TDM	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	623LA2	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	6370A6	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	681Y0M	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	6885NJ	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	693MRX	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	698Y2M	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	6L4LT7	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	6N4J52	UHR BIPOLAR 28X44mm	07613327016895

UH1-44-28	6T3JAY	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	6V19AE	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	6V6DW2	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	6W0A59	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	7307MN	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	758K2A	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	768V6X	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	7D5MN2	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	7L230N	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	7M6YVP	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	7X2PH3	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	7X7XM0	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	8125RN	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	867YNX	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	886V7A	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	8E6V9P	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	8M114H	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	8M4W64	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	8R438X	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	8X4PRJ	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	971KAN	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	987WJ0	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	9H673H	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	9K656K	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	9V0L0L	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	A56J67	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	A85T4Y	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	A86V5E	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	AA7VL4	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	AM3PPW	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	AW6JML	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	DD37HD	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	DH5N6T	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	DL5E4X	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	DR69YE	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	E47VAT	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	EA6WY1	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	EK5N6A	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	EK8KDD	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	EN17M3	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	ER2VAV	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	HW14NW	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	J477D9	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	J63M24	UHR BIPOLAR 28X44mm	07613327016895

UH1-44-28	J644K6	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	JE7K73	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	JT2PYD	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	K47JDR	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	KD174L	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	L548RP	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	LE83EE	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	LH10XA	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	LM74JP	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	LR3KT5	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	LV0PP9	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	M45LKM	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	MLKMWA	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	MLL0X7	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	MLL3T3	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	MLL3X8	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	MLL3YN	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	MLL5V9	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	MLL6E1	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	MLL70V	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	MLLA90	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	MLLDXM	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	MLLECY	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	MLLLKA	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	MLLNDM	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	MLLR0R	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	MLLY2V	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	MT2YK4	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	MT56MX	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	N082LN	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	N35W9W	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	N4476X	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	N60KXY	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	N780A7	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	N80J7M	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	NR5X7N	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	NW7H6M	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	PN3K4T	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	PV87N1	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	R20NVJ	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	RM8HJJ	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	RY4JDM	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	T6297W	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	T73HN0	UHR BIPOLAR 28X44mm	07613327016895

UH1-44-28	TP366X	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	V281XV	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	V64XWE	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	V95MVD	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	VA7769	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	VL3RJJ	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	VP8E88	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	W33XH8	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	W969RX	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	WD03AK	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	WJ4YYD	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	WK2AHL	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	WW6MVP	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	X16M08	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	XM319D	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	XN437D	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	Y703W8	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	YJ012E	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	YV08TH	UHR BIPOLAR 28X44mm	07613327016895
UH1-44-28	YY4TKX	UHR BIPOLAR 28X44mm	07613327016895
UH1-45-22	8Y0KWH	UHR BIPOLAR 22X45MM	07613327016901
UH1-45-22	HR869M	UHR BIPOLAR 22X45MM	07613327016901
UH1-45-22	R7168K	UHR BIPOLAR 22X45MM	07613327016901
UH1-45-22	WR4MT4	UHR BIPOLAR 22X45MM	07613327016901
UH1-45-26	056MNE	UHR BIPOLAR 26X45mm	07613327016925
UH1-45-26	059MNE	UHR BIPOLAR 26X45mm	07613327016925
UH1-45-26	188AWK	UHR BIPOLAR 26X45mm	07613327016925
UH1-45-26	243WTY	UHR BIPOLAR 26X45mm	07613327016925
UH1-45-26	352XY0	UHR BIPOLAR 26X45mm	07613327016925
UH1-45-26	442K1P	UHR BIPOLAR 26X45mm	07613327016925
UH1-45-26	8L58NK	UHR BIPOLAR 26X45mm	07613327016925
UH1-45-26	AJ0H2R	UHR BIPOLAR 26X45mm	07613327016925
UH1-45-26	AN15H5	UHR BIPOLAR 26X45mm	07613327016925
UH1-45-26	HX6PAR	UHR BIPOLAR 26X45mm	07613327016925
UH1-45-26	HY2V81	UHR BIPOLAR 26X45mm	07613327016925
UH1-45-26	L50T23	UHR BIPOLAR 26X45mm	07613327016925
UH1-45-26	MLL0AE	UHR BIPOLAR 26X45mm	07613327016925
UH1-45-26	MLL5J4	UHR BIPOLAR 26X45mm	07613327016925
UH1-45-26	MLL9TV	UHR BIPOLAR 26X45mm	07613327016925
UH1-45-26	MLLKNT	UHR BIPOLAR 26X45mm	07613327016925
UH1-45-26	N88PKJ	UHR BIPOLAR 26X45mm	07613327016925
UH1-45-26	P31A9Y	UHR BIPOLAR 26X45mm	07613327016925
UH1-45-26	RJ8JPR	UHR BIPOLAR 26X45mm	07613327016925
UH1-45-28	1A3013	UHR BIPOLAR 28X45mm	07613327016918

UH1-45-28	360RPN	UHR BIPOLAR 28X45mm	07613327016918
UH1-45-28	486EJ1	UHR BIPOLAR 28X45mm	07613327016918
UH1-45-28	4N85AM	UHR BIPOLAR 28X45mm	07613327016918
UH1-45-28	D6348Y	UHR BIPOLAR 28X45mm	07613327016918
UH1-45-28	EY829J	UHR BIPOLAR 28X45mm	07613327016918
UH1-45-28	MLL1V3	UHR BIPOLAR 28X45mm	07613327016918
UH1-45-28	MLL210	UHR BIPOLAR 28X45mm	07613327016918
UH1-45-28	MLL24J	UHR BIPOLAR 28X45mm	07613327016918
UH1-45-28	MLL44A	UHR BIPOLAR 28X45mm	07613327016918
UH1-45-28	MLL5RM	UHR BIPOLAR 28X45mm	07613327016918
UH1-45-28	MLL94Y	UHR BIPOLAR 28X45mm	07613327016918
UH1-45-28	MLLDXY	UHR BIPOLAR 28X45mm	07613327016918
UH1-45-28	MLLJRN	UHR BIPOLAR 28X45mm	07613327016918
UH1-45-28	MLLVWK	UHR BIPOLAR 28X45mm	07613327016918
UH1-45-28	MLLW3M	UHR BIPOLAR 28X45mm	07613327016918
UH1-45-28	MLLY7V	UHR BIPOLAR 28X45mm	07613327016918
UH1-45-28	X43WJN	UHR BIPOLAR 28X45mm	07613327016918
UH1-46-26	1839Y2	UHR BIPOLAR 26X46mm	07613327016949
UH1-46-26	1M2AX3	UHR BIPOLAR 26X46mm	07613327016949
UH1-46-26	1P2E82	UHR BIPOLAR 26X46mm	07613327016949
UH1-46-26	1V5AN9	UHR BIPOLAR 26X46mm	07613327016949
UH1-46-26	226AWD	UHR BIPOLAR 26X46mm	07613327016949
UH1-46-26	232X83	UHR BIPOLAR 26X46mm	07613327016949
UH1-46-26	2465J7	UHR BIPOLAR 26X46mm	07613327016949
UH1-46-26	2K0P93	UHR BIPOLAR 26X46mm	07613327016949
UH1-46-26	2R3XD6	UHR BIPOLAR 26X46mm	07613327016949
UH1-46-26	3614YL	UHR BIPOLAR 26X46mm	07613327016949
UH1-46-26	394D8M	UHR BIPOLAR 26X46mm	07613327016949
UH1-46-26	3P7091	UHR BIPOLAR 26X46mm	07613327016949
UH1-46-26	4V2PVP	UHR BIPOLAR 26X46mm	07613327016949
UH1-46-26	5389M0	UHR BIPOLAR 26X46mm	07613327016949
UH1-46-26	5N2N75	UHR BIPOLAR 26X46mm	07613327016949
UH1-46-26	9X7YK5	UHR BIPOLAR 26X46mm	07613327016949
UH1-46-26	A35V7K	UHR BIPOLAR 26X46mm	07613327016949
UH1-46-26	DW3J8W	UHR BIPOLAR 26X46mm	07613327016949
UH1-46-26	DY85VA	UHR BIPOLAR 26X46mm	07613327016949
UH1-46-26	E242NY	UHR BIPOLAR 26X46mm	07613327016949
UH1-46-26	H07J22	UHR BIPOLAR 26X46mm	07613327016949
UH1-46-26	JV88X1	UHR BIPOLAR 26X46mm	07613327016949
UH1-46-26	L04PPV	UHR BIPOLAR 26X46mm	07613327016949
UH1-46-26	LY44K5	UHR BIPOLAR 26X46mm	07613327016949
UH1-46-26	MLL50M	UHR BIPOLAR 26X46mm	07613327016949
UH1-46-26	MLL789	UHR BIPOLAR 26X46mm	07613327016949
UH1-46-26	MLLA53	UHR BIPOLAR 26X46mm	07613327016949

UH1-46-26	NE45DD	UHR BIPOLAR 26X46mm	07613327016949
UH1-46-26	V94NTT	UHR BIPOLAR 26X46mm	07613327016949
UH1-46-26	YY5N98	UHR BIPOLAR 26X46mm	07613327016949
UH1-46-28	1804DD	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	29772K	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	2E5PX4	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	3402DJ	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	3P0NDP	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	4H0W7N	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	4J5R3M	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	4K0NDA	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	5818XA	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	5R5H8L	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	6N77PT	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	6T493L	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	7K0L36	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	7L1KRT	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	8015WY	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	802XY1	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	8P7D00	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	905HH2	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	AD0E56	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	D3227P	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	E83XET	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	EV4W50	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	HT21XV	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	JP0TKR	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	L00E8N	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	M25670	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	MLL0Y0	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	MLL1W9	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	MLL46R	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	MLL499	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	MLL602	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	MLL8Y7	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	MLLAYE	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	MLLDY9	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	MLLJWN	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	MLLK40	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	MLLKEH	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	MLLNE3	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	MLLRL2	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	MLLV8T	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	MLLVTD	UHR BIPOLAR 28X46mm	07613327016956

UH1-46-28	MLLXV3	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	MR613J	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	NT30DY	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	P72N17	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	VE411K	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	WH517T	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	XH01PP	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	Y66HVE	UHR BIPOLAR 28X46mm	07613327016956
UH1-46-28	YE5XLP	UHR BIPOLAR 28X46mm	07613327016956
UH1-47-26	1Y2K34	UHR BIPOLAR 26X47mm	07613327016987
UH1-47-26	21197T	UHR BIPOLAR 26X47mm	07613327016987
UH1-47-26	3E0A48	UHR BIPOLAR 26X47mm	07613327016987
UH1-47-26	4E7MKW	UHR BIPOLAR 26X47mm	07613327016987
UH1-47-26	7K1VLW	UHR BIPOLAR 26X47mm	07613327016987
UH1-47-26	A58HN9	UHR BIPOLAR 26X47mm	07613327016987
UH1-47-26	J38RE2	UHR BIPOLAR 26X47mm	07613327016987
UH1-47-26	J60A4E	UHR BIPOLAR 26X47mm	07613327016987
UH1-47-26	LT43X0	UHR BIPOLAR 26X47mm	07613327016987
UH1-47-26	M15HKN	UHR BIPOLAR 26X47mm	07613327016987
UH1-47-26	MK442Y	UHR BIPOLAR 26X47mm	07613327016987
UH1-47-26	MLL2LR	UHR BIPOLAR 26X47mm	07613327016987
UH1-47-26	MLL52N	UHR BIPOLAR 26X47mm	07613327016987
UH1-47-26	MLLRMW	UHR BIPOLAR 26X47mm	07613327016987
UH1-47-26	MLLVRW	UHR BIPOLAR 26X47mm	07613327016987
UH1-47-26	MLLWVH	UHR BIPOLAR 26X47mm	07613327016987
UH1-47-26	XD724J	UHR BIPOLAR 26X47mm	07613327016987
UH1-47-28	2821VP	UHR BIPOLAR 28X47mm	07613327016970
UH1-47-28	2N68HT	UHR BIPOLAR 28X47mm	07613327016970
UH1-47-28	3682T6	UHR BIPOLAR 28X47mm	07613327016970
UH1-47-28	3M0Y7N	UHR BIPOLAR 28X47mm	07613327016970
UH1-47-28	4E58YD	UHR BIPOLAR 28X47mm	07613327016970
UH1-47-28	661DW6	UHR BIPOLAR 28X47mm	07613327016970
UH1-47-28	8K6225	UHR BIPOLAR 28X47mm	07613327016970
UH1-47-28	8M60J8	UHR BIPOLAR 28X47mm	07613327016970
UH1-47-28	9D22DE	UHR BIPOLAR 28X47mm	07613327016970
UH1-47-28	9Y0X3H	UHR BIPOLAR 28X47mm	07613327016970
UH1-47-28	A68HK5	UHR BIPOLAR 28X47mm	07613327016970
UH1-47-28	AL4LYE	UHR BIPOLAR 28X47mm	07613327016970
UH1-47-28	J60DLA	UHR BIPOLAR 28X47mm	07613327016970
UH1-47-28	KM65D9	UHR BIPOLAR 28X47mm	07613327016970
UH1-47-28	LJ3RRT	UHR BIPOLAR 28X47mm	07613327016970
UH1-47-28	MLL8V0	UHR BIPOLAR 28X47mm	07613327016970
UH1-47-28	MLL984	UHR BIPOLAR 28X47mm	07613327016970
UH1-47-28	MLLH9P	UHR BIPOLAR 28X47mm	07613327016970

UH1-47-28	MLLLLD	UHR BIPOLAR 28X47mm	07613327016970
UH1-47-28	MLLND6	UHR BIPOLAR 28X47mm	07613327016970
UH1-47-28	PV1K03	UHR BIPOLAR 28X47mm	07613327016970
UH1-47-28	PW100N	UHR BIPOLAR 28X47mm	07613327016970
UH1-47-28	TD5KRD	UHR BIPOLAR 28X47mm	07613327016970
UH1-47-28	WX0R02	UHR BIPOLAR 28X47mm	07613327016970
UH1-48-26	1J8K9L	UHR BIPOLAR 26X48mm	07613327017007
UH1-48-26	2Y6JK3	UHR BIPOLAR 26X48mm	07613327017007
UH1-48-26	350D6V	UHR BIPOLAR 26X48mm	07613327017007
UH1-48-26	3845V3	UHR BIPOLAR 26X48mm	07613327017007
UH1-48-26	3D6216	UHR BIPOLAR 26X48mm	07613327017007
UH1-48-26	44827N	UHR BIPOLAR 26X48mm	07613327017007
UH1-48-26	580Y68	UHR BIPOLAR 26X48mm	07613327017007
UH1-48-26	726W7P	UHR BIPOLAR 26X48mm	07613327017007
UH1-48-26	7L7TPP	UHR BIPOLAR 26X48mm	07613327017007
UH1-48-26	8E7A1J	UHR BIPOLAR 26X48mm	07613327017007
UH1-48-26	9V8P94	UHR BIPOLAR 26X48mm	07613327017007
UH1-48-26	9Y6HML	UHR BIPOLAR 26X48mm	07613327017007
UH1-48-26	H88E74	UHR BIPOLAR 26X48mm	07613327017007
UH1-48-26	K54MN4	UHR BIPOLAR 26X48mm	07613327017007
UH1-48-26	L783M7	UHR BIPOLAR 26X48mm	07613327017007
UH1-48-26	MLL3E9	UHR BIPOLAR 26X48mm	07613327017007
UH1-48-26	MLL52W	UHR BIPOLAR 26X48mm	07613327017007
UH1-48-26	MP56YN	UHR BIPOLAR 26X48mm	07613327017007
UH1-48-26	N00HDA	UHR BIPOLAR 26X48mm	07613327017007
UH1-48-26	N26TW7	UHR BIPOLAR 26X48mm	07613327017007
UH1-48-26	TR4TT5	UHR BIPOLAR 26X48mm	07613327017007
UH1-48-26	W016KY	UHR BIPOLAR 26X48mm	07613327017007
UH1-48-26	XJ48DN	UHR BIPOLAR 26X48mm	07613327017007
UH1-48-26	YA4YPN	UHR BIPOLAR 26X48mm	07613327017007
UH1-48-28	1831VA	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	1M4NYV	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	1T0L93	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	1W0E1J	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	254N4A	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	294M71	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	2A2RJ0	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	2T45AP	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	3K8VM4	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	3T6PD5	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	3Y8DJX	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	437ANA	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	448V4T	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	4623VM	UHR BIPOLAR 28X48mm	07613327017014

UH1-48-28	484LE9	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	4X0M53	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	5D50WT	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	5L074T	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	5T0PX8	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	6E2MX7	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	6K5Y7V	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	6L43JP	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	7A1HVV	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	8J8DDN	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	975VXR	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	9786P2	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	9D6M48	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	9M6KHT	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	9T2991	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	9X5MJ3	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	A228M1	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	A54PLW	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	A77V1M	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	AV5ANE	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	AX18R8	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	DK0M6J	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	DY3LR8	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	E600KE	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	EY5P9V	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	H558EP	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	H63W11	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	HP74LE	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	J72DLP	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	J76YJ6	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	JK7JY4	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	JL1T01	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	K47WL6	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	K728DH	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	KT3K27	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	L35VMM	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	M15D88	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	ML731E	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	MLL18D	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	MLL562	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	MLL6MH	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	MLL74R	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	MLLM4T	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	MLLMWV	UHR BIPOLAR 28X48mm	07613327017014

UH1-48-28	MLLNTV	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	MLLPV9	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	MLLV38	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	MLLX14	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	N046DW	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	N06WWP	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	N20WDM	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	N75AY1	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	P419XW	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	RP3D0X	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	TD218R	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	V15WA1	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	VH06PA	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	VN3P18	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	VR853H	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	WM8KJT	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	XE2ND8	UHR BIPOLAR 28X48mm	07613327017014
UH1-48-28	YH0N93	UHR BIPOLAR 28X48mm	07613327017014
UH1-49-32	324H9H	UHR BIPOLAR 32X49MM	07613327017069
UH1-51-32	9J6N9E	UHR BIPOLAR 32X51MM	07613327017120
UH1-52-26	5J8PM2	UHR BIPOLAR 26X52mm	07613327017168
UH1-52-28	V02MME	UHR BIPOLAR 28X52mm	07613327017182
UH1-52-32	5D4NMJ	UHR BIPOLAR 32X52MM	07613327017212
UH1-52-32	9T612E	UHR BIPOLAR 32X52MM	07613327017212
UH1-54-32	MNNAAL	UHR BIPOLAR 32X54MM	07613327017267
UH1-55-28	R30MME	UHR BIPOLAR 28X55mm	07613327017281
UH1-58-26	1R8TNH	UHR BIPOLAR 26X58mm	07613327017359
UH1-58-26	414A2L	UHR BIPOLAR 26X58mm	07613327017359
UH1-58-26	4W1HD7	UHR BIPOLAR 26X58mm	07613327017359
UH1-58-26	AV0J7J	UHR BIPOLAR 26X58mm	07613327017359
UH1-58-26	E226XA	UHR BIPOLAR 26X58mm	07613327017359
UH1-58-26	E50DAW	UHR BIPOLAR 26X58mm	07613327017359
UH1-58-26	M31665	UHR BIPOLAR 26X58mm	07613327017359
UH1-58-26	WK1WP2	UHR BIPOLAR 26X58mm	07613327017359
UH1-58-26	X90ML6	UHR BIPOLAR 26X58mm	07613327017359
78-5-0910	30582401	SCORPIO + MOBILE BEARING PCS	07613327366006
78-5-0910	30582601	SCORPIO + MOBILE BEARING PCS	07613327366006
82-3-0318	28208551	NRG KNEE P/S	04546540364319
72-4-0324	23505751	SCORPIO TS TIBIAL INSERT	07613327034356
82-2-0512	23325151	NRG KNEE CR	04546540400611
82-2-0518	23640051	NRG KNEE CR	04546540400635
82-2-0912	26862851	NRG KNEE CR	04546540400758
72-16-0512	19093651	SCORPIOFLEX TOTAL KNEE CR	07613154020515
72-16-0512	27835951	SCORPIOFLEX TOTAL KNEE CR	07613154020515

72-16-0912	17879451	SCORPIOFLEX TOTAL KNEE CR	07613154020652
72-16-0512	17376651	SCORPIOFLEX TOTAL KNEE CR	07613154020515
72-16-0512	23125951	SCORPIOFLEX TOTAL KNEE CR	07613154020515
72-16-0518	25621551	SCORPIOFLEX TOTAL KNEE CR	07613154020539
72-16-1112	22764851	SCORPIOFLEX TOTAL KNEE CR	07613154020737